

# Protocol for preparation of competent cells for transformation

For two transformations

## Materials:

- 0.1 M Calcium Chloride chilled on the ice;
- Overnight bacteria I culture or bacteria I colonies;

## Procedure:

1. Add 20  $\mu$ l of the overnight bacteria I culture or pick a colony to 1 ml of LB antibiotic liquid medium, Incubate at 37 degree in a shaker till the OD600 value reaches 0.4–0.6.
2. Put the tubes on ice to incubate for 5 min.
3. Pellet bacterial cells by 5 min centrifugation at 5000 rpm, discard the supernatant.
4. Resuspend cells in 600  $\mu$ l of ice–chilled 0.1 M Calcium Chloride solution. Incubate on ice for 30 min.
5. Centrifuge for 5 min at 5000 rpm in a microcentrifuge tube , discard the supernatant.
6. Resuspend the pelleted cells in 100  $\mu$ l of ice–chilled 0.1 M Calcium Chloride solution. Incubate on ice.
7. Add 50  $\mu$ l of the prepared cells to each tube containing DNA sample , mix and incubate on ice for 30 min.
8. Transform subsequently as the transformation protocol.

**Note:**

1. Make sure the cells are not left in the centrifuge at ambient temperature for more than 5 min as this will significantly decrease the transformation efficiency.
2. The rpm at centrifugation is not higher than 5000, as a high rpm may cause the lysis of cells.