

# PCR with iPROOF High Fidelity DNA Polymerase (BIORAD)

Nadine and Effie

July 27, 2010

---

<b>Where</b>	On the bench, PCR machine on the bench left to the door	
	<b>Quantity</b>	<b>What</b>
<b>Material needed</b>	for 50 $\mu$ l of PCR mixture 10 $\mu$ l	5x iProof HF buffer
	1 $\mu$ l	of each primer
	1 $\mu$ l	dNTP mix of 10mM (Roche tube in the -20° freezer with purple cap)
	0.5 $\mu$ l	High Fidelity Polymerase
	0.5 $\mu$ l	template
		calculate dH <sub>2</sub> O quantity needed to have 50 $\mu$ l
<b>Steps</b>	<ol style="list-style-type: none"><li>1. Add everything together.</li><li>2. Put in the PCR thermal cycler (PCR machine), Place carefully and use green plastic rack if required.</li><li>3. Program the PCR thermal cycler: <b>1.</b> 98°C for 30s <b>2.</b> 98°C for 7.5 seconds <b>3.</b> 60°C for 20 sec (this depends on the primers, therefore look up the annealing temperature of the primers used) <b>4.</b> 72°C you have to calculate the time required. This polymerase has an efficiency of 15-30s/kb. (calculate for the longest fragment.) <b>6.</b> 72°C for 5 min <b>7.</b> 4°C "forever"</li><li>4. Save the file in the MAIN folder → OK→ SAVE→RUN → Select the block you use → RUN → Change sample volume to 50<math>\mu</math>l, LidT should be 105°C and the box just below ("hotlid") should ALWAYS be checked</li><li>5. You can use STATUS to check the progress of the PCR.</li></ol>	
<b>Warnings</b>	<ul style="list-style-type: none"><li>• Sign in for PCR on the White Board just above the machine (write iGEM and from when to when you will use it.)</li><li>• It is possible to save polymerase and only do a 10 <math>\mu</math>l mix if you only want to check if your insert is there and don't actually want to amplify it.</li></ul>	

---