

# Transformation protocol

## Materials:

- Plasmid samples or ligation product;
- Commercially competent cells;
- LB non-antibiotic liquid medium;
- LB antibiotic agar plates

## Procedure:

1. Get the competent cells from  $-70$  degree, and wait for its fusion.  $30-50 \mu$ l of competent *E.coli* cells for each sample. Put microcentrifuge tubes to chill on ice for at least 2 min.
2. Add 2 – 3  $\mu$ l of each plasmid sample or all the ligation product into the competent cells in the microcentrifuge tubes. Mix and incubate on ice for 30 min.
3. Heat pulse for 90 sec, at  $42$  degree. Put back to ice and incubate for 5 min.
4. Add 200  $\mu$ L LB non-antibiotic liquid medium into each microcentrifuge tube. Shake the microcentrifuge tubes in shaker, at  $37$  degree, for 30 min to recover.
5. Plate 150  $\mu$ L of the liquid medium with transformed cells immediately, on prewarmed LB antibiotic agar plates. Incubate overnight at  $37^{\circ}\text{C}$  for 10–14h.

## Tips:

- All procedures are performed on ice.
- Make sure the cells are not left at ambient temperature for more than 5

min as this will significantly decrease the transformation efficiency.

■ When got out from the shaker, the competent cells may form pellet in the microcentrifuge tubes. You need to resuspend the cells before plating.

#### **References:**

\*Current protocols in molecular biology.