

# 1 DNA double digestion protocol

## Materials:

- DNA sample(s) in water or TE buffer
- 10x digestion buffer
- Restriction enzyme s (EcoRI or SpeI or XbaI or PstI)
- DNA loading buffer (if electrophoresis is subsequent)
- Agarose gel 1.5% (or different depending on expected band sizes)

## Procedure:

1. Test the concentration of the DNA sample(s).

2. Pipet the following into a microfuge tube:

20uL reaction system	50uL reaction system
DNA around 1ug	around 2.5ug
10x Digestion buffer 2uL	5uL
1 <sup>st</sup> Enzyme 1–1.5uL	2.5–4uL
2 <sup>nd</sup> Enzyme 1–1.5uL	2.5–4uL
ddWater Rest of volume	Rest of volume

3. Incubate at recommended temperature (37.0 degrees) for 2 or 4 hours

(0.5~2h for enzymes of NEB, 4h for enzymes of Takara).

4. Take 2 to 5 uL of the digested sample, add loading buffer, and run it on the agarose gel to check the result, or take the entire sample to run to extract a wanted fragment).

## Tips:

1. DNA:

■ For identification of DNA, use 0.4 ug/uL DNA; (or 2uL from a nice DNA mini prep)

■ For cloning, 1ug/uL DNA is enough.

2. Buffer: we' d better use the buffer that comes with the enzyme, which means buffers from other company may cause some abnormal results.

3. Enzyme: the maximum volume that an enzyme can be used is 1/10 of the total reaction volume (example: 2 uL for 20 uL reaction system). If you want to do overnight digestion, add less enzyme (example: 1 uL for 20 uL reaction system). It is necessary to point that too many enzymes will reduce the efficiency of enzyme digestion with glycerol in it.

4. Gel: make sure to run the uncut DNA as a control along with the digested DNA sample(s). And, always run a DNA marker!

**References:**

\*Current protocols in molecular biology (3.1.1–3.1.2)