

### **Cytosol MBP Sample Preparation**

1. Add 50  $\mu$ l of the overnight bacterial culture to 50 ml of LB medium-containing ampicillin (100 g/ml). Cultures were grown at 37 °C until the A600 reached 0.6 and then induced with 0.5 mM IPTG at 30 °C for 3 h.
2. The cells were harvested by centrifugation at 6000 rpm for 10 min. The pelleted cells were suspended in 5 ml buffer A (40 mM Tris-HCl pH 8.0, 100 mM NaCl). Then the cell were lysed by ultrasonic, and the cell lysate was separated by centrifugation at 13000 rpm for 15 min.
3. Add 50  $\mu$ l SDS-PAGE loading buffer into 50  $\mu$ l supernatant, then heat it at 95°C for 10 min.
4. The pellet were suspended in 5 ml buffer A. Add 50  $\mu$ l SDS-PAGE loading buffer into 50  $\mu$ l suspension, then heat it at 95°C for 30 min.

### **Periplasmic MBP Sample Preparation**

1. Add 50  $\mu$ l of the overnight bacterial culture to 50 ml of LB medium-containing Kanamycin (40 g/ml). Cultures were grown at 37 °C until the A600 reached 0.6 and then induced with 0.5 mM IPTG at 20 °C for 6 h.
2. The cell of 1 ml culture were harvested by centrifugation at 6000 rpm for 10 min. The pelleted cells were suspended in 100  $\mu$ l buffer B (1 mg/ml lysozyme, 20%(m/v) sucrose, 30 mM Tris-HCl pH8.0, 1 mM EDTA), ice-bath for 20 min. The mixture were separated by centrifugation at 13000 rpm for 1 min.
3. Add 50  $\mu$ l SDS-PAGE loading buffer into 50  $\mu$ l supernatant, then heat it at 95°C for 10 min.
4. The pellet were suspended in 5 ml buffer A. Add 50  $\mu$ l SDS-PAGE loading buffer into 50  $\mu$ l suspension, then heat it at 95°C for 30 min.