



Lab Diary for Gyrase Binding Site

Date: Tuesday, 20 August 2011 11:00

Topic: Re-suspension of primers and PCR reaction for the 3 GBS

- 1) Carried out re-suspension of primers (forward and reverse primers of mu, pSC101 and pBR322).
- 2) Conducted a PCR of gene MU, pSC101 and pBR 322.

Materials	Volume (μ l)
water	32.5
5X HF buffer	10.0
dNTPs	1.0
Forward primer	2.5
Reverse primer	2.5
Template	1.0
Phusion Enzyme	0.5
Total	50

PCR setting:

Step 1- 1 cycle

98°C for 30 sec

Step 2 – 30 cycle

98 °C for 10 sec

65 °C for 30 sec

72 °C for 15 sec


Step 3- 1 cycle

72 °C for 10 min

Step 4

4 °C – hold

3) Performed gel electrophoresis



Date: Wednesday, 21 August 2011 11:00

Topic: PCR Purification, Digestion, Ligation and Transformation.

1) Carried out PCR purification and recovered 50 μ l of PCR product.

2) Set up restriction digestion with PCR purification product.:

Materials	CM Backbone(μ l)	pSC101 (μ l)	pBR322 (μ l)	MU (μ l)
DNA	20	20	20	20
EcoRI	0.5	0.5	0.5	0.5
PstI	0.5	0.5	0.5	0.5
DpnI	0.5	-	-	-
Buffer 2	5	5	5	5
BSA	0.5	0.5	0.5	0.5
water	23	23.5	23.5	23.5
Total	50	50	50	50


3) Set up the following ligation mixtures:

Material	A1 (μ l)	A2 (μ l)	A3 (μ l)
CM backbone	4.5	4.5	4.5
MU	4.0		
pBR322		4.0	
pSC101			4.0
Ligase	0.5	0.5	0.5
Ligase Buffer	1.0	1.0	1.0
Total	10	10	10

4) Transformation

Organism: TOP 10 *E. Coli*

- 1) Carried out transformation for Gyrase binding sites MU, pSC101 and pBR322.
- 2) Plated the transformed cells on the chloramphenicol agar plates.
- 3) Observed colonies in all the three plates.



Date: Thursday, 22 August 2011 11:00

Topic: Setup Overnight Cultures

- 1) Set up overnight cultures from a single colony from each of the transformed cell plates in selective LB medium. Left the cultures to grow overnight inside the shaker.


Date: Friday, 23 August 2011 11:00

Topic: Preparing Glycerol stock, Miniprep, Digestion, Gel.

- 1) Set up glycerol stock from the overnight culture and stored it at -80°C .
- 2) Mini-prepped the overnight culture into 10 μl volume.
- 3) Digested 7.5 μl of the mini-prep sample with EcoRI and PstI. Ran the digestion mixture on agarose gel by electrophoreses. The bands confirmed the size of the final construct and plasmid backbone.

Material	B1 (MU) (μl)	B2(pBR322) (μl)	B3 (pSB101) (μl)
Miniprep sample	7.5	7.5	7.5
EcoRI	0.5	0.5	0.5
PstI	0.5	0.5	0.5
Buffer2	1.0	1.0	1.0
BSA	0.5	0.5	0.5
Total	10	10	10

4) Performed gel electrophoresis for Digestion and mini-prep samples.



min- miniprep sample
D- Digestion sample

Date: Saturday, 24 August 2011 11:00

Topic: Sequencing.

1) Sent 10 μ l of mini-prep sample to Wolfson Institute for Biomedical Research for sequencing of all three GBS.

Digestion and ligation had failed, the experiment was repeated using new digestion and ligation protocols.

Experiment continued with three GBS (MU, pBR 322, pSC 101) and two backbones (pSB1C3, pSB3T5).

Date: Thursday, 27 August 2011 11:00

Topic: Digestion

1) Set up restriction digestion of pSB1C3:

Material	E13 (μ l)
pSB1C3	4.5
EcoRI	1.0
PstI	1.0
Buffer3	5.0
BSA	0.5

dH2O	38
Total	50


Date: Sunday, 28 August 2011 16:59

Topic: digestion and ligation

1) Set up restriction digestion for MU, pBR322, pSC101 and pSB3T5:

Materials	E17 (μ l)	E18(μ l)	E19(μ l)	E20(μ l)
Mu	20			
pBR322		19		
pSC101			22	
pSB3T5				3.5
EcoRI	1	1	1	1
PstI	1	1	1	1
Buffer 3	5	5	5	5
BSA	0.5	0.5	0.5	0.5
Water	22.5	23.5	20.5	39.0
Total	50	50	50	50

2) carried out gel electrophoresis of digestion product and all were successful.



3) Set up the following ligation:

Materials	M5(μ l)	M6(μ l)	M7(μ l)	M8(μ l)	M9(μ l)	M10(μ l)
Mu (E17)	6	6				
pBR322 (E18)			6	6		
pSC101 (E19)					6	6
pSB1C3 (E13)	2		2		2	
pSB3T5 (E20)		2		2		2
Ligase	1	1	1	1	1	1
Ligase buffer	2	2	2	2	2	2
Water	9	9	9	9	9	9
Total	20	20	20	20	20	20

4) Transformation

Transformed a stock of competent E coli with the ligation mixture M5, M6, M7, M8, M9, M10


Plated M5, M7, M9 on chloramphenicol plates

Plated M6, M8 & M10 on tetracycline plates

Date: Monday, 29 August 2011 18:00

Topic: Overnight cultures

1) Observed colonies in all the six plates.



2) Set up overnight cultures of five colonies, from the chloramphenicol plate and tetracycline plates in selective LB medium. Left the culture to grow overnight inside the shaker.

Date: Tuesday, 30 August 2011 11:00


Topic: Mini-preps and competent cells

1) Set up glycerol stock from the chloramphenicol overnight culture and stored it at -80°C, Forgot to set up glycerol stock from the tetracycline overnight cultures.

2) Mini-prepped the overnight culture into 50 µl volume and measured concentration using nano-drop.

Table showing nano-drop readings.

Sample µl	1 nanodrop (ng)	2 nanodrop (ng)	Mean nanodrop (ng)	Volume of DNA for Digestion	Plasmid	BioBrick
M5 mini-prep (1)	3.3	3.2	3.3	154	Mu/ pSB1C3	Mu
M5 mini-prep (2)	46.4	42.7	44.6	11		
M5 mini-prep (3)	30.7	30.7	30.7	16		
M5 mini-prep (4)	25.6	30.7	28.2	18		
M5 mini-prep (5)	80.0	75.4	77.7	6		
M6 mini-prep (1)	45.9	41.8	43.9	11	Mu/ pSB3T5	
M6 mini-prep (2)	70.4	85.4	77.9	6		
M6 mini-prep (3)	94.0	99.7	96.9	5		



M6 mini-prep (4)	83.7	89.6	86.7	6	pBR322/ pSB1C3	pBR322
M6 mini-prep (5)	72.1	75.5	73.8	7		
M7 mini-prep (1)	36.3	37.2	36.8	14		
M7 mini-prep (2)	37.2	37.8	37.5	13		
M7 mini-prep (3)	102.4	96.5	99.5	5		
M7 mini-prep (4)	34.9	36.6	35.8	14		
M7 mini-prep (5)	38.6	34.1	36.4	14	pBR322/ pSB3T5	
M8 mini-prep (1)	144.3	157.0	150.7	3		
M8 mini-prep (2)	75.5	75.8	75.7	7		
M8 mini-prep (3)	209.0	208.9	209.0	2		
M8 mini-prep (4)	85.7	87.4	86.6	6		
M8 mini-prep (5)	105.8	105.3	105.6	5	pSC101/ pSB1C3	pSC101
M9 mini-prep (1)	168.9	132.9	150.9	3		
M9 mini-prep (2)	15.0	16.0	15.5	32		
M9 mini-prep (3)	54.4	56.0	55.2	9		
M9 mini-prep (4)	42.9	44.8	43.9	11		
M9 mini-prep (5)	97.8	97.8	97.8	5	pSC101/ pSB3T5	
M10 mini-prep (1)	74.3	96.3	85.3	6		
M10 mini-prep (2)	71.2	72.5	71.9	7		
M10 mini-prep (3)	84.2	82.5	83.4	6		
M10 mini-prep (4)	nil	nil	76.3	7		
M10 mini-prep (5)	70.1	59.2	64.7	8		



Date: Wednesday, 31 August 2011 12:00

Topic: Digestion

1) Carried out digestion of all the mini-preps and performed gel electrophoresis

Sample μ l	E22	E23	E24	E25	E26	E27	E28	E29	E30	E31	E32	E33	E34	E35	E36
M5 mini-prep (1)															
M5 mini-prep (2)		11													
M5 mini-prep (3)			16												
M5 mini-prep (4)				18											
M5 mini-prep (5)					6										
M6 mini-prep (1)						11									
M6 mini-prep (2)							6								
M6 mini-prep (3)								5							
M6 mini-prep (4)									6						
M6 mini-prep (5)										7					
M7 mini-prep (1)											14				
M7 mini-prep												13			

2) Performed gel electrophoresis



3) Repeated gel for digested mini-preps for M9, M10.

Date: Thursday, 1 September 2011 11:00

Topic: Sent for sequencing

Sent 10 µl of mini-prep sample to Wolfson Institute for Biomedical Research for sequencing. Sent the following mini-prep sample for sequencing: M5 (3)

Mu/pSB1C3; M7 (3) **Mu/ pSB3T5**; M9 (2) **Mu/ pSB3T5**.

Sequencing Results:

Mu- success

Mu

```

TTGTTTCAGACGCTGAATAACAGCCGGGAGCGTCAGCTTGCTTTTTCCCTCTGGC
TAAACGCCTGGCAGGATCACGTCAGACCGGTCATGAAATCCGGGCGCTGGGG
TTTATGACGGGATCCGTGGCTGAAGCCGTATATCAGCGCACCGGGAATATGCCA
GCCCGTCTGCTTGTAATGAACGGGAAAAGTCTGGCGACCACTGCTGATGCCGC
CCTGAAACCGGAGGATTTACAGCGTCTGCCGTCGCTGATGGCAAACACAGG
CGGTATTGTGGGACAGGGA
    
```

The sequencing data was analysed with 'ApE' software and all the parts for the construct are present.

PBR3222- rev primer wrong.

pSC101- rev primer wrong.

Construct a plasmid with spy device with MU GBS.

Date: Friday, 02/September 2011 11:00


Topic: Digestion and ligation for MU and SPY.

1) Set up restriction digestion for Spy and mu

Spy was restricted with XbaI and PstI, MU was digested with ECOR1 and SpeI

Materials	E52	E53
SPY-GFP (mini2)	1.5	
Mu		13
ECORI		1
XbaI	1	
PstI	1	
SpeI		1
Buffer	5	5
BSA	0.5	0.5
Water	41.0	29.5
Total	50	50

2) Carried out gel electrophoresis of digestion product and all were successful.





3) Set up the following ligation mixtures

Materials	M12
SPY-GFP (E52)	5
Mu(E53)	5
pSB1A3	3
Ligase	1
Ligase Buffer	2
Water	4
Total	20

4) Transformation

Carried out transformation of competent cell stock with Spy/MU miniprep sample. Plated cells for Spy/MU on ampiciline plates.

Date: Saturday, 03/September 2011 11:00

Topic: Overnight cultures for mu/spy, spy and overnight shake flask cultures of 3 GBS for 1D chloroquine gel.

1) Seed stock: Added 10 ul of spy/GFP glycerol stock from 2009 to 10 ml LB medium falcon tubes.

Incubated at 37°C for 5 hours at 250 rpm

Added 10ul of kanamycin to LB medium and left the cell to grow overnight.

Seed stock: Select 1 colony from the M12 (MU/SPY) plate and set up overnight culture.

Set up overnight shake flask cultures for 3 GBS on pSB1C3 and control of **psB1C3 (spy with RFP insert)**

Added 50 ml LB medium with no antibiotic in each flask

Set to grow at 37 °C at 175 rpm.

Date: Monday, 05/September 2011 11:00

Topic: Miniprep for the 3 GBS shake flask and setup chloroquine gel

1) Prepared 3 mini-prep samples from each shake flask culture (50ul each)

2) Measured nano drop concentration.


3) Loaded a fixed mass (600ng) of DNA from each sample into a 1D chloroquine gel.

Date: Tuesday, 06/September 2011 11:00

Topic: visualised chloroquine gel for GBS

1) stained, destained and visualised chloroquine gel for 3 GBS

The 1D chloroquine agarose gel electrophoresis



Date: Wednesday, 07/September 2011 11:00


Topic: Shake flask Setup for MU/SPY and SPY and OD reading

- 1) Setup shake flask (37°C at 175 rpm for inoculation)
- 2) Collected 2 ml samples every 2 hours for 24 hours.
- 3) 1ml will be used to get OD reading and the other 1 ml sample will be centrifuged, pelleted and stored at -20 °C.

OD Results and Graph

Time (hrs)	OD (S)	OD (MS)
<i>Inoculum</i>	1.61	1.47
0	0.11	0.09
2	0.63	0.71
4	2.05	1.53
6	1.96	1.85
8	2.45	1.89
10	2.06	1.75

12	2.3	1.87
14	1.98	2.8
16	1.77	2.53
18	1.94	1.62
20	1.32	1.32
22	1.37	1.27
24	1.64	1.38



Date: Thursday, 08/September 2011 11:00

Topic: GFP Scan

- 1) Re-suspend the shake flask samples in 400 ul of water.
- 2) Loaded 200ul on to the 96-well plate for GFP Scan.
- 3) Samples were taken from the shake flask every 2 hours to detect the amount of fluorescent given by the cells.

GFP Scan parameters

SAFIRE II; Serial number: 12904200032; Firmware: V 2.10 12/2007		
Safire2; XFLUOR4SAFIREII Version: V 4.62n		
Date:		10/09/2011
Time:		14:27
Measurement mode:		Fluorescence Top
Excitation wavelength:		483 nm
Emission wavelength:		525 nm



Excitation bandwidth:		20 nm
Emission bandwidth:		20 nm
Gain (Manual):		62
Number of reads:		1
FlashMode:		High speed
Integration time:		20 μ s
Lag time:		0 μ s
Plate definition file:		GRE96ft.pdf
Plate with cover		
Part of the plate:		A1 - E12
Z-Position (Manual):		10200 μ m
Target Temperature:		22 $^{\circ}$ C
Current Temperature:		30.5 $^{\circ}$ C

Result – GFP SCAN

Time (hours)	2	4	6	8	10	12	14	16	18	20	22	24
Spy	613	162	215	258	252	280	247	325	328	364	295	3589
		1	0	1	1	9	0	0	2	2	4	
Mu/Spy	858	396	571	650	694	672	716	100	803	711	697	7426
		4	3	1	6	2	7	61	7	7	7	
Water	20	18	17	19	17	18	19	19	19	16	17	19
Spy	593	160	213	256	250	279	245	323	326	362	293	3570
		3	3	2	4	1	1	1	3	6	7	
Mu/Spy	838	394	569	648	692	670	714	100	801	710	696	7407
		6	6	2	9	4	8	42	8	1	0	

