

Protocol for preparation of competent cells for transformation

For two transformations

Materials:

- 0.1 M Calcium Chloride chilled on the ice;
- Overnight bacterial culture or bacterial colonies;

Procedure:

1. Add 20 μ l of the overnight bacterial culture or pick a colony to 1 ml of LB antibiotic liquid medium, Incubate at 37 degree in a shaker till the OD600 value reaches 0.4-0.6.
2. Put the tubes on ice to incubate for 5 min.
3. Pellet bacterial cells by 4 min centrifugation at 4000 rpm, discard the supernatant.
4. Resuspend cells in 600 μ l of ice-chilled 0.1 M Calcium Chloride solution. Incubate on ice for 30 min.
5. Centrifuge for 4 min at 4000 rpm in a microcentrifuge tube, discard the supernatant.
6. Resuspend the pelleted cells in 100 μ l of ice-chilled 0.1 M Calcium Chloride solution. Incubate on ice.
7. Add 50 μ l of the prepared cells to each tube containing DNA sample, mix and incubate on ice for 30 min.
8. Transform subsequently as the transformation protocol.

Note:

1. Make sure the cells are not left in the centrifuge at ambient temperature for more than 5 min as this will significantly decrease the transformation efficiency.
2. The rpm at centrifugation is not higher than 4000, as a high rpm may cause the lysis of cells.

3. Competent cells prepared with this protocol are suitable for direct use only.
Freezing down and storage at -70°C is not recommended.
4. The culture can be kept at 4 degree for one week and used for preparation of competent cells, but culture stored longer than 10 ten days is not suitable for competent cells.

References:

*Current protocols in molecular biology