

# NucleoSpin Plasmid QuickPure : Isolation of high-copy plasmid DNA from E. coli

This protocol is for a preparation of up to 15µg of high-copy plasmid or DNA using the Macherey-Nagel Nucleospin plasmid QuickPure Kit.

## Procedure

**1.** Use 1-5mL ( we used 1.5mL ) of a saturated E.Coli LB culture, pellet cells in a standard benchtop microcentrifuge for 30 s at 11,000 x g. Discard the supernatant and remove as much of the liquid as possible.

**2.** Add 250 µL Buffer A1. Resuspend the cell pellet completely by vortexing or pipetting up and down. Make sure no cell clumps remain before addition of Buffer A2!

**3.** Check Buffer A2 for precipitated SDS prior to use. If a white precipitate is visible, warm the buffer for several minutes at 30 – 40 °C until precipitate is dissolved completely. Cool buffer down to room temperature (18 – 25 °C).

**4.** Add 250 µL Buffer A2. Mix gently by inverting the tube 6 – 8 times. Do not vortex to avoid shearing of genomic DNA. Incubate at room temperature for up to 5 min or until lysate appears clear.

**5.** Add 300 µL Buffer A3. Mix thoroughly by inverting the tube 6 – 8 times. Do not vortex to avoid shearing of genomic DNA!

**6.** Centrifuge for 5 min at 11,000 x g at room temperature. Repeat this step in case the supernatant is not clear!

**7.** Place a NucleoSpin Plasmid Column in a Collection Tube (2 mL) and decant the supernatant from step 3 or pipette a maximum

of 750  $\mu\text{L}$  of the supernatant onto the column. Centrifuge for 1 min at 11,000 x g. Discard flow-through and place the NucleoSpin Plasmid Column back into the collection tube.

**8.** Add 600  $\mu\text{L}$  Buffer A4 (supplemented with ethanol ). Centrifuge for 1 min at 11,000 x g. Discard flow-through and place the NucleoSpin Plasmid Column back into the empty collection tube.

**9.** Centrifuge for 2 min at 11,000 x g and discard the collection tube.

**10.** Place the NucleoSpin Plasmid Column in a 1.5 mL tube and add 50  $\mu\text{L}$  Buffer AE. Incubate for 1 min at room temperature. Centrifuge for 1 min at 11,000 x g.

This protocol is extracted from "Plasmid DNA purification user manual" by **Macherey-Nagel**

