

Fermentas digestion

This protocol aims at cutting plasmids on specific restriction sites using restriction enzymes.

Procedure

- 1.** In an eppendorf tube, add the following solutions in any order :
 - 250 ng DNA (if the concentration is unknown, use 5 μ L)
 - 2.5 μ L (1X) or 5 μ L (2X) 10X tango buffer (this depends on the enzymes used)
 - water to get a volume of 25 μ L after addition of the enzymes
- 2.** Add 1 μ L of each one of the desired enzymes.
- 3.** Incubate at 37°C for 1h30
- 4.** Incubate at 70°C for 10 mn to inactivate the restriction enzymes.

