

Fluorescence Microscopy – Smear Culture Technique (Fixed Bacteria)

- 4mL of LB broth with Ampicillin (the antibiotic resistance of the pSB1A3 backbone, into which the constructs were cloned) was inoculated from colonies off a master plate (or a scrape off a frozen glycerol stock) of positive clones and grown overnight at 37°C in a shaking incubator.
- 1mL of each bacterial culture was then used to inoculate 25mL of LB broth with Ampicillin and grown for \pm 2 hours in a shaking incubator at 37°C until the cells reached mid-log phase.
- 2 x 5mL of each sample culture was aliquoted into separate 50mL falcon tubes (to have enough aeration during induction).
- One tube was left untreated; 150 μ L of 15mM theophylline was added to the other (to make a concentration of 1.5mM theophylline – to induce activation of the theophylline riboswitches).
- These 5mL aliquots were placed in a shaking incubator at 37°C for 90 min each to allow for activation.
- 1mL of each culture was then aliquoted into a 1.5mL Eppendorf tube and centrifuged for 10min at 4000rpm.
- Supernatant was poured off.
- The bacterial pellet was re-suspended in 50 μ L of LB broth.
- 5 μ L of these re-suspended cultures were then mixed with 5 μ L of an aqueous mounting medium on normal glass microscope slides.
- This mixture was smeared across the glass slide with a long cover slip, before placing the cover slip on top of this smear.
- Slide is ready for viewing under a fluorescent microscope using a 100X objective.