

This protocol is based on the Qiagen Miniprep instruction manual.
This document is version 1.00: last updated 8.12.11.

∅∅∅∅∅ **DNA Purification (PCR/Digest)**

This protocol will purify DNA based on the QIAGEN PCR purification system. DNA must be dissolved in a binding buffer, then run through the column. This protocol can operate for the purification of both digests and PCR, though PCR purification is its primary function.

Solutions

PCR/Digest solution
Buffer PB
Buffer PE
Buffer EB

Materials

100,1000µl Pipettes
Centrifuge tubes
Spin columns with flow-through containers

*You will also need access to a
Centrifuge*

Procedure

1. Add **5x(volume of DNA solution)µL Buffer PB** to each centrifuge tube. Mix by vortexing.
2. Add the contents of the tube to a **spin column** with a **1000µL pipette**.
3. Place the column in a **flow-through container**. Centrifuge at **13 000 rpm** for **1 minute**. Decant the flow-through.
4. Add **750µL Buffer PE** with a **1000µL pipette**. Centrifuge at **13 000 rpm** for **1 minute**. Decant the flow-through.
5. Centrifuge at **13 000 rpm** for **1 minute**. This will remove any residual buffer in the column. It is not necessary to decant the flow-through, as the column collection tube will be discarded.
6. Remove the column from its collector and insert into a fresh, labeled **1.5mL centrifuge tube**.
7. Add **30-50µL Buffer EB** to the center of the column's filter with a **100µL pipette**. Ensure the buffer is absorbed by the filter.
8. Centrifuge at **13 000 rpm** for **5 minutes**. The flowthrough will contain the desired DNA.