

Adapted from a UCLA protocol.  
This document is version 1.02. Last updated: 6.12.11.

## ⊘⊘⊘⊘ **Making LB Plates**

This protocol will generate plates containing LB nutrient broth to grow bacteria, agar, and your choice of antibiotic. While generating the initial LB is does not require strong sterile technique as it will be sterilized in the autoclave, a strongly bacteria-free environment will be necessary for pouring the plates.

After agar LB is autoclaved, the liquid should be handled very carefully-any introduced bacterial will contaminate the stock.

This protocol will produce 25 plates. The external protocol on *autoclaving* will be used in addition to this one. A 1 hour gap separates procedures *I* and *II*, where *I* prepares the agar solution and *II* fills the plates.

### *Compounds*

dH<sub>2</sub>O

LB agar powder

Your favorite antibiotic

### *Materials*

1000mL graduated cylinder

500mL Erlenmeyer flask

10mL Autopipettes

Pipettes capable of 250μL, tips

Weighing boats

Petri dishes

Permanent markers

Masking tape

*You will also need access to an*

Autoclave

Refrigerator

Scale

### *External protocols*

Autoclaving

### *Protocol I*

1. Mass **10g LB agar powder** in a **weighing boat**.
2. Pinching the sides of the boat, pour the entire contents of the boat into an empty **500mL Erlenmeyer flask**.
3. Fill a **1000mL graduated cylinder** with **250mL dH<sub>2</sub>O**.
4. Empty the contents of the graduated cylinder into the Erlenmeyer flask.

⇒ Autoclave in accordance with the **autoclaving protocol**.

### *Protocol II*

1. Let the autoclaved **LB agar solution** sit on a counter until cool enough to hold without a glove.
2. **Ensuring the work area is sterile**, remove the sterile **Petri dishes** from the bag and stack in 5-high piles.
3. Carefully remove the cover on the **LB agar solution**.

4. Add **250 $\mu$ L** of **your favorite antibiotic** to the LB solution.
5. Transfer **10mL LB agar solution** to one of the plates with a **10mL autopipette**. Carefully lift the lid of the plate just high enough to fit the pipette under. Swirl each plate to distribute the agar. Once the pipette has been emptied, set the plate off to the side, spreading them out over the counter.
6. Repeat **step 5** with the other **24 plates**. Wait **20 minutes** or until the agar is solid.
7. Flip each plate over so the agar side of the plate is oriented upwards to prevent condensation. Stack and store in a sterile plastic bag, making sure to label the bag.

Despite being treated with antibiotic, each plate can still be contaminated and should be handled carefully and in a sterile manner.