

PCR protocol

1. Pipet the following into a microfuge tube:

reagent	50µ l reaction system	final concentration
2× EsTaq MasterMix	25µ l	1×
Forward Primer, 10µ M	2µ l	0.4µ M
Reverse Primer, 10µ M	2µ l	0.4µ M
Template DNA	<1µ g	<1µ g/reaction
RNase-Free Water	up to 50µ l	

The table below provides standard reaction conditions for PCR.

Mg ²⁺	KCl	dNTPs	Primers	DNA polymerase	Template
DNA					
1.5 mM	50 mM	200 µM	1 µM	1-5 units	1 pg to 1 µg

2. cycling condition:

procedure	temperature	time
initial denaturation	94°C	2min
denaturation	94°C	30s
primer annealling	55-65°C	30s
extending	72°C	1min/Kb
final extending	72°C	2min
cooling down	4°C	10min