



Clean deletions

pMAD see Arnaud et al., Appl Environ Microbiol 2004, 70:6887
(<http://www.ncbi.nlm.nih.gov/pubmed/15528558>)

1. Cloning
 - a. Amplify up and down fragments (~1 kb upstream and downstream from the fragment to be deleted)
 - b. Join both fragments (see joining PCR in PCR methods)
 - c. Clone joined fragments in pMAD (in *E. coli* DH5 α)
2. Transform the sequenced plasmid into *E. coli* BMH7118 and isolate it again
3. Transformation in *B. subtilis*
 - a. Transform the plasmid isolated from *E. coli* BMH7118 (~1,5 – 2 μ g) in *B. subtilis* W168
 - b. Select the transformants at 30°C on plates with X-gal (100 μ g/ml) and MLS selection
 - c. Pick blue colonies
4. 1st temperature shift
 - a. Prepare an overnight culture with MLS selection at 30°C
 - b. Inoculate 10 ml LB with MLS at OD₆₀₀ 0,1 and incubate at 30°C for 2 hours
 - c. Shift the temperature to 42°C and incubate for 6 hours
 - d. Plate out dilutions (10⁻² to 10⁻⁵) on plates with X-gal and MLS and incubate at 42°C
5. Colony PCR - check the integration of pMAD in the genome



Possibilities after recombination:



1&2 pMAD check fwd and rev primers (253 & 254)

6. 2nd temperature shift
 - a. Inoculate the positive cones in LB without selection and incubate for 6 hours at 30°C
 - b. Shift the temperature to 42°C and incubate for 3 hours
 - c. Plate out dilutions (10^{-2} to 10^{-5}) on plates with X-gal (without selection) and incubate at 42°C
 - d. Pick white colonies and check the deletion with colony PCR
 - e. Check if the mutants are MLS sensitive!