

This protocol is based on a protocol by Knight.
This document is version 1.00. Last updated: 8.12.11.

∅∅∅∅∅ **Phusion PCR**

PCR is a rapid method of amplifying low concentrations of DNA into high concentrations capable of being viewed or used.

This protocol will utilize NEB's *Phusion* polymerase, a high power and high accuracy polymerase capable of amplifying large DNA fragments at very high rates.

******Phusion* should only be used when for cloning or producing stocks when the DNA of interest is known to work with PCR! Do not use *Phusion* for colony PCR or routine assays unless absolutely necessary, as it is five times as expensive as standard *taq* polymerase!**

Solutions

HF or GC buffer (depending on DNA %GC content)

dNTP solution

Forward primer solution, concentrated

Reverse primer solution, concentrated

Phusion polymerase solution

ddH₂O

Template solution

Materials

10, 50 μ L pipette

PCR tubes

1.5mL centrifuge tubes.

You will also need access to a

PCR machine

Tabletop quick centrifuge

Procedure

1. Generate a PCR master mix by adding these reagents to a **1.5mL centrifuge tube**.
 \Rightarrow *This reaction will produce a 100 μ L mix; 50 μ L reactions (or less) can be easily achieved by dividing these volumes proportionally.*

62 μ L **diH₂O** or **ddH₂O**

2.0 μ L **dNTP solution**

5.0 μ L **Forward primer solution**

5.0 μ L **Reverse primer solution**

1.0 μ L ***Phusion* polymerase solution**

2. Properly divide the solution among your samples and add it to each tube in a **clean PCR tube strip**.
3. Add 5/(# of samples) μ L **template solution** to each tube. Mix by uptaking and expelling liquid with your pipette.
4. Centrifuge the tubes in a **tabletop quick centrifuge** for **20 seconds**.
5. Design a program for your PCR reaction. If using VF2 and VR primers, run the samples on the following PCR program:

A **95°C for 15 mins**

B **94°C for 15 seconds**

C **56°C for 20 seconds**

D **68°C for 30 seconds** *per kb expected fragment size*

Repeat B-D 35 or so times.

E **68°C for 20 minutes**

F **4°C indefinitely**

6. View the results of the reaction using the **gel electrophoresis** protocol.