

This protocol is based on a protocol by Knight.  
This document is version 1.04. Last updated: 8.25.11.

## ∅∅∅∅∅ **Colony PCR**

Colony PCR is a fast and efficient method of screening ligated colonies for properly ligated inserts. As many ligations rely on the random chance of a part's proper ligation, it becomes important to screen colonies bearing ligated plasmids for these proper clones. Colony PCR suspends raw bacterial colonies in water or LB, then adds the solution itself to the PCR mixture; the high temperature of the PCR reaction will lyse the cells and free the DNA.

While it is possible to screen by growing overnight, digesting, and running a subsequent gel, this method takes far more time and resource; by the same token, however, colony PCR can be highly unreliable due the wildly uncontrollable nature of what is essentially contaminating a PCR mixture with the complete contents of a bacterial cell.

Before starting, decide how many colonies you intend to pick and screen.

### *Solutions*

5Prime master mix  
diH<sub>2</sub>O or ddH<sub>2</sub>O  
Bacterial colonies (on plate)  
Forward, reverse primer solution

### *Materials*

10, 50 $\mu$ L pipette  
PCR tubes  
Colony picker

*You will also need access to a  
PCR machine*

### *Procedure*

1. Add **20 $\mu$ L H<sub>2</sub>O** to tubes in a PCR tube strip. One tube should exist for every colony screen desired.
2. Using a **sterile colony picker**, pick a colony and rub the picker around inside the PCR tube until a cloudy-ish solution forms. This should be done for all colonies you intend to screen in the 20 $\mu$ L. Sterilize the implement in between picking. This will be your **colony template solution**.
3. Generate a PCR master mix.  
*⇒ Multiply this mixture by the number of samples you intend to screen and run a subsequent PCR reaction on. This is necessary, as it can be very difficult to accurately measure 0.2 $\mu$ L of polymerase in glycerol. It also can be helpful to add 3-5 $\mu$ L extra water to the master mix (or make 10-20 $\mu$ L extra), in case you run out of mix for your last sample due to pipetting error.*

**4 $\mu$ L diH<sub>2</sub>O or ddH<sub>2</sub>O**  
**0.25 $\mu$ L Forward primer solution**  
**0.25 $\mu$ L Reverse primer solution**  
**5 $\mu$ L 5Prime master mix**

4. Add **9.5 $\mu$ L PCR master mix** to a **second clean PCR tube strip** according to your sample.
5. Add **0.5 $\mu$ L colony template solution** to each tube containing PCR mix.

6. Design a program for your PCR reaction. If using VF2 and VR primers, run the samples on the following PCR program:  
⇒ *Alternatively, use the PCR or similar program already entered into the PCR machine. The current standard program is "MSBTVEVR."*

A **95°C for 15 mins**

B **94°C for 30 seconds**

C **56°C for 30 seconds**

D **68°C for 1 minute** *per kb expected fragment size*

Repeat B-D 35 or so times.

E **68°C for 20 minutes**

F **4°C indefinitely**

7. View the results of the reaction using the **gel electrophoresis** protocol.