

## Lab plan for construct of Pveg+RBS+YbdN+BMP-2+BBa\_0015

### Final construct:

**Pveg promoter + RBS + Signal peptide YbdN + Mature mouse BMP-2+ BBa+0015  
double terminator in pSB1AK3**

### Lab plan

#### *Construction of Pveg +RBS+YbdN+BMP-2 in pSB1C3*

##### **1. Double Digestion of YbdN+BMP-2 in pSB1C3**

Double digest construct YbdN+BMP-2 in pSB1C3 with XbaI and PstI

##### **2. Gel Electrophoresis of digestion product**

- Prepare 2% agarose gel
- Expected band size: 2070 b.p. and 439 b.p.

##### **3. Gel Purification of YbdN+BMP-2**

- Cut the gel that shows 439 b.p. band.
- Purify gel using Favorgen gel purification kit.

##### **4. Double Digestion of Pveg+RBS in pSB1C3**

Double digest construct Pveg+RBS in pSB1C3 with SpeI and PstI

##### **5. Dephosphorylation of digested Pveg+RBS in pSB1C3**

Into the digestion product of Pveg+RBS in pSB1C3, add 2  $\mu$ l of Antarctic phosphatase buffer and 1  $\mu$ l Antarctic phosphatase.

##### **6. DNA purification of digested Pveg+RBS in pSB1C3**

DNA purification using Favorgen PCR cleanup kit.

##### **7. Ligation of Pveg+RBS in pSB1C3 and YbdN+BMP-2**

- Mix the insert: YbdN+BMP-2 and backbone: Pveg+RBS in pSB1C3 with ratio of backbone:insert equals to 1:4.
- Into the mixture also add T4 ligase buffer and T4 ligase
- In another tube, prepare the mixture using backbone only without insert as negative control.

## **8. Transformation of Pveg+RBS+YbdN+BMP-2 in pSB1C3 into *E. coli***

- a. Take out the competent cell from -80 freezer.
- b. Cool down the cell in ice for 10 minutes.
- c. Add all ligation products into the cells. Leave it in ice for 10 minutes.
- d. Put the tubes into 42°C water bath for exactly 90 seconds.
- e. Put them back into ice for 2 minutes.
- f. Add 1ml LB in all cells. Incubate the tubes in 37°C for 1 hour.
- g. Spin down the cells. Transfer 100  $\mu$ l of each tube of cells into separate Chloramphenicol (25ng/  $\mu$ l) plates. Spread the plate. Incubate in 37°C overnight.
- h. Observe the colonies formed on next day.

## **9. Inoculation of Pveg+RBS+YbdN+BMP-2 in pSB1C3 in *E. coli***

- a. Prepare Chloramphenicol (25ng/  $\mu$ l) LB solution.
- b. Add 5 ml of LB solution into sterile falcon tubes.
- c. Pick colony from transformed plate and dip into LB solution in falcon tube.
- d. Incubate at 37°C overnight

## **Construction of Pveg +RBS+YbdN+BMP-2+BBa\_0015 in pSB1AK3**

### **1. Plasmid extraction of Pveg+RBS+YbdN+BMP-2 in pSB1C3**

Extract the plasmid using Favorgen plasmid extraction miniprep kit.

### **2. Double Digestion of Pveg+RBS+YbdN+BMP-2 in pSB1C3**

Double Digest construct Pveg+RBS+YbdN+BMP-2 in pSB1C3 with EcoRI and SpeI

### **3. Gel electrophoresis of digested Pveg+RBS+YbdN+BMP-2 in pSB1C3**

- a. Prepare 2% agarose gel
- b. Expected band size: 2070 b.p. and 570 b.p

### **4. Gel Purification of digested Pveg+RBS+YbdN+BMP-2**

- a. Cut the gel that shows 570 b.p. bands
- b. Purify the DNA from the gel using Favorgen gel purification kit.

### **5. Double Digestion of BBa\_0015 in PSB1AK3**

Double digest biobrick BBA\_0015 in PSB1AK3 with EcoRI and XbaI

## 6. Dephosphorylation of Bba\_0015 in pSB1AK3

Into the digestion product of Bba\_0015 in pSB1AK3, add 2  $\mu$ l of Antarctic phosphatase buffer and 1  $\mu$ l Antarctic phosphatase.

## 7. DNA purification of digested BBa\_0015 in PSB1AK3

DNA purification using Favorgen PCR cleanup kit.

## 8. Ligation of Pveg+RBS YbdN+BMP-2 with BBa\_0015 in PSB1AK3

- Mix the insert: Pveg+RBS+YbdN+BMP-2 and backbone: BBa\_0015 in pSB1AK3 with ratio of backbone:insert equals to 1:4.
- Into the mixture also add T4 ligase buffer and T4 ligase
- In another tube, prepare the mixture using backbone only without insert as negative control.

## 9. Transformation of Pveg+RBS+YbdN+BMP-2+BBa\_0015 in pSB1AK3 into E.coli

- Take out the competent cell from -80 freezer.
- Cool down the cell in ice for 10 minutes.
- Add all ligation products into the cells. Leave it in ice for 10 minutes.
- Put the tubes into 42°C water bath for exactly 90 seconds.
- Put them back into ice for 2 minutes.
- Spread the onto Ampicilin (150 ng/ $\mu$ l) plate. Incubate in 37°C overnight.
- Observe the colonies formed on next day.

## Construct Pveg+RBS+YbdN+BMP-2+BBa\_0015 screening

### 1. Colony PCR of Colonies from plate Pveg+RBS+YbdN+BMP-2+BBa\_0015 in pSB1AK3 in E.coli

Prepare the following mixture:

Sample	Each PCR Tube
ddH <sub>2</sub> O	16
10X thermopol PCR buffer	2
Forward Primer (10 $\mu$ M)	0.5
Reverse Primer (10 $\mu$ M)	0.5
dNTP (10mM)	0.5
Taq polymerase	0.5
Total Volume ( $\mu$ l)	20

Forward Primer sequences:

5'- GATCATTCTAGAGAAAGGAGGTTGTTTGCATGGTG-3'

Reverse Primer sequences:

5'-ATGATCACTAGTATTATTAAACGACACCCGCAGCCC-3'

Thermocycle temperature setup

Cycle number	Denaturation	Annealing	Polymerization
Initial denaturation 1 cycle	10 min at 95 °C	-	-
30 cycles	30 sec at 95 °C	30 sec at 64°C	1 min at 68 °C
Final extension 1 cycle	-	-	5 min at 68 °C

## 2. Gel electrophoresis of colony PCR product

- Prepare 2% agarose gel
- Expected band size: 439 b.p.