

β-Galactosidase Assay for *E. coli* (Miller, 1972)

Example of culture preparation

- Inoculate 10 ml LB medium 1:100 or to an OD₆₀₀ 0.1 with a fresh overnight culture carrying a promoter-*lacZ*-fusion and incubate on a shaker at 37°C
- At OD₆₀₀ 0.4-0.8 split the culture to 3 ml samples, induce one sample with e. g. an antibiotic, leave one sample as an uninduced control (or whatever conditions you need)
- After induction (for example, 30 min), measure OD₆₀₀ and harvest 2 ml of cell culture by centrifugation
- Store cell pellets at -20°C or continue directly with the assay

β-Galactosidase Assay

- Resuspend the cell pellet in 2 ml Z-buffer
- In a 2 ml eppendorf cup, prepare three different dilutions with Z-buffer (for example, 0, 1:5 and 1:10, final volume 1 ml), use 1 ml Z-buffer as reference
- Add 25 µl SDS and 50 µl chloroform (from here on, work under fume hood) and mix by vortexing
- Incubate 5 min at room temperature
- Add 200 µl ONPG, mix well and record time (=t₀)
- Incubate at room temperature until the sample turns yellow
- Stop the reaction by adding 500 µl Na₂CO₃, mix well and record time (=t_s)
- If the samples do not turn yellow, stop the reaction after 60 min
- Centrifuge (7 min, 13000 rpm, RT)
- Measure OD₄₂₀ of the supernatant, use a cuvette with everything but the cells as blank
- Calculate promoter activity according to the formula:

$$MillerUnits = \frac{OD_{420} * V * 1000}{t * OD_{600}}$$

t time of reaction in min (T_s - T₀)
 V dilution factor

Solutions

Na ₂ CO ₃	1 M
ONPG (2-nitrophenyl-β-D-galactopyranoside)	4 mg/ml in Z-buffer
SDS	0.1% (w/v)
Chloroform	

Z-buffer (pH 7.0)	Na ₂ HPO ₄ * 2 H ₂ O	60 mM	10.68 g
	NaH ₂ PO ₄ * H ₂ O	40 mM	5.52 g
	KCl	10 mM	0.75 g
	MgSO ₄	1 mM	0.24 g
	H ₂ O		ad 1000 ml

In the original protocol, the Z-buffer contains 100 µg/ml Chloramphenicol. I did the assay without Chloramphenicol and it worked.

Protocol generously provided by the lab
Prof. Thorsten Mascher
Großhadernerstr. 2-4
82152 Planegg-Martinsried
www.syntheticmicrobe.bio.lmu.de