

## Construction of Pveg promoter + spoVG RBS + lytC + linker + RPMrel + GFP + double terminator (K733007)

### STAGE 1

#### 1) PCR amplification of [consensus RBS + GFP + double terminator] region using BBa\_E0840 as the template.

##### *Forward primer design*

5' – [17bp prefix region, including XbaI restriction site] [11bp B. subtilis consensus RBS] [22bp overlap with GFP including 3bp spacer preceding GFP start codon] – 3'

##### *Forward primer sequence*

5' – CGCGGCCGCTTCTAGAGAAAGGAGGTGTTAGATGCGTAAAGGAGAAGAAC – 3' (50bp)

##### *Reverse primer design*

5' – [18bp overlap with pSB1A2 plasmid downstream of suffix] – 3'

##### *Reverse primer sequence*

5' – TACCGCCTTTGAGTGAGC – 3' (18bp)

##### *PCR reaction set-up*

Reagents	Quantity per reaction
ddH <sub>2</sub> O	Add up to total 20μL volume
Phusion® HF Buffer (5X)	4μL
DNA Template - BBa_E0840	10ng
Forward Primer (10μM)	0.5μL
Reverse Primer (10μM)	0.5μL
dNTP (10μM)	0.5μL
Phusion® Polymerase	0.5μL

##### *PCR program*

Step Number	Step Name	Temperature (°C)	Duration (minutes:seconds)
1	Initial Denaturation	95	05:00
2	Denaturation	95	00:30
3	Annealing	56	00:30
4	Extension	72	01:00
5	Goto step 2 for 30 cycles	-	-
6	Final Extension	72	05:00
7	Storage	8	∞

#### 2) Reaction clean-up of above PCR.

Purification of the PCR products using the materials and protocols packaged in the FavorPrep™ GEL/PCR Purification Kit.

### 3) Digestion of [consensus RBS + GFP + double terminator] PCR product and pSB1C3.

*Digestion reaction set-up*

Reagents	Quantity per reaction	
	PCR Product	pSB1C3
ddH <sub>2</sub> O	Add up to total 20μL volume	Add up to total 20μL volume
10X NEBuffer 4	2μL	2μL
DNA	1μg	1μg
XbaI	0.5μL	0.5μL
PstI-HF	0.5μL	0.5μL

Reaction incubation at 37°C for 2 hours.

### 4) Gel purification of digested pSB1C3.

Run the digestion product in 1% agarose gel for an appropriate time. Identify the band resolved at 2070bp and cut it from the gel. This is the pSB1C3 backbone. Purify the gel cutting using the materials and protocols packaged in the FavorPrep™ GEL/PCR Purification Kit.

### 5) Reaction clean-up of PCR product digestion reaction.

Purify the digestion products using the materials and protocols packaged in the FavorPrep™ GEL/PCR Purification Kit.

### 6) Ligation of PCR product into standard backbone pSB1C3.

*Ligation reaction set-up*

Reagents	Quantity per reaction	
	Ligation	(-) Control
ddH <sub>2</sub> O	Add up to total 10μL volume	Volume of insert, then add up to total 10μL volume
10X T4 Ligase Buffer	1μL	1μL
Insert (mass ratio value)	3	-
Backbone (mass ratio value)	1	1
T4 Ligase	0.5μL	0.5μL

Reaction incubation at room temperature for 2 hours.

### 7) Transformation of above ligation product.

We transformed using competent *Escherichia coli* DH10B cells. We are also selecting for chloramphenicol resistance. The following steps take these factors into account.

- Add whole ligation product to 100μL competent cells.
- Place mixture on ice for 10 minutes.
- Heat shock mixture at 42°C for 90 seconds.
- Place mixture on ice for 2 minutes.
- Add 1ml LB to the mixture.

- Recover the cells by incubation at 37°C for 1 hour.
- Centrifuge mixture at 16.1rcf for 90 seconds.
- Remove majority of the supernatant.
- Resuspend pellet in an appropriate volume of LB.
- Transfer suspension to 25µg/ml chloramphenicol LB plate and spread.
- Incubate at 37°C overnight.

### 8) Plasmid extraction of successfully transformed colonies.

Extract [consensus RBS + GFP + double terminator] construct in pSB1C3 using the materials and protocols packaged in the FavorPrep™ Plasmid DNA Extraction Mini Kit.

## STAGE 2

### 1) PCR amplification of [Pveg promoter + spoVG RBS + lytC + linker + RPMrel] region using BBa\_K316037 as the template.

*Forward primer design*

5' – [6bp cap] [20bp overlap with standard prefix] – 3'

*Forward primer sequence*

5' – GATCATGAATTCGCGGCCGCTTCTAG – 3' (26bp)

*Reverse primer design*

5' - [8bp cap] [7bp SpeI restriction site] [6bp reverse-complementary double stop codon] [27bp reverse-complementary sequence of codon-optimized RPMrel] [15bp reverse-complementary overlap with linker] - 3'

*Reverse primer sequence*

5' - GTTCTTCACTAGTATTATTAACACATCGGGCGATCTTCGATCGGACAGGCCGCGGCTTTCGC - 3' (63bp)

*PCR reaction set-up*

Reagents	Quantity per reaction
ddH <sub>2</sub> O	Add up to total 20µL volume
Phusion® HF Buffer (5X)	4µL
DNA Template - BBa_K316037	100ng
Forward Primer (10µM)	0.5µL
Reverse Primer (10µM)	0.5µL
dNTP (10µM)	0.5µL
Phusion® Polymerase	0.5µL

*PCR program*

Step Number	Step Name	Temperature (°C)	Duration (minutes:seconds)
1	Initial Denaturation	95	05:00
2	Denaturation	95	00:30

3	Annealing	65	00:30
4	Extension	72	01:30
5	Goto step 2 for 25 cycles	-	-
6	Final Extension	72	05:00
7	Storage	8	∞

## 2) Reaction clean-up of above PCR.

Purify the PCR products using the materials and protocols packaged in the FavorPrep™ GEL/PCR Purification Kit.

## 3) Digestion of [Pveg + spoVG RBS + lytC + linker + RPMrel] PCR product and [consensus RBS + GFP + double terminator] construct in pSB1C3.

*Digestion reaction set-up*

Reagents	Quantity per reaction	
	PCR Product	pSB1C3 Construct
ddH <sub>2</sub> O	Add up to total 20μL volume	Add up to total 20μL volume
10X NEBuffer 2	2μL	2μL
DNA	1μg	1μg
EcoRI	0.5μL	0.5μL
XbaI	-	0.5μL
SpeI	0.5μL	-

Reaction incubation at 37°C for 2 hours.

## 4) Reaction clean-up of PCR product and pSB1C3 construct digestion reaction.

Purify the digestion products using the materials and protocols packaged in the FavorPrep™ GEL/PCR Purification Kit.

## 6) Ligation of PCR product into standard backbone pSB1C3.

*Ligation reaction set-up*

Reagents	Quantity per reaction	
	Ligation	(-) Control
ddH <sub>2</sub> O	Add up to total 10μL volume	Volume of insert, then add up to total 10μL volume
10X T4 Ligase Buffer	1μL	1μL
Insert (mass ratio value)	3	-
Backbone (mass ratio value)	1	1
T4 Ligase	0.5μL	0.5μL

Reaction incubation at room temperature for 2 hours.

## 7) Transformation of above ligation product.

Again consider the following steps suitable for transformation of *E. coli* DH10B competent cells for selection with chloramphenicol.

- Add whole ligation product to 100µL competent cells.
- Place mixture on ice for 10 minutes.
- Heat shock mixture at 42°C for 90 seconds.
- Place mixture on ice for 2 minutes.
- Add 1ml LB to the mixture.
- Recover the cells by incubation at 37°C for 1 hour.
- Centrifuge mixture at 16.1rcf for 90 seconds.
- Remove majority of the supernatant.
- Resuspend pellet in an appropriate volume of LB.
- Transfer suspension to 25µg/ml chloramphenicol LB plate and spread.
- Incubate at 37°C overnight.

### **8) Plasmid extraction of successfully transformed colonies.**

Extract [Pveg promoter + spoVG RBS + lytC + linker + RPMrel + consensus RBS + GFP + double terminator] construct in pSB1C3 using the materials and protocols packaged in the FavorPrep™ Plasmid DNA Extraction Mini Kit.