

## **Overnight culture with one antibiotic**

*DO NOT FORGET to perform a glycerol stock from the performed overnight culture! See bottom of the page for instructions*

### **Procedure:**

The amount of LB should be "6X+2" ml, where X is the number of the samples with the given antibiotic resistance, and 1 ml for the negative control

- Label "6X+2" number of 50 ml Falcon tubes with the proper sample name, clone number and negative control
- Put "6X+2" ml of LB into the tube of negative control
- Add antibiotic into the LB of the negative control
  - Since we have the stock with the right, add "6X+2"  $\mu$ l of antibiotic will be good
  - Mix the antibiotic with LB by inverting the Falcon tube several times
- Take up 1 ml of LB with antibiotic into each sample tube
- Pick one colony of each sample clone with pipette tip, and transfer into the LB solution
- Make sure everything is sterile and CHANGE pipette tip between sample. You can leave the pipette tip in the O/N culture.
- Grow the samples and the negative control overnight in 37°C in a shaker

Expected result: no growth in negative control, full growth in all the sample tubes

Repeat the process for every other antibiotic needed.

## **Glycerol stock**

Take the overnight culture from the previous description, transfer 600  $\mu$ l of culture sample into 2 ml tubes. Add 400  $\mu$ l of 50% glycerol to obtain a final volume of 1 ml and 20% glycerol concentration.

Store the labeled 2 ml tubes in -80°C freezer indefinitely.