

iGEM2013 – Microbiology – BMB – SDU	
Project type: USER cloning	Creation date: 6/8-13
Project title: USER cloning of pSB1C3 and AmilCP	Written by: MHK and SF
Sub project:	Performed by: MHK and SF

1. SOPs in use

SOP0006_v01_USER_cloning

2. Purpose

To compare the colonies of E. coli containing just AmilCP with the colonies containing AmilCP and DXS in the pSB1C3 plasmid.

3. Overview

Day	SOPs	Persons	Experiments
1	SOP0006_v01	MHK, SF SF	USER PCR of pSB1C3 and AmilCP USER PCR of Amil
2	SOP0014_V01	SF, MHK	Gelpurification of ON PRC of Amil and Plasmid
3	SOP0007_v01 SOP0009_v01	SF SF	Excision and ligation of PCR products in USER cloning TSB Trnsformation

4. Materials required

Materials in use

Name	Components (Concentrations)	Manufacturer / Cat. #	Room	Safety considerations
Primer 7				
Primer 13				
Primer 52				
Primer 53				
Grøn 24				
Grøn 26				

5. Other comments

6. Experiment history

Date (YY.MM.DD)	SOPs	Alterations to SOPs and remarks to experiments
13.08.06	SOP0006_v01	<p>The PCR program was altered from the SOP version to the program used for succesful plasmid PCR:</p> <p>95 deg 2 min 95 deg 30 sec* 48 deg 30 sec* 68 deg 2:30 min* 72 deg 10 min *37 cycles</p> <p>Amil Template: grøn 24 Primers: 13+52</p> <p>Plasmid Template: Grøn 26 Primers: 7+53</p> <p>50 µL was loaded in each well on the gel.</p>

13.08.06	SOP0006_v01	<p>The PCR program:</p> <p>95 deg 2 min 95 deg 30 sec* 38 deg 30 sec* 72 deg 1:00 min* 72 deg 10 min *30 cycles</p> <p>amil Template: grøn 24 Primers: 13+52</p> <p>ON</p>
13.08.06	SOP0006_v01	<p>The PCR program was altered from the SOP version to the program used for succesful plasmid PCR:</p> <p>95 deg 2 min 95 deg 30 sec* 48 deg 30 sec* 68 deg 2:30 min* 72 deg 10 min *37 cycles</p> <p>Plasmid band 1-3 Template: Grøn 26 Primers: 7+53</p> <p>Plasmid band 4-6 Template: Green 106 Primers: 7+53</p> <p>ON</p>

13.08.08	SOP0007_v01	<p>The reaction mixtures were:</p> <p>10 fmol plasmid : 20 fmol amil 2 µL buffer 1 µL USER enzyme water to 20 µL</p> <p>and:</p> <p>10 fmol plasmid : 50 fmol amil The rest as above</p> <p>Control:</p> <p>10 fmol plasmid : 0 fmol amil The rest as above.</p> <p>Incubated for 30 min at 37 deg. C. Incubated for 30 min at 25 deg. C.</p>
13.08.08	SOP0009_v01	<p>0,55 µL culture was used for ech sample. 10 µL of each sample from above was used. 1 µL RFP (on plasmid) was used. One sample was pure bacteria.</p> <p>55 min Phenotypical expression.</p>

7. Sample specification


Sample name	Sample content	From	Used for / Saved where
Green 106	pSB1C3 USER	USER PCR from first PCR	iGEM fridge
Green 107	pSB1C3 USER	USER PCR from second PCR	iGEM fridge
Green 108	pSB1C3 USER	USER PCR from second PCR	iGEM fridge
Green 109	AmilCP USER	USER PCR from second PCR	iGEM fridge
Green 110	AmilCP USER	USER PCR from second PCR	iGEM fridge

8. Remarks on setup

9. Results and conclusions

06.08.13

PCR was run on 1% agarose gel.




Well 2-4= amil, well 5-7=plasmid. Upper band in well 5 was cut out and purified.

Nanodrop:

Green 106: 3,3 ng/ μ L

13.08.07

PCR Amil:




the two bands were cut out and purified.

Nanodrop:

Grøn 109: 26,4

Grøn 110: 82,4.

PCR Plasmid:



The slight bands in well 2-4 was purified together. The band in well 5 was also purified.

Nanodrop:

Grøn 107: 9,3 ng/ μ L

Grøn 108: 5,0 ng/ μ L

10. Appendices