

<b>iGEM2013 – Microbiology – BMB – SDU</b>	
<b>Title:</b>  <b>SOP number:</b> SOP0015_v01  <b>Version number:</b> 01	<b>Date issued:</b> 2013.06.19  <b>Review date:</b>  <b>Written by:</b> ASF

### 1. Purpose

To ligate pieces of DNA

### 2. Area of application

Cloning

### 3. Apparatus and equipment

Apparatus/equipment	Location (Room number)	Check points	Cr
Vortex		•	
Pipettes (p20, p10)		•	
		•	
		•	
		•	
		•	
		•	

#### 4. Materials and reagents – their shelf life and risk labelling

Name	Components	Supplier / Cat. #	Room (hallway storage)	Sa
Purple pipette tips		Contact lab-manager	Micro storage	
Green pipette tips		Contact lab-manager	Micro storage	
Eppendorftubes		Contact lab-manager	Micro storage	
Distilled water		Contact lab-manager	Micro storage	
Ligasebuffer		Agilent Technologies	Freezer at 1. Floor	
Ligase				
DNA piece 1				
DNA piece 2				

#### 5. QC – Quality Control

#### 6. List of other SOPs relevant to this SOP

#### 7. Environmental conditions required

#### 8. Procedure

1. Prepare the ligation mixture and mix by pipetting up and down
2. Leave the mixture overnight at 16°C
- 2a. If there is no time leave the ligation solution at 22.5°C for 30mins. Then denature the ligase at 65°C for 10min.
3. Use ligation solution for transformations

Reagents	Volume
10x T4 DNA ligase buffer	2 µL
T4 DNA ligase (add last!)	1 µL
PCR product (cut) of each brick which is to be ligated – or 1 part plasmid and 5 part bricks	5 µL or 10 fmol Plasmid, 0, 10 and 20 fmol PCR
H2O	to reach a total volume of 20µL

--	--

## 9. Waste handling

Chemical name	Concentration	Type of waste (C, Z...)	Remarks

## 10. Time consumption

## 11. Scheme of development

Date / Initials	Version No.	Description of changes
13.06.19 / ASF	01	The SOP has been written
13.06.26 /PRA	01	The SOP has been approved

## 12. Appendixes

