

## iGEM EXPERIMENT 3 – CONSTRUCTION OF CiChE STRAINS

### Strategy

Transform plasmids constructed in experiment 2 in the *E. coli* strain constructed in experiment 1

### SOP

Make *E. coli* DendA DE3 Ch1.ccdA-Pmb1GFP-CmFRT electrocompetent and transform with plasmids pSB#X#T7ccdB

1. Incubate *E. coli* MG1655  $\Delta$ endA DE3 Ch1.ccdA-Pmb1GFP-CmFRT in 3ml LB+Cm o/n @ 37 °C
2. Add around 500  $\mu$ l of the overnight culture to 10 ml LB+Cm (in 50ml falcon tube) in order to have an OD<sub>660 nm</sub> between 0.05 and 0.1. Also add 500  $\mu$ l 1mM L-arabinose solution.
3. Incubate at 37 °C till OD 0.6
4. Place falcon tube on ice during 30 min.

**From now on it's important to keep everything cold! Use chilled tips and recipients!**

5. Centrifuge 4 min at 1250rcf (4 °C)
6. Resuspend the pellet in 45 ml cold MQ-water and centrifuge again during 5 min at 1250rcf (4 °C)
7. Resuspend the pellet in 1 ml cold MQ-water and transfer into a cold eppendorf
8. Centrifuge at 12000g during 20 - 30 seconds (4 °C)
9. Resuspend the pellet in 50  $\mu$ l cold MQ-water by moving the pipet point around (do not pipet up and down!)
10. Transfer 50  $\mu$ l of the suspension into a new cold eppendorf, what 's left over will be used as the control
11. Add 3  $\mu$ l of the plasmid DNA (100ng-200ng) and mix gently by turning:
  - a. pSB3T5-T7ccdB
  - b. pSB4A5-T7ccdB
  - c. pSB6A1-T7ccdB
12. Transform the strain by electroporation with the Gene Pulser:
  - a. Transfer the sample to the gene Pulser Cuvette using a cold tip
  - b. Tap gently to make sure that the sample is evenly distributed between the sides of the Cuvette. Also take care that everything is on the bottom of the cuvette.
  - c. Insert the cuvette into the pulse slider and close slide [pay attention to the correct orientation of the cuvette]

- d. Pulse once by pressing the two red button until the beep
- 13. Add immediately 1 ml SOC to the cells when the Gene Pulser beeps. Remove the cuvette from the pulse chamber and invert several times to mix. Transfer into a new culture tube.
- 14. Also add 1 ml SOC to what is left of the washed cells (= control, see above)
- 15. Incubate the transformed cells and the control during 2 h at 37 °C on a shaker
- 16. Transfer 50 µl of the culture + 50 µl SOC on a LB+Tet/Amp (340 µg/ml) plate.
- 17. Inoculate single colony in 3 ml LB+Tet/Amp and incubate o/n @ 37 °C
- 18. Archive strains
  - a. Create new strain number in sMEMO strain list
  - b. Add 0.5 ml o/n culture to 0.5 ml sterile 70% glycerol
  - c. Mix thoroughly and store @ -80°C