

## iGEM-Group – Microbiology – BMB – SDU

**Title:** Agarose gel DNA

**SOP number:** SOP0020\_v01

**Version number:** 01

**Date issued:** 2012.10.30

**Review date:** 2013.12.01

**Written by:** Michelle Madelung

### 1. Purpose

---

To visualise or purify DNA run on an agarose gel

### 2. Area of application

---

This procedure is valid for all DNA samples

### 3. Apparatus and equipment

---

Apparatus/equipment	Location (Room number)	Check points	Criteria for approval/rejection
Heating Cupboard	V18-403b-2	• Preheated	58 °C
Microwave oven	V18-403b-2	• Max strength	Never use metal in the microwave
Fume hood	V18-403b-2		
Power supply	V18-403b-2		
UV chamber (UVIDOC/TEC)	V18-403b-2		Use only when DNA is not to be used further
Transilluminator (UV light table)	V18-403b-2		Use for DNA samples to be purified
Erlenmeyer flask (500 ml) - glass			
Measuring cylinder (500 ml) - plastic			
Gel casting tray and comb	V18-403b-2		
Electrophoresis chamber	V18-403b-2		
Bull's eye spirit level	V18-403b-2		Centre the bobble in the bull's eye

#### 4. Materials and reagents – their shelf life and risk labelling

---

Name	Components (Concentrations)	Manufacturer / Cat. #	Room	Safety considerations
Tris-acetate-EDTA (TEA buffer)	0.4 M Trizma base 0.01 M EDTA / Titriplex III 0.2 M Sodium acetate trihydrate Acetic acid ≥ 99.8 %	Sigma / T1503-1kg MERCK / 1.08418.1000 Millipore Sigma-Aldrich / 33209-1L	10x Anne Mette	cas # 77-86-1   cas # 64-19-7
Demineralized water				
Agarose	SeaKem LE Agarose	Lonza	V18-403b-2	
Ethidium bromide	0.07 %	AppliChem	V18-403b-2	Can be carcinogenic – use gloves at all times
Photo paper		Mitsubishi electric	V18-403b-2	
Weighing paper		Contact lab-manager	BMB-storage	
Tesa Universal Tape		Tesa – 19 mm	BMB-storage	

## 5. QC – Quality Control

---

% agarose gel	Range of separation (bp)
0.3	5000 – 60000
0.5	1000 – 20000
0.8	800 – 10000
1.0	400 – 8000
1.2	300 – 7000
1.5	200 – 4000
2.0	100 – 3000

Size of comb	Amount of sample (μl)
Small thin	10
Small thick	20
Big thin	25-30
Big thick	50

## 6. List of other SOPs relevant to this SOP

---

## 7. Environmental conditions required

---

Ethidium Bromide can be carcinogenic and should therefore be handled with care. Always wear gloves in the agarose gel room and if something gets on the gloves change into a new pair.

## 8. Procedure

---

- 8.1 Look in table in paragraph 5 to identify percentage of gel to be used
- 8.2 Add appropriate amount of agarose to 500 ml Erlenmeyer flask
- 8.3 Add 300 ml buffer, mark water line on flask
- 8.4 Place a piece of weighing paper over the opening (NEVER METAL)
- 8.5 Place the flask in the microwave oven for 4-5 min
- 8.6 Check that all agarose has been dissolved or heat extra
- 8.7 Add demineralized water to the line marked before to replace the water that has evaporated
- 8.8 Mix by swirling the liquid in the bottle
- 8.9 Let the agarose cool to 60 °C
- 8.10 Add Ethidium bromide (1 drop per 60 ml) while swirling in fume hood
- 8.11 On the flask write your initials, percentage gel, the date and +/- Ethidium bromide
- 8.12 Place the flask in a heating cupboard (58 °C) until use
- 8.13 Attach tape at the open ends of the plastic agarose forms and insert a comb to form the wells in the agarose gel
- 8.14 Place the form on an even surface, check by using a bull's eye spirit level
- 8.15 Pour agarose into the form and let it set
- 8.16 Remove comb and tape and place the form in an agarose tub
- 8.17 Add buffer until the gel is submerged (this buffer can be used for three runs)
- 8.18 To DNA samples add loading buffer 1 µL per 5 µL sample, mix
- 8.19 When all samples are ready, load samples onto gel and add appropriate ladder to the first and last well if possible (paragraph 12)
- 8.20 Run the gel for 30 minutes at 75 V and unlimited amps
- 8.21 Take a photo of the gel in UV chamber (UVIDOC/TEC) or if the DNA is to be used further look at the gel on Transilluminator (UV light table) – run longer if necessary

## 9. Waste handling

---

Chemical name	Concentration	Type of waste (C, Z...)	Remarks
Used buffer		EtBr waste	
Agarose gel		GMO yellow waste	
Once use plastic		GMO yellow waste	

## 10. Time consumption

---


- Total-time 2 hours.
- Hands-on-time 1 hour.

11. Scheme of development

Date / Initials	Version No.	Description of changes
12.10.30 / MM	01	The SOP has been written
13.01.17 / MM & TK	01	The SOP has been approved

12. Appendixes


Fermentas GeneRuler™ 1 kb DNA Ladder Plus



Fermentas GeneRuler™ DNA Ladder Mix



Fermentas GeneRuler™ 50bp DNA Ladder



Fermentas GeneRuler™ 100bp DNA Ladder

