

JMJ-Group – Microbiology – BMB – SDU

Title: PCR with MyTaq

Date issued: 2013.06.18

SOP number: SOP0011_v01

Review date: 2013.06.18

Version number: 01

Written by: SIS & TJK

1. Purpose

To amplify wanted fragment of template via MyTaq PCR.

2. Area of application

This procedure is valid for all *E. coli* strains

3. Apparatus and equipment

Apparatus/equipment	Location (Room number)	Check points	Criteria for approval/rejection
PCR machine	Laboratory (class 1) – V18-403b-2	<ul style="list-style-type: none"> Check program 	Appropriate PCR program
Pipettes (p100,20,10)			
PCR tube rotator	Laboratory (class 1) – V18-403b-2		

4. Materials and reagents – their shelf life and risk labelling

Name	Components (Concentrations)	Manufacturer/Cat. #	Room	Safety considerations
Water	Demineralised milli-Q autoclaved water	Milli-Q water purification system (Millipore)	RT	
MyTaq TM HS Red Mix	http://www.bioline.com/documents/product_inserts/MyTaq%E2%84%A2%20HS%20Red%20Mix.pdf#zoom=130	Bioline	V18-405a-2	
Reverse primer	Made specific to the template	Sigma-Aldrich		
Forward primer	Made specific to the template	Sigma-Aldrich		
PCR tubes		Eppendorf	Micro storage	
1.5 ml tubes		Contact lab-manager	BMB storage	
Green pipette tips		Contact lab-manager	Micro storage	
Purple pipette tips		Contact lab-manager	Micro storage	

5. QC – Quality Control

For more than 2 samples; a premix of primers, MyTaq, primers and water is mixed and aliquoted before adding template.

6. List of other SOPs relevant to this SOP

iGEM2013_SOP007_v01_excision and ligation of PCR products in USER cloning

iGEM2013_SOP008_v01_preparing primers

7. Environmental conditions required

8. Procedure

- 8.1 Mix primers, water and MyTaq as described under PCR set up; paragraph 12 with the template
- 8.2 Be sure that all the samples are in the bottom of the PCR tubes by spinning on PCR tube rotator
- 8.3 Place in PCR machine
- 8.4 Start the appropriate PCR program for MyTaq or design one yourself, see paragraph 12
- 8.5 Keep at 4-5 °C until use, if more than 2 days waiting time place in -20°C

9. Waste handling

Chemical name	Concentration	Type of waste (C, Z...)	Remarks
Once used plastic		GMO yellow waste	

10. Time consumption

- Total-time 2 hours
- Hands-on-time 0.5 hour

11. Scheme of development

Date / Initials	Version No.	Description of changes
13.06.18 / SIS & TJK	01	The SOP has been written
13.06.18 / PRA	01	The SOP has been approved

12. Appendixes

PCR set up, 1 sample:

Notes: Smaller volume than recommended by Bioline

Template:	variable
Primers (10 μ M):	0.5 μ l each
MyTaq HS Red Mix, 2x:	5 μ l
Water (sterile):	3,5 μ l
Total:	10 μ l

PCR cycling conditions:

Step 1: Initial denaturation:	95 °C	2 min
Step 2: Denaturation:	95 °C	15 sec
Step 3: Annealing:	55 °C	15 sec (depending on the primer sequences, 2-5 °C below the lowest T _m of the primers)
Step 4: Extension/Elongation:	72 °C	10-30 sec (30 sec. pr. kilo bp)
Step 5: Repeat step 2-4:		30 times
Step 6: Extra elongation:	72 °C	2.5 min
Step 7: Keep the samples cold	4 °C	until the samples is removed