

Transformation of *E. coli* with plasmid DNA

Grow *E. coli* culture overnight (18-19 hours). Dilute fresh culture 20 times and incubate in shaker for 2 hours to the OD=0.3 – 0.4 (for the cotransformation use LB with the antibiotics).

1. Take 2 ml of culture to the centrifuge tube.
2. Keep on ice for 10 minutes.
3. Centrifuge for 3 min at 6 500 g at +4°C. Discard supernatant.
4. Resuspend the cell pellet in 0.8 ml of cold 0.1 M CaCl₂.
5. Keep the cells on ice for 15 min.
6. Centrifuge for 2 min at 6 500 g at +4°C. Discard supernatant.
7. Resuspend the cell pellet in 100 µL of 0.1 M CaCl₂.
8. Keep the tubes on ice for 4-5 min.
9. Add 10 µL (50ng) of plasmid DNA to the cells.
10. Keep on ice for 15 min.
11. Heat shock (put the tubes in water bath heated to 42 °C for 2 min).
12. Immediately put the tubes on ice and keep for 30 seconds.
13. Add 1 ml of LB to the cells.
14. Incubate at 37 °C for 1-2 hours (more than 1 hours for the cells, resistant to 2 or more antibiotics).
15. Plant the cells on plates with LB and corresponding antibiotics.