

NOVEMBER

Week 1

Thursday 05/11/13

Biobrick stocktake of 2013 iGEM Macquarie_Australia parts: 11/11/13

ChIG - sufficient plasmid stock

DVR1 - sufficient plasmid stock

ChIM - sufficient plasmid stock

ChII2 - need more plasmid stock

POR- need more plasmid stock

YCF - need more plasmid stock

Plasto - need more plasmid stock

Gun4- need more plasmid stock

CTH1 - need more plasmid stock

ChID- need more plasmid stock. Question whether the 2013 part is really the reported sequence - something appears to be missing.

Send all for re-sequencing to verify DNA sequence as per registry entries.

DVR1 re-transformation:

Was re-done using gibson assembly and then transformed.

NOVEMBER

Week 2

Tuesday 12/11/13

Sequencing Results: all parts except ChID were correct.

ChID Fix

ChID is missing 50 bp. Strategy to correct is to use Apal and Mlul restriction enzymes to cut out 50bp from clone of ChID in pET vector from Willows group and re-insert into our BioBrick vector.

Apal and Mlul were used in a single digest according to manufacturer's instructions and as per ligation protocol on methods wiki. Fragments run on 1% agarose and gel purified. However, digestions were incomplete as viewed on agarose gel. Need to do separate digests for next attempt.

Increase stocks

Did plasmid preps to get more of: ChII1; ChII2; YCF54; ChIP, DVR1; POR

ChIH Biobrick correction

Attempt to make ChIH (BBa_K1080001) using combination of gblocks and PCR products, as designed by Macquarie_Australia 2013 iGEM team.

Assembly strategy is: G-Block -1 (470bp) + PCR-1 (304bp) + G-Block-2 (499bp) + G-Block-3 (499bp) + PCR-2 (984bp) + G-Block-4 (500bp) + G-Block-5 (481bp) + PCR-3 (673bp)

Double restriction enzyme digest was carried out to combine PCR1 and Gblock2. After the two sections were ligated and extended, straight PCR was done. The PCR worked as judged by agarose gel.

November

Week 4

Tuesday

26/11/13

Composite parts Assembly

Biobrick (BB) ChII1 is combined with ChII2 biobrick in AMP backbone. Method is via 3A assembly. Use 500ng of each part and insert into 500ng of amp backbone. Ligation for 16oC for 30 mins then 80oC for 20 mins. Leave plates over weekend at room temperature.

Growth on plates : 1 colony on low plate, hundreds on high plate.

Assembly of ChIH

PCR of individual fragments from ChIH: 29/11/13

1. G1 - G1F + G1R
2. G2 - G2F + G2R
3. G3 - G3F + G4R
4. G4 - G4F+ G4R
5. G5 - G5F + G5R
6. G6 - G6F + G6R
7. PCR1+ G2- H1F+ G2R
8. (PCR1 + G2) + G1
9. G3 + PCR2- G3F+ H2R
10. G5 + G6- G5F +G6R

ChIH G1 and G4 fragments were amplified but the rest of the bands were not clear and were likely to have failed.

Friday 29/11/13

ChID: another attempt at PCR using fragments: 29/11/13:

G block 1


G block 4

G4 + (G5-G6)

G1 + PCR1

G2 + (G3+PCR2)

Extremely faint bands seen for ChID amplification. G1 and G4 appear to work but bands are very faint on agarose gel.



Faint to no bands viewed for G2 & G3+PCR2. Re-attempt necessary.


DECEMBER

Week 1

Friday: 06/12/13

Digestion & Ligation of ChIM gene of lac promoter into CAM backbone

ChIM in AMP backbone vector and lac in backbone were digested using iGEM restriction digestion protocol EcoRI and Pst1 restriction enzymes.



Small amount of growth seen on ChIM, lacA & lacB indicating that they were successfully incorporated into the DHS← cells. Sequencing to confirm required.

Transformation of Kanamycin Resistant backbone

We need more of the kanamycin biobrick. Transformation of kanamycin backbone into *E. coli* cells to produce large amounts of KAN backbone for future ligations.


Monday
09/12/13

PCR reaction for ChIH and ChID

The overall aim of the week was to build ChIH fragment and PCR ChID. Using the standard PCR protocol, G1+H1, G2 (G3+H2), G4 (G5+G6), ChID (2) and ChID (3) were run.

The result showed another G1+PCR1 failure. It was also suggested however to use BioBrick primers. Distinct bands for G2+(G3/PCR2) and G4+(G5/G6) were present and proved correct. This assumption was made that these results were correct.

ChID 2 and 3 showed a band present at approximately 1500 bp which was also assumed to be correct in relation to the actual size of 1681 bp.



The next step was to rePCR G1+PCR1 with BBF + HR2 and BBvF + HR2, gel extraction of G2+(G3/H2) and G4+(G5/G6), ChID 2 and 3.

Tuesday
10/12/13


The fragments to be PCR'd and the primers are presented on the following table;

Fragments to PCR	Primers
G1	BBF+G1R
G1+P1	BBF+H1R
G2+(G3-P2)	G2F+P2R
ChII1	BBVF2+BBVR
ChII2	BBVF2+BBVR
ChID	BBVF2+BBVR

The standard PCR method was adopted to run the reaction.

Results of the PCR reaction;

All but ChII 1 and 2 failed



Continued ChIH construction

At this stage, the ChIH gene construct was continued;
 G1-P1-G2-G3-P2-G4-G5-G6


The ChID gene was cut from the gel and extracted with another attempt to PCR.

To test for protein expression, the successful 3A gene was combined with lac creating a composite.

Continuing the construct of ChIH, a PCR reaction was performed to identify the successful or unsuccessful attempt in the composite build in addition to DVR1 identification.

Gene fragment	Primers
P1+G2	H1F+G2R
(G3+P2)+(G4-G5-G6)	GBF+G6R
DVR1	BBVF2+BBVR

The PCR reaction screening the attempting to construct ChIH failed.




Digest of DVR1

The next step was the insertion of DVR1 into the plasmid vector. The plasmid vector and the plasmid containing the gene of interest were ligated with EcoR1 and Pst1. The gene was introduced into the vector by means of 1 vector to 3 insert to maximise insertion efficiency.

ChIH construction by Gibson assembly

The failure of the construction of the ChIH gene subjected the attempt in the construction of the gene using Gibson assembly. The provided gel image proved the construction also failed.



Tuesday
07/01/14

ChIH construct PCR


In the attempt to yield a positive result in the construction of ChIH, each P1+G2, (G3+P2) and (G4-G5-G6) were PCR'd separately in the attempt to successfully join the individual components.

Gene fragment	Primer
P1-G2	P1F + G2R
G3-P2	G3F + P2R
G4-G5-G6	G4F + G6R
(P1-G2) + (G3-P2)	P1F + P2R
(G3-P2) + (G4-G5-G6)	G3F + G6R
ChID	ChID F + ChID R

Thursday
09/01/14

Gel analysis of PCR gel of ChIH constructs and ChID

The results obtained would indicate the band to extract for Gibson assembly. The gel image showed positive results



The marked were cut out and stored for gel extraction.

To compare the sizes, 25-500ng of plasmid were digested with and without lac. The expected size was approximately 200 bp.

The amplification of ChID was faint indicating an issue with the construction of the gene.

Friday
10/01/14


ChIH screening

The bands on gel corresponding to ChIH were extracted to screen for the correct sizes. The result of the gel extraction showed low concentration indicating poor construction of gene.

ChID, ChII and Plastocyanin screening

Digests were performed with enzymes EcoR1 and Pst1 to comment on the sizes of the inserts including ChID, ChII1 and Plastocyanin. These were also run against the corresponding components including lac.

Results showed that the success of ChII1 and plastocyanin + lac. ChID and ChII1/2 however showed negative results.




Saturday

11/01/14

Continued construction of ChIH PCR

Gene fragment	Primers
P1-G2	F1 + G2R
G3-P2	G3F + P2R
G4-G5-G6	G4F + G6R
CHID	NF2 + NR2

The results obtained from the gel yielded a successful result for P1-G2, responsible for the construction of CHIH and negative results for the remaining samples on the gel.



JANUARY

Thursday


30/1/14

PCR: It is thought that the excess template in the previous PCR may have been responsible for the failure of PCR amplification. Template dilutions of 1/10 and 1/100 were tested by running another pcr.

The PCRs carried out were:

- ChID (new template), diluted
- ChIH G3 -P2 PCR template
- ChIH gel run + extracted template

Results:




ChID and ChIH G3-P2 PCR template did not work, however, ChIH G3 -P2 PCR template was successful.

PCR for ChID blocks:

G4 + (G5-G6) X3 = G4F + G6R

P1- G2 (from the original templates) x3 = P1F +G2F.

Results:




G4 + (G5-G6) and P1- G2 PCRs worked. G4 + (G5-G6) showing a band of 1700 bp in length and P1-G2 showing 800 bp in length.

PCR continuation:

(P1-G2) + (G3-P2)

G3-P2 + G4-G5-G6

Results:



None of the PCRs from the ChID blocks worked.

FEBRUARY

Week 1

Monday

3/2/2014

Digestion of DVR1 was run.

Following the digest, a ligase reaction was conducted and transformation performed. Note, the concentration of DVR part in comparison to the concentration of the plasmid was 1.5 times more. Plates incubated overnight.

Protein expression of lac+plasto & lac+ChlI

Expression of protein via lac promotor using 2uL of IPTG was done for each sample to amplify protein expression. SDS-PAGE was run according to methods.

Tuesday

4/2/14

SDS PAGE attempt #2

Here we conducted a second SDS PAGE for lac plasto and lac ChlI.


Lane order: 1-4 lac plasto, 5 is the ladder, 6-10 lac ChlI

Expression of proteins was not visible by eye. No image of the gel was recorded. We think we need to do mass spec (MALDI/TOF/TOF) to identify proteins in bands. Discuss with APAF (Australian Proteomic Analysis Facility) at Macquarie University to ask if they can help us with performing mass spec.

Testing new ligase: new ligase purchased as concerns were that our ligase was old and the reason ligations were not successful

ChIH was digested with E+P restriction enzymes as per methods. To ligate, the ligation mixture comprised of 8.5uL DNA, 0.5uL ligase and 1uL of buffer.

Results:



The gel showed inconclusive results and requires further clarification. Further testing methods such as re-inserting the biobrick into another vector and growing it on a plate with the second vector antibiotic have been suggested.

PCR:

Fragment 1 - (P1-G2) + (G3-P2) = P1F, P2R

Fragment 2- (G3-P2) + (G4-G5-G6) = G3F, G6R

Fragment 3- (P1-G2) + (G3-P2) + (G4-G5-G6) = P1F, G6R.

For such large fragments, the preliminary melting step was completed twice prior to the addition of the primers because of the long fragments. The rest of the process was continued on the regular loop as in other PCR protocol.


Wednesday

5/2/14

Western Blot:

Western blot for ChII1 and plasto were carried out.

Results: The plasto lanes did not show any expression, however ChII showed good expression in lanes 2,4 and 5.



Thursday 6/2/14: Gibson Assembly of ChIH:

G1: 3uL

P1-G2 exosap: 0.5uL

G3-P2 exosap: 4.8uL

G4-G5-G6 exosap: 1.0uL


Cam vector: 29uL with Gibson mix: 12.3uL

or AMP vector 1.4uL with Gibson mix: 10.8 uL

- Plated out

PCR of Gibson Assembly product for ChIH: Standed PCR x2 using BBF/ BBR/ BBVF2, BBVR.

Results:



February

week 2

Wednesday

12/2/14

Nanodrop of ChIH fragments

ChIH Fragments:

P1-G2:

- A= 19.3 ng/ml
- B= 23.4 ng/m
- C= 18.2 ng/m

G3-P2:

- A = 27.8 A= 141ng/ml
- B= 9.9 ng/ml
- C= 47.5 ng/ml

G4-G5-G6:

- A= 141ng/ml
- B= 24.6 ng/ml
- C= 62.5 ng/ml

ChID Fragments: D1, D2, D3

Thursday

13/2/14

Gel Electrophoresis for ChIH and ChID fragments:

Top Gel Lane order: 1- Ladder, 2- G1, 3-5 - P1-G2, 6-8- G3-P2, 9-11- G4-G5-G6, 12-14- ChID

Bottom Gel Lane Order: 1+ 2- Ladder, 3-5- P1-G2, 6-8- G3-P2, 9-11 G4-G5-G6

Results: ChIH fragments appear not have been digested. P1-G2 and G4-G5-G6 didn't have plasmids on the gel so they did not digest. ChID has digested with MLU and partially APAI. ChID plasmids from lanes 12 and 13 were added together for further re-digestion.

New Digestion using E+P from previous gel electrophoresis for ChIH:

Lane Order: 1- G1, 2- G3-P2 **A**, 3- G3-P2 **B**, 4- G3-P2 **C**, 5- G4-G5-G6 **A**, 6- G4-G5-G6 **C**, 7- Plasto Control

New Digest for ChID with APAI: The ChID being redigested is a combination of lanes 12 and 13 from the previous electrophoresis gel.

ChID Ligation : ChID was ligased and transformed into E.coli. 2 colonies grew on the 300uL plate and 1 colony on the 30uL plate.


Friday 14/2/14

Plasmid nanodrop

Lac	Gun	A	33.8 ng/ μ l
		B	17.8 ng/ μ l

Lac	ChIM	A	33.5 ng/ μ l
		B	26.8 ng/ μ l
Lac	ChIP	A	32.2 ng/ μ l
		B	14.4 ng/ μ l
Lac	ChI2	A	18.5 ng/ μ l
		B	15.2 ng/ μ l
Lac	POR	A	35.1 ng/ μ l
		B	24.5 ng/ μ l
Lac	ChIG	A	8.3 ng/ μ l
		B	13.2 ng/ μ l
Lac	YCF54	A	29.3 ng/ μ l
		B	33.9 ng/ μ l
Lac	CTH1	A	59.1 ng/ μ l
		B	51.6 ng/ μ l

ChIH re-digest:



Digests checking biobricks:

With/without lac	Biobrick	Fragment	Expected weight	Measured
Lac	Gun	A	930	~930
		B		~930
Lac	ChIM	A	873	~1050
		B		~1050
Lac	ChIP	A	1299	~1500
		B		~1500
Lac	ChI2	A	1212	~1400
		B		~1250
Lac	POR	A	1067	~1250

		B		~1250
Lac	ChIG	A	1050	~1250
		B		~1250
Lac	YCF54	A	471	~650
		B		~650
Lac	CTH1	A	1152	~1350
		B		~1350
	DVR1	A		~1106
		B		~1106
		C		~1106
		D		~1106

Results: Majority look as though they match the expected band length. All digests excluding DVR1 include lac.

AUGUST (SEMESTER 2)

WEEK 1

Thursday

07/08/14

Off and running with the whole team of 12 Biomolecular Major students. We sat through a full day of learning about iGEM; we had a discussion of project goals and aims; as well as a refresher course on how the Chlorophyll pathway works. Roles were assumed by our wiki-chiefs and those interested in gaining sponsorship and promotional roles were also filled. The wet lab group started to discuss the plan for the first wet lab next week as well as looking over the protocols.

WET LAB WEEK 2

Thursday

14/08/14

Project Name: After brainstorming many names, we decided on "The Green Machine" as our title and the slogan "Follow the biobrick road" as our theme to carry throughout our wiki page.

Stocks: Made many stocks, plates and buffers as per methods.

Nanodrop: leant how to use the nanodrop to quantitate all parts

Gene Info: We did another stock-take of parts & checked that we had enough to perform ligations to assemble or planned three Operons. We discussed the strategy for how we were going to make each of the three Operons. The parts we require to assemble our pathway are as follows:

ChID - 2240bp

ChII1 - 1202bp

ChII2 - 1298bp

GUN4 - 782bp

ChIH - 4207bp

CTH1 - 1382bp

YCF54 - 556bp

Plasto - 410bp

ChIM - 959bp

POR - 1154bp

DVR1 - 1193bp

ChIP - 1385bp
ChIG - 11366bp

ChID BioBrick Correction

As before, previous BioBrick (BB) from 2013 had a 50bp deletion/error within the ChID (900bp) from using a single restriction digest. Our Aim is to excise the entire ChID gene using Apa1 and Mlu1 restriction enzymes and to insert a complete ChID gene into a BB. We attempted this experiment again but did restriction enzyme digestions with the two enzymes separately to improve efficiency of cutting.

To prevent re-joining after digestion of our cut vector, we treated our samples with Alkaline Phosphatase (Fast A.P.). The DNA was then run on a 1% Agarose gel. 5 bands were identified and using a 1Kb ladder a complete ChID (900bp) band was found and excised for ligation.

WET LAB WEEK 3

Thursday

21/08/14

More plasmid prep was done, the following 6 genes were inserted into an Ampicillin backbone. Cells were grown to extract more plasmid for stocks.

ChI1
ChI2
YCF54
ChIP
DVR1
POR

Composite Part Assembly: Trouble-shooting with the BioBrick assembly protocol, we found that if we ligated in a particular way then the plasmid linearises itself and then cannot be cut for the making of composite parts. We then started working on forming test composite parts. Our stocktake was also completed.

Transformations: Electroporation does not appear to be working well. We changed to heat shock to transform our cells. There may be a problem with our electro-competent cells. Made more electro-competent cells to test. Also made chemical competent cells for heat-shock transformation

DRY LAB WEEK3

Thursday

21/08/14

Decided on Outreach ideas. FINALLY. Online reality contest "So You Think You Can Synthesise". Had discussions of framework for competition, making a trailer.

Met up with MQ Media Team later during the week

Started to put together the sponsorship package

WET LAB WEEK 4

Thursday

28/08/14

Digests

Digests of GUN4+ChII2 & ChID to check results from last week

Competent cells

Electroporation does not seem to work and create viable competent cells ergo we shall stick to the heat shock methodology for further preps. More cells were made and used for plasmid preps, BB's and composite parts.

Composite Parts

4:1 insert – vector ratio for Fast AP ligation steps to produce:

AMP backbone

CTH1 + YCF54

CTH1 + Plasto

ChIP + ChIG

CAM backbone

ChID

KAN backbone

GUN4 + ChII2

GUN4 + ChII1

ChII1 + GUN4

WET LAB WEEK 5

Thursday

04/09/14

A busy week!

Composite parts that had growth were digested with X & P and run on 1% Agarose gel to check insert size.


PCR was also performed with BioBrick Forward and BioBrick Reverse primers to see if we could confirm correct assembly of composite parts.

further composite parts were assembled on the AMP backbone


ChIM + YCF54

POR + DVR1

POR + ChIP



Composite part RE digest images from plates that had growth



Friday

05/09/14

Liquid cultures of composite part transformants

Plasmid preps

RE digest of each part – into CAM BB

Competent cell prep: both chemical and electro-competent cells were made

New composite parts made:

/ChIM+YCF54/

/POR+DVR1/

/POR+ChIP/

WET LAB WEEK 6

Wednesday 10/09/14

Ran PCR products from last week

ChID

CTH1 + YCF54

CTH1 + Plasto

ChII2 + GUN4

GUN4 + ChII2

Plasmid preps done

Thursday 11/10/14

Sequencing

Plasmids were positive, sent to Macrogen for sequencing:

CTH1 + YCF54

ChID

Gun4 + ChII1

CTH1 + Plasto

Composite Checks


These composites were cut as a single digest and a double digest, and run on an agarose gel. Sizes were compared,

POR + ChIP


POR + DVR1

ChIM + YCF54

CTH1 + Plasto



Composite part screenings



Composite Parts

New composite parts were made, and built upon, transformed and plated out:

/CTH1+YCF54/ + Plasto

/CTH1+YCF54/ + ChIM

ChII2 + ChII1

ChII2 + GUN4

Open Day

Saturday 13/09/14

Set up chromatography reactions in preparation for open day on Saturday.

Fluorescent plates drawn, grown, ready to go for Saturday.

Plasmid prep done for previous composite parts. CHIH has finally worked

WET LAB WEEK 7

Monday

15/09/14

Nanodrops of plasmid preps.

single and double digests of every second plasmid (a, c, e, f) run gel. Send for sequencing if successful.

Wednesday

17/09/14

gels re-labelled

transformations

CTH1 + YCFS4 + Plasto <- ChIM

GUN4 + ChID + CHII2

Thursday

18/09/14


New Composites

- POR + ChIP + ChIG
- POR + DVR1 + ChIG
- POR + DVR1 + ChIP
- ChID


All transformed and plated out.

ChID Fix

Apa1 & Mlu1 digests with the backbone being treated with Fast AP in Mlu1 digest reaction. Ligations was as usual. Then transformed and plated out onto CAM plates. Gel digests were run to resolve the 50bp difference (850-900). The resolution was seen.



^ Colony screen apal mluI ChID 850 and 900 + ChIH



^ pET Apal MluI Digest Gel

ChI H PCR reaction was also run

WET LAB – MIDSEM BREAK W1


Monday

22/09/14

Liquid Cultures from Thursday plates.

- POR + ChIP + ChIG
- POR + DVR1 + ChIG
- POR + DVR1 + ChIP
- ChID

Restriction Enzyme Digest of CTH1 + YCF54 + Plasto + ChIM + ChID composite part.




Gel run of ^ composite part and ChID to compare with pET to confirm to presence of the 50bp.

Tuesday

23/09/14

EcoRI + X/P digests for plasmids from Mon.

ChID A + M individual digests, compared against pET with same digest.



Re-screen CTH1 + YCFS4 + Plasto + ChIM colonies and grow in liquid culture

**Wednesday
24/09/14**

Recheck: ChIH in KAN and CAM both were not okay.
CTH...ChIM: not okay therefore re screen plates from liquid cultures.

Composite part-> POR + DVR1 + ChIP + ChIG
Transform: ChII2 + GUN4 + ChID} new composite and plate out.

Rescreen plated samples of CTHI + Plasto & ChIM + YCFS4 all in AMP and put into liquid culture.

**Thursday
25/09/14**

Liquid culture of ChII2 + GUN4 + ChID

Plasmid prep done in afternoon

POR ...ChIG transformed and plated out
CTHI ...ChIM gel resolved and excited for friday purification.
ChII2 + GUN4, CTH1 + YCF54 and POR + DVR1 into CAM backbones.

prepare ChII2 + GUN4, CTH1 + YCF54, ChID for sequencing
Plasmid prep of CTHI + Plasto, ChIM + YC5S4.

**Friday
26/09/14**

CTHI ..ChIM gel band purified
Liquid cultures of new composite parts (ChII2 + GUN4 + ChID POR..ChIG)

Plasmid prep of CTHI + plasto (A-D) and ChIM + YCF54 (A-D)
Restriction enzyme screen plasmid prep from thursday.

Transformation of GEL extracted plasmids (linear [total of 10 for transformation] + circular + ligation)
*1ul ligase and 4.5 ul ligase buffer.
37oC for 1 hour and 80oC for 20mins.

→ 5ul linear plasmid → 50ul competent cells (chemical)
→ 10ul circular plasmid → 50ul competent cells

Gel Run for ChII2 + GUN4 + ChID, ChIH1 + Plasto & ChIM + YCF54 cuts.

LANE ORDER:

1.1....13.2- ChII2 + GUN4 + ChID - ChIH1 + Plasto - ChIM + YCF54

*note: The wells for the last 3 lanes did not accept much of the sample. The sample would float to the surface when being inserted.

The transformants were plated out, the digest resolution wasn't clear and needs to be re run on monday.
liquid cultures of Sunday transformants 28/9 (CTH1...ChIM gel purified CTH1 + YCF54 (CAM), ChII2 + GUN4 (CAM)).

WET LAB – MIDSEM BREAK W2

Monday

29/09/14

plasmid preps of:

ChIM CAM

CTH1..ChIM gel purified

CTH1..ChIM re-screen

POR..ChIG

POR + DVR AMP

CTH1 + YCF CAM

ChII2 + GUN4 CAM


Re Digest and Gel which all results were good:

ChII2 + GUN4 + chID

CTH1 + plasto AMP

ChIM + YCF54 AMP

Gel:



ChIH linearized and gel purified.

Digest + gel: CTH1..ChIM and POR...ChIG

CTH..ChIM plasmid prep kept to H, I, J

POR ..ChIG plasmid prep kept to 1:1, 1:3, 1:6

ChIH upper and lower bands excised & gel purified

Tuesday
30/09/14

Re-transformed and plated out:

CTH..ChIM x 3
POR..ChIG x 3
ChIH x 2
CTH..+ POR.. mixed x 2

Gel run on CAM transformants and re-screen of:

POR + DVR, ChII2 + GUN4, CTH + YCF

Composite part created, transformed and plated out: ChII2 + GUN4 + ChID + ChII1
ChIM cut into CAM Backbone, transformed and plated out. Tested for registry

PCR reactions run on final constructs and checked for

F & R from ends:

CTH...ChIM
POR...chIG

Wednesday (1/10/14)

PCR:

1ul BB
5ul buffer (x10)
1ul F primer
1ul R primer
1ul dNTP
0.25 ul Taq
0.75 ul H₂O (to Ebel 50ul)

Master Mix:

50ul buffer
10ul dNTP
2.5ul Taq
407.5ul H₂O

----- X -----

C1 + C2 → CTH1 + ChIM (H)
(CTF + BBR) (CMR + BBR)

C3 + C4 → CTH1 + ChIM (I)
(CTF + BBR) (CMR + BBR)

C5 + C6 → CTH1 + ChIM (J)
(CTF + BBR) (CMR + BBR)

P1 → POR A P2 → PORB P3 → PORC

-----X-----

Setup for functional assays and plasmid preps of ChIM & ChII1 + ChII2

Thursday

2/10/14

Cyclase assay (100ul)
10uM MPE - 8ul
1mM NADP - 10ul
10mM G-P-P - 10ul
Assay buff 1 x (50mM tricine, 2mM MgCl₂, 1mM DTT 10% glycerol, pH 8.0) - 50ul 2x

0.5ul of G-6-P-dehydrogenase
total volume = 78ul allowing 22ul for other additions.

- #1 Blank (water)
- #2 22ul CTH1 part 1
- #3 11ul CTH1 part 2
- #4 11ul CTH

Bradford functional assays were done on induced cell pellets
10% Glycerol + 5mM Tricine NaOH ph8.0 + 2mM MgCl₂ + 1mM DTT
Results:


		5µl	2.5 dilution
POR - ChIG	4	~ 1.25	1
	3	~ 1.75	1.5
	2	~ 2	1.5
	1	~ 2	1.5
CTH - ChIM	1	~ 1	1
	2	~ 1	1

Friday
3/10/14


Protein weight estimates:

- ChIM - 30440 Da
- CTH1 - 43873.3 Da
- YCFS4 - 17073.7 Da
- Plasto - 10339 Da
- Chli 1 - 39952 Da
- GUN4 - 2450.6 Da
- ChID - 76420.1 Da
- POR - 41871 Da
- DVR1 - 37034 Da
- ChIP - 47011 Da
- ChIG - 36880 Da

Protein gels run of two complete composites.



^CTH1-ChIM composite



^POR-ChIG composite

The CTH1-ChIM composite showed good separation of products. The POR-ChIG composite had separation but only three parts were able to be easily identified. As Annotated the selected bands were cut-out for in-gel digestion & analysis by MALDI- TOF/TOF.

WET LAB WEEK 8

Wednesday


8/10/14

Plasmid preps were done on the Chli1 - ChliD composite parts

Nanodrops:

Sample	nucleic acid conc (ng/ μ l)	260/280
Chli1 + ChliD	179.6	1.92
Chli1 + ChliD	296.3	1.87
Chli1 + ChliD	261.7	1.93
Chli1 + ChliD	648	1.90
Chli1 + ChliD	118.7	2.07
Chli1 + ChliD	179.1	1.79
Chli1 + ChliD	553.7	1.93
Chli1 + ChliD	321.6	1.92

All 8 were then digested and run on a gel, 5-10mL liquid cultures were also prepared and left to incubate overnight.



5-10mL liquid cultures were made of the following parts in an AMP backbone with a *lac* promoter. Chli2, YCFS4, ChIP, POR, GUN4, ChIM, ChIG, CTH1, CTH1-ChIM, POR-ChIG for large scale growth (50mL) for functional assays.

Thursday (9/10/14)

New composite Chli1 + ChID + GUN4 transformed and plated out. Intermediates and final parts prepped and sent for sequencing to confirm that the inserts are what they're supposed to be. Glycerol stocks were made of the current intermediates and final composites. Final parts in an AMP backbone were induced ($OD_{600} = 0.4-5$) for functional assays.

WET LAB WEEK 9


**Monday
13/10/14**

Plasmid prep of Chli1+ChID+GUN4, All cultures were screened and prepped for sequencing.


Nanodrops:

Sample	nucleic acid conc (ng/ μ l)	260/280
Chli1+ChID+GUN4	539.2	1.86
Chli1+ChID+GUN4	315.4	1.91
Chli1+ChID+GUN4	493.7	1.90
Chli1+ChID+GUN4	164.0	1.91
Chli1+ChID+GUN4	217.5	1.90
Chli1+ChID+GUN4	447.0	1.90
Chli1+ChID+GUN4	499.2	1.90

All composites re-run on gels for results page, SDS-PAGE gels and MS/MS prep. POR-ChIG and Chli1-GUN4 parts were re-transformed for lysate harvesting.



^ Chli - ChIM Final gels



^ POR-ChIG Final gel

Tuesday
14/10/14


POR assays were run with GUN4 acting as a negative control. We expected to see a major peak at $\sim 630\text{nm}$ and a secondary peak at $\sim 670\text{nm}$, but on both runs (a 20minute and 1.5hours) the second peak was still not visible. Cyclase assays were also run mirroring experimentation from 2/10/14 on CTH1 and POR, these were run overnight with an expected peak at $\sim 590\text{nm}$ showing presence of the Mg - protoporphyrin intermediate. The Mg is apparent but no intermediates have been generated.

Wednesday
15/10/14


Liquid cultures were made from plate cultures from the previous day (POR-ChIG, CTH1, Chli1-ChID-GUN4, POR). The liquid cultures from the POR-ChIG & Chli1+ChID+GUN4 composites were then induced for growth as 50mL cultures for French Pressing. The resulting proteins were then run on SDS-PAGE gels and used for functional assays; gels were

destained and bands cut for in-gel digestion and MALDI-TOF/TOF.

**Thursday
16/10/14**



^ CTH1 - ChIM Final gel



^ POR-ChIG Final Gel

Friday (17/10/14)

LAB BOOK - PROTOCOLS

PLASMID PREP

1. Centrifuge @13200 rpm for 10min 1.8mL of overnight cultures in 2mL eppendorfs
2. Discard supernatant and add 1.9mL of culture and centrifuge again @13200 rpm for 10min
3. Resuspend pelleted cells in P1(250µl)
4. add 250µL of P2 invert 4-6 times (turns homogenous blue)
5. Add 350µL of N3 & mix by inverting (turns colourless)
6. Centrifuge @ 10min 13000 rpm to create a pellet
7. Transfer supernatant in QIA Prep Spin Column by pipetting
8. Centrifuge 30-60sec - discard flow through
9. Wash Qia Prep Spin Column with 0.5mL of PB
10. Centrifuge for 30-60 seconds - discard flow through
11. Wash spin column by adding 0.75 mL PE buffer
12. Centrifuge for 30-60 seconds -> discard flow through and centrifuge for a further 1min
13. Place QIAprep Column in clean eppendorf
14. Elute DNA add 50µl of water, stand for 1min, centrifuge for 1min
15. QIA prep -> Spin miniprep buffer

BIOBRICK ASSEMBLY

Digest

1. Prepare the following
 - 1.1. Upstream Plasmid Prep (U)
 - 1.2. Downstream Plasmid Prep (D)
 - 1.3. Destination backbone (P)
 - 1.4. NE Buffer 2
 - 1.5. BSA
 - 1.6. 3 PCR tubes - labelled (U/D/P)
2. Add 250ng of each part to its respective tube
 - 2.1. adjust total volume to 21.25µl with water
3. Add 2.5µL of NE Buffer 2 to each tube
4. Add 0.25µL of BSA to each tube
5. RE digest - total volume ~ 50µl
 - 5.1. (U) 0.5µL EcoR1 - HF + 0.5µL Spe1
 - 5.2. (D) 1µL Xba1 + 1µL Pst1
 - 5.3. (P) 1µL EcoR1-HF + 1µL Pst1
6. Mix & Spin Down
7. Incubate @ 37°C ~ 15min then @ 80°C for 20min
 - 7.1. [OPTIONAL] Run a gel with 20µL of each
 - 7.2. [OPTIONAL] Store @ -20°C

Ligation

1. Prepare the following
 - 1.1. 10x T4 DNA ligase Buffer
 - 1.2. T4 DNA Ligase
 - 1.3. PCR tube (L)
2. Add 11µL of water to the L tube
3. Add 2µL from each digest tube (U/D/P)
4. Add 2µL of 10x Buffer
5. Add 1µL of T4 Ligase total volume = 20µL
6. Incubate @ room temperature for 10min
7. Incubate @ 80°C for 20min
8. Store @ -20°C or transform.

FAST AP LIGATION

HEAT-SHOCK TRANSFORMATION

1. Add 2 μ L ligation reaction to chemically competent cells and mix
2. Incubate on ice for 30min
3. Heat shock @ 42°C for 30sec
 - 3.1. Store on ice for 2min
4. Add 950 μ L of SOC medium
5. Incubate @ 37°C for 1hr with shaking
6. Plate out 5-200 μ L onto LB plates w/ antibiotic required & IPTG if necessary
7. Incubate overnight @ 37°C

ELECTROPORATION TRANSFORMATION

1. Add 1 μ L ligation reaction to a fresh PCR tube
2. Add 50 μ L electroporated cells
3. pipette mix into an electric cuvette and pulse
4. Add 1 μ L warmed LB broth to cuvette and mix gently
5. Quickly transfer contents to a fresh tube
6. Incubate @ 37°C for 1hr.