

Transformation

Protocol:

1. Take competent cells out of -80°C. Thaw the competent cells (100 µl) **on ice** for 5 minutes.
2. Add 5 µl of ligation product to the competent cells and mix gently by stirring the mixture.
3. Incubate the cells on ice for 20 minutes.
4. Heat shock at 42°C for 30 seconds.
5. Put the tubes back on ice for 5 min.
6. Add 200 µl of LB medium into the tube.
7. Incubate the mixture in 37°C shaking incubator for 60 minutes. Shake vigorously (250 rpm).
8. Plate 2 µl, 20 µl and the rest of each cell culture onto three LB agar plate containing appropriate antibiotic respectively*.

* The concentrations corresponding to the antibiotics used

Antibiotics	Concentration (µg/ml)
Ampicillin	100
Chloramphenicol	34

9. Incubate the plates at 37°C overnight.
10. Pick colonies.