

iGEM TU/e 2014
Biomedical Engineering



Eindhoven University of Technology
Room: Ceres 0.04
Den Dolech 2, 5612 AZ Eindhoven
The Netherlands
Tel. no. +31 50 247 55 59
2014.igem.org/Team:TU_Eindhoven

FACS - DBCO-PEG₄-5/6-TAMRA

Table of contents

Title		
FACS - DBCO-PEG4-5/6-TAMRA	1	Stock solutions
	2	Preparation of FACS samples

1 Stock solutions



DBCO-PEG₄-5/6-TAMRA, molecular weight: 936.06 g/mol; λ_{abs} 545 nm; λ_{em} 565 nm

- DBCO-PEG₄-5/6-TAMRA in DMSO (5 mM)
- DBCO-PEG₄-5/6-TAMRA in DMSO (833 μM)
- Buffer: PBS

2 Preparation of FACS samples

- Prepare following tubes:

Tube	[DBCO]	Buffer volume to add	Cells ($10^9/\text{mL}$)	DBCO volume to add (μL)		DBCO/tag ratio
				5 mM	833 μM	
1	0	150 μL	50 μL			
2	30 μM	0 μL	200 μL	1.21		1,1E+04
3	5 μM	0 μL	200 μL		1.21	1,8E+03
4	30 μM	100 μL	100 μL	1.21		2,2E+04
5	5 μM	100 μL	100 μL		1.21	3,6E+03
6	30 μM	150 μL	50 μL	1.21		4,4E+04
7	5 μM	150 μL	50 μL		1.21	7,3E+03

- React DBCO tubes for 1h or 6h a in shaking block at 4°C and 500 rpm
- Prepare FACS samples after 1 and 6 hours:
 - Spin down the cells for 5 min at 13,400 rpm
 - Resuspend with 1 mL ice cold PBS
 - Spin down the cells for 5 min at 13,400 rpm and put on ice until FACS
 - Right before FACS: resuspend with 200 μL PBS