

iGEM TU/e 2014
Biomedical Engineering

Eindhoven University of Technology
Room: Ceres 0.04
Den Dolech 2, 5612 AZ Eindhoven
The Netherlands
Tel. no. +31 50 247 55 59
2014.igem.org/Team:TU_Eindhoven

FACS - Antibody Titration

Table of contents

Title			
FACS - Antibody Titration	1	Stock solutions	3
	2	Preparation of FACS samples	3

1 Stock solutions

- HA Anti Mouse Monoclonal Antibodies in PBS (4.52 µM)
- HA Anti Mouse Monoclonal Antibodies in PBS (1 µM)
- HA Anti Mouse Monoclonal Antibodies in PBS (200 nM)
- DBCO-PEG₄-5/6-TAMRA (5 mM)
- Buffers: PBS (DBCO) or PBS-0.1%BSA (antibody)

2 Preparation of FACS samples

- Prepare following tubes:

Tube	[DBCO]	Cells (10^9)	DBCO volume to add (µL)		DBCO/tag ratio
			5 mM		
1	0	200 µL			
2	30 µM	200 µL	1.21		182.2

- React DBCO tubes for 1h in shaking block at 4°C and 500 rpm
- Prepare FACS samples:
 - Spin down the cells for 5 min at 13,400 rpm
 - Discard the supernatant and resuspend with 1 mL ice cold PBS
 - Spin down the cells for 5 min at 13,400 rpm
 - Discard the supernatant and put the pellets on ice until FACS
 - Right before FACS: resuspend with 200 µL ice cold PBS
- Prepare following tubes:

Tube	[Ab]	PBS-0.1%BSA	Cells (10^7)	Antibody volume to add (µL)			Ab/tag ratio
				4.52 µM	1 µM	200 nM	
1		150 µL	50 µL				
2	1 µM	50 µL	50 µL	28.41			1546.6
3	316 nM	150 µL	50 µL	15.03			818.2
4	100 nM	150 µL	50 µL	4.53			246.6
5	31.6 nM	500 µL	50 µL		17.95		216.3
6	10 nM	500 µL	50 µL		5.56		67.0
7	3.16 nM	1 mL	50 µL			16.87	40.6
8	1 nM	1 mL	50 µL			5.28	12.7
9	316 pM	10 mL	50 µL			15.92	38.3
10	100 pM	10 mL	50 µL			5.03	12.1
11	31.6 pM	40 mL	50 µL			6.33	15.2
12	10 pM	40 mL	50 µL			2.00	4.8

- Prepare FACS samples:
 - Spin down the cells for 5 min at 13,400 rpm
 - Discard the supernatant and resuspend with 1 mL ice cold PBS-0.1%BSA
 - Spin down the cells for 5 min at 13,400 rpm
 - Discard the supernatant and put the pellets on ice until FACS
 - Right before FACS: resuspend with 200 µL ice cold PBS-0.1%BSA