


High Performance Liquid Chromatography (HPLC)

- The HPLC was used for substrate analytics. We were able to detect glucose, xylose, succinate and fumarate in the media.
- Sample preparation:
 - Centrifuge 1.5 ml of bacteria culture for 5 minutes at 14.000 rpm
 - Transfer supernatant to a fresh tube. Avoid to resuspend the cells
 - Add 2 µl of sodium azide (100g/ L) to each sample.
 - Store the samples at -20 °C until you analyse them by HPLC
- HPLC measurement:
 - 750 µl of each sample are transferred into HPLC vials with a septum
 - The samples are determined in duplicates together with two sets standards á three concentrations
 - The concentration of the standards depends on the maximum expected carbon source concentration in the media
- Sample run:
 - It was used the *VA 300/7.7 NUCLEOGEL SUGAR 810 H* column for sample separation
 - Eluent was 3.5 mM H₂SO₄
 - The column was heated at 65 °C and the injected sample volume was 20 µl
 - Detection was achieved by conductivity measurement



- For performing enzyme assays with purified enzymes or crude cell extracts we were using following HPLC-MC system.
- HPLC: LaChromUltra (Hitachi Europe, United Kingdom)
- Analysis: microTOF-Q hybrid quadrupole / time of flight mass spectrometer (Bruker Daltonics, Bremen, Germany)
- Ionization: Electrospray ionization (ESI)
- Column: SeQuant ZICpHILLIC column (150 by 2.1 mm, Merck KGaA, Darmstadt, Germany)
- Solutions: 10 mM Ammonium bicarbonate solution (eluent A), acetonitrile (eluent B)
- Injection volume: 2 μl
- Flow rate: 150 $\mu\text{l min}^{-1}$
- Gradient [t_{min} , % of eluent B]: t_0 : 80%, t_{30} : 10%, t_{35} : 10%, t_{40} : 80%, t_{60} : 80%
- Negative ionization mode
- Internal mass calibration: Formate (0.1 M), in 50% (vol/vol) isopropanol

