


**iGEM TU/e 2014**  
Biomedical Engineering

Eindhoven University of Technology  
Room: Ceres 0.04  
Den Dolech 2, 5612 AZ Eindhoven  
The Netherlands  
Tel. no. +31 50 247 55 59  
2014.igem.org/Team:TU\_Eindhoven

## Monitoring SPAAC with UV-VIS

General protocols from the dry lab



## Table of contents

<b>Title</b>	<b>1</b>	<b>Reagents</b>	<b>3</b>
Monitoring SPAAC with UV-VIS	2	Settings	3
	3	Monitoring SPAAC reaction	3

## 1 Reagents

- PBS with pH 7.2
- Unnatural amino acid p-azido-L-phenylalanine (pAzF) in DMSO (5%) and PBS dissolved (47.6  $\mu\text{M}$ )
- DBCO-PEG<sub>4</sub>-NHS ester solved in PBS (142.9  $\mu\text{M}$ )

## 2 Settings

Software: Spectramanager

Measurement type: Spectrum measurement

Measuring absorption spectrum	
Photometric mode	Abs
Response	Medium
UV/Vis bandwidth	1.0 nm
Scan Speed	40 nm/min
Start	340 nm
End	190 nm
Data interval	1 nm
Vertical size	Auto
Scan mode	Continuous
Accumulation/cycle	Accumulation
No. Of cycles	192
Time between cycles	300 s
Temperature	20 °C
Stirrer	250 rpm

## 3 Monitoring SPAAC reaction

- Prepare blank sample by filling cuvettes with 2.5 mL PBS
- Start blank measurement
- Prepare sample in cuvette containing 1.75 mL DBCO-PEG<sub>4</sub>-NHS ester solution and 0.75 mL pAzF solution (100 : 33.3  $\mu\text{M}$ )
- Start sample measurement