

Testing the incorporation of acetylated thymidine

Reaction Protocol:

1. Prepare MOPS buffers
2. Prepare 7x 1.5mL tubes
3. Prepare Dilutions/Aliquots
 - a. Buffer: 10X - check pH of each
 - b. Primer: 1 μ L {200pmol/ μ L} in 9 μ L water
 - c. **Ac_dTTP**: 2 μ L {10mM}, no dilution (make ten aliquots of 5 μ L)
 - d. dTTP: 0.5 μ L {100mM} in 9.5 μ L water
 - e. TdT: 5 μ L {20U/ μ L} in 5 μ L
4. Add components to all "-a" tubes (plus controls) as follows:


	ddH ₂ O	buffer	CoCl ₂	primer	monomer		enzyme	EDTA	dye
C	8 μ L	65 1 μ L		1 μ L				4 μ L	2 μ L
M-65-a	5 μ L	65 1 μ L	1 μ L	1 μ L	Ac 1 μ L		1 μ L	4 μ L	2 μ L
M-65-b					dT 1 μ L			4 μ L	2 μ L
M-79-a	5 μ L	79 1 μ L	1 μ L	1 μ L	Ac 1 μ L		1 μ L	4 μ L	2 μ L
M-79-b					dT 1 μ L			4 μ L	2 μ L
ccccc	5 μ L	65 1 μ L	1 μ L	1 μ L	dT 2 μ L		1 μ L	4 μ L	2 μ L
ccccc					dT 1 μ L			4 μ L	2 μ L

5. Get 100mM EDTA ready before adding TdT.
6. React all "-a" tubes (and controls) at 37°.
7. At each "-a" tube's 10min time point, take 5 μ L and move it into its respective "-b" tube.
8. Add additional unmodified nucleotide to respective "-b" tubes.
9. Before loading into gel, add 2 μ L 6X DNA Loading Dye

Assay Conditions (for each reaction tube):

component	moles	volume	source
TdT Rxn Buffer	-	1 μ L	10X buffer
CoCl ₂	2.5nmol	1 μ L	2.5mM CoCl ₂ solution
primer (5'-FAM-dT ₁₅)	20pmol	1 μ L	10x dilution of 1 μ L from 1 μ g/ μ L (200pmol/ μ L) tube
Ac_dTTP	10nmol	1 μ L	10mM solution
TdT	10Units	1 μ L	2x dilution of 1.5 μ L from 20U/ μ L
water (to complete rxn vol.)	-	5 μ L	
after reaction:			
EDTA (to stop reaction)	0.4 μ mol	4 μ L	100mM solution
DNA Loading Dye	-	2 μ L	6X solution with Bromophenol Blue

First attempt:



Second attempt:

