

iGEM TU/e 2014

Biomedical Engineering

Eindhoven University of Technology Room: Ceres 0.04 Den Dolech 2, 5612 AZ Eindhoven The Netherlands Tel. no. +31 50 247 55 59 2014.igem.org/Team:TU_Eindhoven

Preparative steps: culture media, agar plates, antibiotics and glycerol



Table of contents

| Title | / 1 | Preparation of Culture Media | 3 |
|--------------------|-------|--|---|
| Preparative steps | / - | • | _ |
| r reparative steps | / 1.1 | Luria-Bertani (LB) medium | 3 |
| | 1.2 | 2YT medium (for protein expression) | 3 |
| | 2 | Preparation of Agar Plates | 3 |
| | 2.1 | LB-agar | 3 |
| | 2.2 | Pouring agar plates | 3 |
| | 3 | Preparation Antibiotics | 4 |
| | 3.1 | Kanamycin | 4 |
| / | 3.2 | Chloramphenicol (for protein expression) | 4 |
| / | 4 | Preparation of Glycerol Stocks | 4 |

1 Preparation of Culture Media

1.1 Luria-Bertani (LB) medium

For 1L

- 10 g peptone
- 10 g NaCl
- 5 g yeast extract
- Add 1L H₂O

1.2 2YT medium (for protein expression)

For 1L

- 16 g peptone
- 5 g NaCl
- 10 g yeast extract
- Add 1L H₂O

2 Preparation of Agar Plates

2.1 LB-agar

For 1L

- 10 g peptone
- 10 g NaCl
- 5 g yeast extract
- 15 g bacto-agar
- Add 1L H₂O

Both media and agar need to be autoclaved (sterilization).

2.2 Pouring agar plates

- After autoclaving the LB-agar at 121 °C for 20 minutes, let the agar cool down to ~50°C.
 Make sure the agar does not start solidifying
- Add kanamycin and /or chloramphenicol to the liquid LB-agar and slowly mix.
- Pour the LB-agar in the petri dishes until the bottom is well covered. Work near the Bunsen burner flame
- Close the lid after filling the plate. Let the agar solidify for ~1 hour on the bench
- Transfer the plates to a bag, in which they should be placed upside down
- Store the plates in the fridge (4°C)

3 Preparation Antibiotics

3.1 Kanamycin

30 mg/mL stock solution in H_2O . Final concentration in cultures and LB-agar: 30 μ g/mL Filter stock solution using 0.2 μ m filter.

3.2 Chloramphenicol (for protein expression)

25 mg/mL stock solution in ethanol. Final concentration in cultures and LB-agar: 25 μ g/mL Filter stock solution using 0.2 μ m filter.

4 Preparation of Glycerol Stocks

- Prepare a 50% v/v glycerol solution in dH₂O and mix well
- Take 500 µL of desired bacterial broth to form stock.
- Add 500 μL of 50% glycerol stock to the broth and mix well
- Snap-freeze using liquid nitrogen and store at -80°C