

Gas-Chromatography (GC) – Mass Spectrometry (MS)

- Modified from Kondo et al. (2012)
- Sample preparation
 - Supernatant from cultivation (100 μl) was extracted with 900 μl of GC standard-grade acetone containing 0.1 % 2-butanol as an internal standard
 - 150 μl were used for GC-MS analysis
 - o GC-MS analysis was realized with the following settings

AS 2000 (Autosampler)

sample volume 1 µl

Trace GC 2000

column	Stabilwax® (Crossbond Carbowax - PEG)		
	(30 m, 0.25 mm internal diameter, 0.25 μm film		

thickness)

Carries Gas Helium (1 ml / min)

		rate [°C]	temp [°C]	hold time [minutes]		
Oven temperature	initial		60	3.00		
	ramp 1	10.0	100	0.00		
	ramp 2	120.0	280	01.02.00		
Injector temperature	250°C					
Split flow	20 ml /	20 ml / min				
MS transfer line	250°C					

Polaris Q (Ion trap MS)

Ion source 220°C Mass range 25-100 start time 2,5 min

- With these settings we generated a calibration line with the following concentrations for further analysis and quantification of our samples.
 - 0,5 % isobutanol
 - 0,1 % isobutanol
 - 0,05 % isobutanol
 - 0,01 % isobutanol 0,001 % isobutanol

