

iGEM TU/e 2014
Biomedical Engineering

Eindhoven University of Technology
Room: Ceres 0.04
Den Dolech 2, 5612 AZ Eindhoven
The Netherlands
Tel. no. +31 50 247 55 59
2014.igem.org/Team:TU_Eindhoven

Colony PCR

Table of contents

Title			
Colony PCR			
	1	Colony picking	3
	2	PCR	3
	3	Gel electrophoresis	3

1 Colony picking

- Using a pipette, pick a number (4-6) colonies and put them into an epp with 15 µL H₂O.

2 PCR

- Construct the PCR mixture in the following way:

Component	Quantity/mass/final concentration	Volume (µL)
DNA (from bacteria)	Pipette tip in 15 µL H ₂ O	1
2x KAPA2G mix	1x	12.5
Primer FW	0.5 µM (stock: 10 µM)	1.25
Primer RV	0.5 µM (stock: 10 µM)	1.25
H ₂ O		9
<i>Total</i>		25

Total mastermix

Component	Quantity/mass/final concentration	Volume (µL)
DNA (from bacteria)	Pipette tip in 15 µL H ₂ O	1 for each separate PCR mixture
2x KAPA2G mix	1x	62.5
Primer FW	0.5 µM (stock: 10 µM)	6.25
Primer RV	0.5 µM (stock: 10 µM)	6.25
H ₂ O		35
<i>Total</i>		125

- Run PCR in the PCR machine with the following protocol:

Step	Temperature (°C)	Time (sec)	Cycles
Denaturation	95	180 (3 min)	1
Denaturation	95	15	35
Annealing	55-60	15	
Extension	72	20 sec/kb	
Final extension	72	600 (10 min)	1
Cooling	4	Hold	1

3 Gel electrophoresis

- Once the colony PCR has been completed the samples (20 µL; first mix with 4 µL loading dye) should be loaded on to an electrophoresis gel and run.

After colony PCR, perform plasmid amplification of the vectors that are correct.