

June 23, 2014

Estrogen Receptor Intein Construction: PCR

Purpose: To synthesize the estrogen receptor intein

Table 1. PCR Protocol to Synthesize Estrogen Receptor

Reagent	Volume (μ l)	Volume (μ l)
Water	12.4	5.4
5X HF PCR Buffer	4	4
dNTPs	0.4	0.4
Primer 1 (10 uM Nterm DNA)	1	5
Primer 2 (10 uM Cterm DNA)	1	5
Phusion Taq enzyme	0.4	0.4
	Total: 20.4	~20.4

Table 2. PCR Conditions for First Round

Temperature (°C)	Time
98	2 minutes
98	10 seconds
60	15 seconds
72	1.5 minutes
72	1 minute
4	Hold

- Repeat for 35 cycles
- Take 2.5 μ l DNA from this PCR reaction for the next PCR

Table 3. PCR Protocol to Synthesize Estrogen Receptor

Reagent	Volume (µl)
Water	31
5X HF PCR Buffer	10
dNTPs	1
Primer 1 (10 uM Nterm DNA)	2.5
Primer 2 (10 uM Cterm DNA)	2.5
DNA	2.5
Phusion Taq enzyme	0.5
	Total: 50

- Same Conditions
- Take 2 - 3 µl DNA and run in 1 % agarose gel

Results: No bands present

Overnight:

- Cultures for Top10 Cells expressing fluorescent proteins

June 24, 2014

PCR to Synthesize Estrogen Responsive Intein

- Repeat protocol from previous day (tables 1, 2, and 3 from June 23, 2014)

PCR to Prepare Fluorescent Protein DNA for Submission to iGEM Registry

Table 1. Protocol for PCR Cocktail

Reagent	Volume (μ l)
Water	155
5X HF PCR Buffer	50
dNTPs	5
Primer 1 (Prefix)	12.5
Primer 2 (Suffix)	12.5
Phusion Taq enzyme	2.5
	Total: 237.5

- Add 47.5 μ l of PCR Cocktail to each tube (5 tubes)
- Add 2.5 μ l of DNA to corresponding tube
- PCR Conditions in Table 2. (June 23, 2014)

Overnight:

- Streak out MACH cells with fluorescent proteins
- Culture of MACH cells containing WT Lac

June 25, 2014

PCR for RFP with His tag

Table 1. Protocol for Single PCR Reaction

Reagent	Volume (μl)
Water	31
5X HF PCR Buffer	10
dNTPs	1
Primer 1 (10 uM Prefix)	2.5
Primer 2 (10 uM Suffix)	2.5
DNA	2.5
Phusion Taq enzyme	0.5
	Total: 50

- PCR Conditions in Table 2. (June 23, 2014)

PCR for Killer Red and Super Nova with BioBricks Prefix and Suffix

Table 2. Protocol for Single PCR Reaction

Reagent	Volume (μl)
Water	31
5X HF PCR Buffer	10
dNTPs	1
Primer 1 (10 uM Prefix)	2.5
Primer 2 (10 uM Suffix)	2.5
DNA (Killer Red or Super Nova respectively)	2.5
Phusion Taq enzyme	0.5
	Total: 50

- PCR Conditions in Table 2. (June 23, 2014)

PCR for Killer Red and Super Nova with Ribosome Binding Site

Table 2. Protocol for Single PCR Reaction

Reagent	Volume (μ l)
Water	31
5X HF PCR Buffer	10
dNTPs	1
Primer 1 (10 uM XbaPreRBS)	2.5
Primer 2 (10 uM PstSuf6His)	2.5
DNA (Killer Red or Super Nova respectively)	2.5
Phusion Taq enzyme	0.5
	Total: 50

- PCR Conditions in Table 2. (June 23, 2014)

Miniprep of WT Lac Plasmid

Protocol for GeneJet Plasmid Miniprep:

<http://www.thermoscientificbio.com/uploadedfiles/resources/k0502-product-information.pdf>

Yield = 225.5 ng/ μ l

Digestion of Wild Type Lac Vector and Fluorescent Proteins for Ligation

Table 3. Protocol for Digestion of WT Lac Plasmid

Reagent	Amount (μ l)
10X FastDigest Buffer	2
Plasmid DNA	8
Restriction Enzyme PstI	2
Restriction Enzyme SpeI	2
Sterile Water	6
Total Volume	20

- Digest at 37 °C for 30 - 45 minutes
- Gel Photo (Name: 062514b)

Table 4. Protocol for Digestion of Fluorescent Proteins, Killer Red and Super Nova Plasmids

Reagent	Amount (µl)
10X FastDigest Buffer	2
Plasmid DNA	14
Restriction Enzyme XbaI	2
Restriction Enzyme SpeI	2
Total Volume	20

- Digest at 37 °C for 30 - 45 minutes
- Gel Photo (Name: 062514c)
 - All fluorescent proteins and Killer Red (prefix/suffix) showed bands
 - None for Killer Red (RBS) or Super Nova
- Protocol: <http://www.thermoscientificbio.com/uploadedFiles/Resources/fast-digestion-dna.pdf>

Ligation of Fluorescent Proteins into WT Lac and BioBrick Vectors

Table 5. Ligation Reaction of Fluorescent Proteins

Reagent	Amount (µl)
Vector DNA (WT Lac or BioBrick Vector)	2
Insert DNA (Fluorescent Proteins)	6
Ligation Buffer	1
T4 Ligase Enzyme	1
Total Volume	10

- Ligate for 10 minutes at room temperature

Transformation of pSB3K3 and WT Lac Plasmid into MACH Cells

Purpose: To be used in cloning

- 50 µl cells for each ligation reaction
- Plate 400 µl transformants on LB + CAM plates
- Incubate at 37 °C overnight

iGEM Transformation Protocol: <http://parts.igem.org/Help:Protocols/Transformation>

June 26, 2014

Overnight:

- Cultures of the transformants (FP's in BBKR and WT Lac Backbones)