

iGEM TU/e 2014
Biomedical Engineering

Eindhoven University of Technology
Room: Ceres 0.04
Den Dolech 2, 5612 AZ Eindhoven
The Netherlands
Tel. no. +31 50 247 55 59
2014.igem.org/Team:TU_Eindhoven

Creating circular RCA template

Table of contents

Title			
Creating circular RCA template	1	Creating circular DNA	3
	2	Removal of linear DNA from mixture	3

1 Creating circular DNA

Make sure that single strand DNA template being used is 3'-hydroxylated and 5'-phosphorylated.

- Combine the following reaction components in a PCR tube:

Component	Quantity/mass/final concentration	Volume (μL)
H_2O	To 200 μL	$x \mu\text{L}$
10x CircLigase II buffer	1x	2 μL
ssDNA	0.5 μM	$y \mu\text{L}$
50 mM MnCl ₂	2.5 mM	1 μL
CircLigase II (200 U)	10 U/ μL	2 μL
<i>Total</i>		20 μL

- Run the following program in a thermo cycler

Step	Temperature (°C)	Time (min)
Ligase activity	60	60
Heat inactivation	80	10
Cooling	4	Hold

2 Removal of linear DNA from mixture

- Add 1000 U (10 μL) of Exonuclease III to the reaction mixture
- Add 100 U (5 μL) of Exonuclease I to the reaction mixture
- Program the thermo cycler as below

Step	Temperature (°C)	Time (min)
Digestion	37	60
Heat inactivation	80	20
Cooling	4	Hold