

**iGEM TU/e 2014**  
Biomedical Engineering

Eindhoven University of Technology  
Room: Ceres 0.04  
Den Dolech 2, 5612 AZ Eindhoven  
The Netherlands  
Tel. no. +31 50 247 55 59  
[2014.igem.org/Team:TU\\_Eindhoven](http://2014.igem.org/Team:TU_Eindhoven)

## Cell viability assay

This protocol is designed to check the viability of the bacteria after a click reaction with a specific DBCO functionalized molecule.

## Table of contents

Title			
Cell viability assay			
	1	<b>Culturing &amp; Protein Expression</b>	3
	2	<b>Prepare stock solutions</b>	3
	3	<b>Preparing reaction samples</b>	3
	4	<b>Preparing agar plates</b>	3
	5	<b>Diluting the samples</b>	3
	6	<b>Plating the bacteria</b>	4
	7	<b>Analyzing</b>	5

# **1      Culturing & Protein Expression**

- Culture the bacteria and induce protein expression according to the protocol Protein Expression.

# **2      Prepare stock solutions**

- 5 mM DBCO-TAMRA
- 5 mM PEG 3350
- Buffer: PBS

# **3      Preparing reaction samples**

- Prepare following tubes:

Tube name	[DBCO & PEG no DBCO]	Cells [10 <sup>9</sup> ]	DBCO volume to add ( $\mu$ L) (5 mM)	PEG (no DBCO)
DBCO 1	30 $\mu$ M	200 $\mu$ L	1.21	0
DBCO 2	30 $\mu$ M	200 $\mu$ L	1.21	0
PEG 1	30 $\mu$ M	200 $\mu$ L	0	1.21
PEG 2	30 $\mu$ M	200 $\mu$ L	0	1.21

- React the tubes for 2h at 500 RPM in shaking block at 4°C. Make sure you cover it with aluminum foil.

# **4      Preparing agar plates**

- Take 240 mL of autoclaved agar and add 240  $\mu$ L chloramphenicol (25  $\mu$ g/mL) and 240  $\mu$ L kanamycine (30  $\mu$ g/mL).
- Pour 12 agar plates.
- Put the agar plates into the 37 °C incubator 15 minutes before plating the cells.

# **5      Diluting the samples**

- Dilute the reacted samples according to the following scheme.

Tube	PBS	Sample name:	Volume sample	Concentration:
Dilution DBCO 1_1	999 µL	DBCO 1	1 µL	10^6 cells / mL
Dilution DBCO 1_2	495 µL	Dilution DBCO 1_1	5 µL	10^4 cells / mL
Dilution DBCO 1_3	180 µL	Dilution DBCO 1_2	20 µL	10^3 cells / mL
Dilution DBCO 1_4	180 µL	Dilution DBCO 1_3	20 µL	10^2 cells / mL
Dilution DBCO 2_1	999 µL	DBCO 2	1 µL	10^6 cells / mL
Dilution DBCO 2_2	495 µL	Dilution DBCO 2_1	5 µL	10^4 cells / mL
Dilution DBCO 2_3	180 µL	Dilution DBCO 2_2	20 µL	10^3 cells / mL
Dilution DBCO 2_4	180 µL	Dilution DBCO 2_3	20 µL	10^2 cells / mL
Dilution PEG 1_1	999 µL	PEG 1	1 µL	10^6 cells / mL
Dilution PEG 1_2	495 µL	Dilution PEG 1_1	5 µL	10^4 cells / mL
Dilution PEG 1_3	180 µL	Dilution PEG 1_2	20 µL	10^3 cells / mL
Dilution PEG 1_4	180 µL	Dilution PEG 1_3	20 µL	10^2 cells / mL
Dilution PEG 2_1	999 µL	PEG 2	1 µL	10^6 cells / mL
Dilution PEG 2_2	495 µL	Dilution PEG 1_1	5 µL	10^4 cells / mL
Dilution PEG 2_3	180 µL	Dilution PEG 1_2	20 µL	10^3 cells / mL
Dilution PEG 2_4	180 µL	Dilution PEG 1_3	20 µL	10^2 cells / mL

## 6 Plating the bacteria

- Plate the following samples on the prepared agar plates.

Tube	Volume for plating	Number of bacteria
Dilution DBCO 1_2	100 µL	1000
Dilution DBCO 1_3	100 µL	100
Dilution DBCO 1_4	100 µL	10
Dilution DBCO 2_2	100 µL	1000
Dilution DBCO 2_3	100 µL	100
Dilution DBCO 2_4	100 µL	10
Dilution PEG 1_2	100 µL	1000
Dilution PEG 1_3	100 µL	100
Dilution PEG 1_4	100 µL	10
Dilution PEG 2_2	100 µL	1000
Dilution PEG 2_3	100 µL	100
Dilution PEG 2_4	100 µL	10

- Incubate the plates for ~16 hours on 37 °C.

## **7      Analyzing**

- The following day take the plates out of the incubator.
- Count the number of colonies and compare them with the original amount of bacteria plated.