

iGEM TU/e 2014

Biomedical Engineering

Eindhoven University of Technology Room: Ceres 0.04 Den Dolech 2, 5612 AZ Eindhoven The Netherlands Tel. no. +31 50 247 55 59 2014.igem.org/Team:TU_Eindhoven

Site Directed Mutagenesis



Table of contents

Title Site Directed Mutagenesis	1	PCR	3
	2	Dpn I Digestion of the Amplification Products	3
,	7	Transformation of YI 10-Gold Ultracompetent Cells	3

1 PCR

Component	Quantity/mass/final concentration	Volume (µL)
dH₂O	Fill up to 25 μL	
10× QuikChange Lightning	1x	2.5
Multi reaction buffer		
Plasmid DNA	100 ng (stock ng/μL)	
Forward primer	100 ng (10 μM stock)	
dNTP mix	1x	1
QuikChange Lightning Multi	1x	1
enzyme blend		
Total		25

Segment	Cycles	Temperature	Time
1	1	95°C	2 min
2	30	95°C	20 sec
		55°C	30
		65°C	30 seconds/kb of
			plasmid length -> 3 min
3	1	65°C	5 min

2 Dpn I Digestion of the Amplification Products

- Add 1 µL of Dpn I restriction enzyme directly to each amplification reaction.
- Gently and thoroughly mix each reaction mixture by pipetting the solution up and down several times. Spin down the reaction mixtures in a microcentrifuge for 1 minute, then immediately incubate each reaction at 37°C for 5 minutes to digest the parental (nonmutated) ds-DNA.

3 Transformation of XL10-Gold Ultracompetent Cells

- Gently thaw the XL10-Gold ultracompetent cells on ice. For each mutagenesis reaction, aliquot 45 μL of the ultracompetent cells to an Eppendorf tube
- Add 2 μ L of the β -ME mix provided with the kit to the 45 μ L of cells. Using an alternative source of β -ME may reduce transformation efficiency.
- Swirl the contents of the tube gently. Incubate the cells on ice for 10 minutes, swirling gently every 2 minutes
- Transfer 4 μ L of the Dpn I-treated DNA from each mutagenesis reaction to a separate aliquot of the ultracompetent cells

- Swirl the transformation reactions gently to mix, then incubate the reactions on ice for 30 minutes
- Heat-pulse the tubes in a 42°C water bath for 30 seconds exactly
- Incubate the tubes on ice for 2 minutes
- Add 0.5 ml SOC medium to each tube and incubate the tubes at 37°C for 1 hour with shaking at 225–250 rpm
- Plate 100 μL of each transformation reaction onto agar plates
- Incubate the transformation plates at 37°C for >16 hours