## **RPU Measurement Method**

- 1. Grow cell culture overnight (Incubate 37°C and shake for 15 hours) in falcon tubes with M9 minimal medium
- 2. Prepare M9 medium with desired volume in a 96 Deep Well plate. The positive control for the experiment is BBa\_I20260, which contains constitutive promoter J23100, and GFP generator.
- 3. Take out 20-30  $\mu$ l of overnight cell culture and mix it with M9 medium in the 96 Deep Well plate
- 4. Incubate in 37°C and shake for 3-4 hours.
- 5. Take out 200ul of cells from the 96 well plates, and put it on a micro test plate 96 well flat bottom.
- 6. Measure the OD<sub>600</sub> and green fluorescence protein expression of the cells plate reader.
- 7. Measure every 30 minutes span, for 2 3 hours. Make sure that the cells are in the mid log phase ( $OD_{600} = 0.3-0.5$ ). In between measurements, keep incubating the cells in 37°C while shaking.
- 8. Calculate the RPU (Relative Promoter Unit) based on the paper "Measuring the activity of BioBrick promoters using an in vivo reference standard" by Kelly et al.

$$RPU = \frac{S_{cell,\phi}^{SS}}{S_{cell,J23101}^{SS}} = \frac{(dF_{\phi}/dt)/ABS_{\phi}}{(dF_{J23101}/dt)/ABS_{J23101}}$$