



# Cohesive-end Assembly Cloning

This protocol was used to perform all the Biobricks assemblies.

Insert calculation:

$$ng\ of\ insert = \frac{ng\ of\ vector \times kb\ size\ of\ insert}{kb\ size\ of\ vector} \times insert:vector\ molar\ ratio$$

## Materials

- Insert 1 previously digested with *EcoRI* and *SpeI*
- Insert 2 previously digested with *XbaI* and *PstI*
- pSB1C3 previously digested with *EcoRI* and *PstI*
- 10X T4 DNA Ligase Buffer
- T4 DNA Ligase (5U/μL) (Thermo Scientific)
- Water, nuclease free
- 0.5 mL microtube

## Apparatus

- Incubator

## Method

1. After thawing and vortexing the solutions (except the enzyme solution), prepare the following mixture in a 0.5 sterile microtube. Keep the microtube on ice.

Insert 1 digest with <i>EcoRI</i> and <i>SpeI</i> (vector:insert)	1:3 to 1:6
Insert 2 digest with <i>XbaI</i> and <i>PstI</i> (vector:insert)	1:3 to 1:6
pSB1C3 digested with <i>EcoRI</i> and <i>PstI</i>	50 ng
10X T4 DNA Ligase Buffer	2 μL
T4 DNA Ligase (5 U/μl)	1 μL (5U)
Water, nuclease free	to 20 μL

2. Spin the mixture.
3. Incubate overnight at 4°C.
4. Transform the ligation in *E. coli* DH5α. (see the Transformation in *E. coli* Protocol) [\[http://2014.igem.org/wiki/images/a/a5/Transformation\\_in\\_Escherichia\\_coli\\_DH5.pdf\]](http://2014.igem.org/wiki/images/a/a5/Transformation_in_Escherichia_coli_DH5.pdf)
5. Plate out the suspension on a LB agar plate containing chloramphenicol (25 μg/mL) and incubate at 37° for 12 hours.
6. Select the white colonies.