

NEB BioBrick Assembly Kit

- We used the BioBrick Assembly Kit to assemble an upstream part with a downstream part into destination plasmid.
- Digestion Protocoll:
 - Digest upstream part with **EcoRI-HF** and **SpeI**
 - Digest downstream part with **XbaI** and **PstI**
 - Digest destination plasmid with **EcoRI-HF** and **PstI**
 - 500 ng part DNA
 - 1 µl of each enzyme
 - 5 µl 10x NEBuffer 2.1
 - to 50 µl H₂O
- Incubate all three restriction digest reactions at 37 °C for 10 minutes and then heat inactivate at 80 °C for 20 minutes.
- Dephosphorylation Protocoll
 - Add 1 µl of AP (Antarctic phosphatase) and 5 µl of 10x AP reaction buffer to digested destination plasmid, incubate for 1 h at 37 °C.
- Ligation Protocoll
 - Ligate the upstream and downstream parts into the digested destination plasmid.
 - 2 µl of each part
 - 2 µl 10x T4 DNA Ligase Buffer
 - 1 µl T4 DNA Ligase
 - 11 µl H₂O
 - Incubate at room temperature for 10 minutes and then heat inactivate at 80 °C for 20 minutes.
- Transform 1-2 µl of the ligation product into 50 µl of competent E. coli cells.

