

Antibiotic stocks protocol

Materials

- Chloramphenicol; stock 1000x 34 mg/ml (323.132 g/mol)
- Kanamycin; stock 1000x 10 mg/ml (484.499 g/mol)
- Ampicillin; stock 1000x 100 mg/ml (349.41 g/mol)
- 99,9% pure EtOH
- dH₂O

Method

1. According to desired amount of antibiotics, calculate the amount in grams that has to be weighed
2. For example, make 11.5mL of Chloramphenicol: $34 \text{ mg/ml} * 11.5\text{mL} = 391\text{mg} / 1000 = 0.391 \text{ gram}$
3. Chloramphenicol has to be diluted in EtOH
4. Dilute Ampicillin and Kanamycin in dH₂O
5. Divide in aliquots of 500 μ L
6. Save in freezer -20°C