


Gibson Assembly

- Modified from [Gibson et al. \(2009\)](#)
- This assembly method is an isothermal, single-reaction method for assembling multiple overlapping DNA molecules. By coordinating the activity of a 5' exonuclease, a DNA polymerase and a DNA ligase two adjacent DNA fragments with complementary terminal sequence overlaps can be joined into a covalently sealed molecule, without the use of any restriction endonuclease.
- Preparation of DNA molecules for *in vitro* recombination
- Generate the complementary sequence overlaps by PCR using the [Phusion DNA-polymerase](#). If necessary add 5 M Betain in the reaction mix by reducing the amount of H₂O to decrease the number of false PCR products.
- Identify the PCR products of interest by gel electrophoresis with known DNA standards.
- Extract the PCR products from the gel by cutting out the DNA fragments and clean them up by using a commercial [PCR clean-up kit](#).
- *in vitro* recombination
- Assembly mixture:
 - 320 µl 5x isothermal reaction buffer
 - 0.64 µl of 10 U ml⁻¹ T5 exonuclease (for DNA molecules overlapping by greater than 150 bp add 3.2 µl of 10 U ml⁻¹ T5 exonuclease)
 - 20 µl of 2 U ml⁻¹ [Phusion DNA polymerase](#)
 - 160 µl of 40 U ml⁻¹ taq DNA ligase
 - add ddH₂O water up to a final volume of 1.2 ml
- aliquote 15 µl of the reagent-enzyme mix and store it at -20 °C
- Thaw 15 µl assembly mixture aliquot and keep it on ice until use.
- Add 5 µl of the purified DNA molecules in equimolar amounts (between 10 and 100 ng of each DNA fragment).
- Incubate the resulting mixture at 50 °C for 15 to 60 min, with 60 min being optimal.
- Transformation ([via heat shock](#) or [via electroporation](#)) without cleaning up the assembly product.



60