

Protocol for a Routine *Taq* PCR

Introduction

All components should be mixed and spun down prior to pipetting. These recommendations serve as a starting point; in order to maximize amplification the reaction conditions may require optimization (see [Taq DNA Polymerase Guidelines for PCR Optimization](#) protocol).

Protocol

1. Prepare the following 50 μ l reaction in a 0.5 ml PCR tube on ice:

** Due to the difficulties in pipetting small volumes of enzyme, Taq DNA Polymerase can be diluted in Diluent F ([NEB #B8006S](#)) or 1X reaction buffer. For example, 1 μ l of Taq DNA Polymerase is mixed with 4 μ l of diluent and 1 μ l of that mixture is added to the reaction. Enzyme diluted in Diluent F can be stored at -20°C for future use.*

Component	25 μ l reaction	50 μ l reaction	Final Concentration
10X ThermoPol or Standard Taq Reaction Buffer	2.5 μ l	5 μ l	1X
10 mM dNTPs	0.5 μ l	1 μ l	200 μ M
10 μ M Forward Primer	0.5 μ l	1 μ l	0.2 μ M (0.05–1 μ M)
10 μ M Reverse Primer	0.5 μ l	1 μ l	0.2 μ M (0.05–1 μ M)
Template DNA	variable	variable	<1,000 ng
Taq DNA Polymerase*	0.125 μ l	0.25 μ l	1.25 units/50 μ l PCR
Nuclease-free water	to 25 μ l	to 50 μ l	

2. Gently mix the reaction and spin down in microcentrifuge.
If the thermocycler does not have a heated cover, add one drop of mineral oil to the reaction tube to prevent evaporation.
3. Cycling Conditions for a Routine PCR:

STEP	TEMP	TIME
Initial Denaturation	95°C	30 seconds
30 Cycles	95°C 45–68°C 68°C	15–30 seconds 15–60 seconds 1 minute/kb
Final Extension	68°C	5 minutes
Hold	4–10°C	