

Testing the effect of pH on acetyl hydrolysis

Reaction Protocol:

1. Prepare six 1.5mL tubes and label them: C, 6.5, 7.0, 7.5, 7.9, NEB
2. Prepare four MOPS 10X buffers (10mL), one at each pH to be tested.
3. Prepare Dilutions/Aliquots
 - a. Primer: 0.5 μ L {200pmol/ μ L} in 4.5 μ L water
 - b. **Ac_dTTP**: 4 μ L {10mM}, no dilution
4. Add the following to all tubes (control only gets primer, buffer, and water):
 - a. 1 μ L Rxn Buffer (from 10X of respective pH)
 - i. add instead 1 μ L of NEB's TdT buffer to the NEB tube
 - b. 1 μ L CoCl₂ (from 2.5mM stock)
 - c. 1 μ L Primer dilution
 - d. 1 μ L **Ac_dTTP** aliquot
 - e. 5 μ L water
5. Warm up raw ingredients in tubes before adding TdT (optional)
6. Get 100mM EDTA ready (for stopping the reaction)
7. Prepare TdT aliquot (at fridge, then move to incubator room for Step 8)
 - a. 1.5 μ L from stock (20U/ μ L) into 1.5 μ L sterile H₂O
8. Add 1 μ L from TdT dilution to the three rxn tubes and start timing immediately
 - a. alternatively, use 0.5 μ L each from stock TdT
9. Leave all four tubes in 37° incubator for exactly 10 minutes
10. At each tube's timepoint, remove from incubator and add 4 μ L EDTA
11. Before loading into gel, add 2 μ L 6X DNA Loading Dye

Assay Conditions (for each reaction tube):

component	moles	volume	source
TdT Rxn Buffer	-	1 μ L	10X buffer
CoCl ₂	2.5nmol	1 μ L	2.5mM CoCl ₂ solution
primer (5'-FAM-dT ₁₅)	5pmol	1 μ L	20x dilution of 0.5 μ L from 1 μ g/ μ L (200pmol/ μ L) tube
Ac_dTTP	10nmol	1 μ L	10mM solution
TdT	10Units	1 μ L	2x dilution of 1.5 μ L from 20U/ μ L
water (to complete rxn vol.)	-	5 μ L	
after reaction:			
EDTA (to stop reaction)	0.4 μ mol	4 μ L	100mM solution
DNA Loading Dye	-	2 μ L	6X solution with Bromophenol Blue

