

PCR Reaction of Hsp60 Promoter June 25, 2014

Purpose: Amplification of Hsp60 Promoter for insertion into biobrick

Reaction Reagents:

1. mCherry Fushion hsp60 promoter template
2. Primers
  - a. Hsp60BBfwd
  - b. Hsp60rev comps
  - c. RBS\_BBrev comps
3. PCR Buffer
4. Double distilled water
5. DNA Polymerase
6. 1kb DNA Ladder

Primer Numbering and Concentration:

1. Hsp60BBfwd: 29.1 nmol/100ul
2. Hsp60rev comps: 26.2 nmol/100ul
3. RBS\_BBrev comps 27.7 nmol/100ul

Final Concentration 100 uM of DNA

Reagent Table 1:

Reagent	3x (ul)
Clear Color Buffer	12
dNTP's	1.2
Primers: 1 + 3	6
Template DNA	3
DNA Polymerase	0.6
ddH <sub>2</sub> O	37.2
Total Volume	60

Reagent Table 2:

Reagent	3x (ul)
Clear Color Buffer	12
dNTP's	1.2
Primers: 1 + 2	6
Template DNA	3
DNA Polymerase	0.6
ddH <sub>2</sub> O	37.2
Total Volume	60

Loading Buffer Concentration Table:


x= total loading volume

$$x = \text{total loading volume}$$

$$x - \frac{1}{5}x = 5$$

Reaction loaded (ul)	Loading Buffer (ul)
4	1
5	1.25
8	2
10	2.5

Result of PCR with Hsp60 June 25th, 2014  
 Friday, June 27, 2014  
 Result of PCR with Hsp60



Gel Analysis:  
 Loading Table

Lane	Sample	Amount Loaded (ul)
1	1kb DNA Ladder	2
2	(1) Hsp60BBfwd (3) RBS_BBrev	6.25
3	Control, no primers	6.25
4	Ladder	6.25
5	(1) Hsp60BBfwd (2) Hsp60BBrev	6.25

6	Control, no primers	6.25
---	---------------------	------

It seems that too much template DNA was used hence the bands did not travel far

## Liquid Culture with DAM- E. coli in LB

Thursday, June 26, 2014

1. Set up 3 liquid cultures with 4 ml of LB each and 4 ul of Ampicillin (1:1000 Concentration)
2. Incubated in 37°C shaker overnight

## Mini-Prep DAM- Negative Liquid Culture

Friday, June 27, 2014

### Buffer Preparation:

Resuspension Solution: add provided RNase A1 store at 4°C for up to 6 months

Wash Solution: Add 35 ml of 96% EtOH to 20 ml wash solution. Store at room temperature (RT); we have 100% EtOH

$$C_1 V_1 = C_2 V_2$$

$$(1) V_1 = 0.96(35ml)$$

$$V_1 = 33.6 \text{ ml of } 100\% \text{ EtOH} + 1.4 \text{ ml } ddH_2O$$

1. Pipet 1.5 mL of culture into microcentrifuge tube and spin for 2 minutes
2. Resuspend pellet in 250 mL of Resuspension solution
3. Add 250 mL of Lysis solution and mix by inverting tube
4. Add 360 mL of neutralization and mix immediately by inversion 4-6 times
5. Centrifuge for 5 minutes to pellet cell debris
6. Transfer supernatant to spin column
7. Centrifuge for 1 minute and discard flow through
8. Add 500 mL of wash solution and centrifuge for 30-60 seconds and discard flow through
9. Repeat wash using 500 mL of wash solution
10. Transfer spin column to fresh 1.5 mL tubes and add 50 mL of elution buffer to center of column
11. Incubate at RT for 2 minutes centrifuge for 2 minutes
12. Store flow through (purified plasmid DNA) at -20 degree Celsius.

Result:

Tube	Concentration (ng/uL)
1	29.1
2	29.5
3	18.4
4	25.9
5	18.5
6	16.8
7	23.9
8	15.6

## PCR of Hsp60 with and without RBS

Friday, June 27, 2014

Making a 1:13 dilution of our template DNA which is 140 ng/uL

Components	Amount (ul)
deionized water	12
Template DNA	1

Reagent Table 1 (w/o RBS):

Reagent	3x (ul)
Clear Color Buffer	12
dNTP's	1.2
Primers: 1 + 3	3
Template DNA	3
DNA Polymerase	0.6
ddH <sub>2</sub> O	37.2
Total Volume	57

Reagent Table 2 (w/ RBS):

Reagent	3x (ul)
Clear Color Buffer	12
dNTP's	1.2
Primers: 1 + 2	3
Template DNA	3
DNA Polymerase	0.6
ddH <sub>2</sub> O	37.2
Total Volume	57

Primer Numbering and Concentration:

1. Hsp60BBfwd: 29.1 nmol/100ul
2. Hsp60rev comps: 26.2 nmol/100ul
3. RBS\_BBrev comps 27.7 nmol/100ul

Final Concentration 100 uM of DNA

1.25 uL of loading dye + 5 uL PCR reaction  
2ul 1 kb DNA ladder

Started gel at 4:26 pm  
Stopped gel at 4:55 pm

Results:

There were no PCR Products

## PCR with Hsp60

Monday, June 30, 2014

Reagent Table 1 (w/o RBS):

Reagent	2.5x (ul)
Clear Color Buffer	25
dNTP's	2.5
Primers: 1 + 3	12 (6 each)
Template DNA	6
Phusion	1.25
ddH <sub>2</sub> O	77.5
Total Volume	125

Reagent Table 2 (w/ RBS):

Reagent	2.5x (ul)
Clear Color Buffer	12
dNTP's	1.2
Primers: 1 + 2	12 (6 each)
Template DNA	6
Phusion	1.25
ddH <sub>2</sub> O	77.5
Total Volume	125

### PCR Tube Labeling:

1 = fwd primer + rev primer RBS

2 = fwd primer + rev primer Hsp60 (without RBS)

### New PCR Condition:

Temperature ( °C)	Time
98	30 secs
98	10 secs
64	15 secs
72	15 secs
72	5 mins
4	hold

Some housekeeping tasks

Monday, June 30, 2014

0.8% Agarose Gel:

25 mL 1x TBE buffer

0.2 g Agarose

1uL EtBr

100 mL LB and autoclaved:

1g tryptone

0.5g yeast extract

1g NaCl

Loaded 1 kb of volume 1 ul and then 6.25 ul of samples (w/RBS in lane 2 and w/o RBS lane 3)

Gel started at 3:30 pm

Gel stopped at 4:40 pm

Made 10, 5 mL LB+AMP, E. coli dam- pBRES36a culture at 6 pm



Mini-Prep

Tuesday, July 1, 2014

Mini-prepped 50 mL of E. coli dam- pBRES36a in LB+AMP

Obtained:

450 ul of 17.6 ng/ul of DNA

50 ul of 14.6 ng/ul of DNA

## Primer Walking and PCR Purification

Wednesday, July 2, 2014

### Primer Walking:

Belle submitted the plasmid (pBRES36a) to GeneWiz for Primer Walking

### PCR Purification:

Purified 45 mL of both RBS and no RBS PCR products using new SV Wizard Kit from Promega

Added 75 mL of 95% EtOH to membrane wash solution

### Final Concentration:

Tube Labeling	PCR Product	Concentration (ng/uL)
1	RBS	17.0
2	No RBS	20.9

## Restriction Digest of PCR Product with XbaI and PstI

Thursday, July 3, 2014

### Result

1. 8.6 ug/mL for RBS PCR product
2. 11.5 ug/mL for no RBS PCR product

## Mini-Prep

Tuesday, July 8, 2014

Mini-prepped 5x5 mL of E. coli dam- pBRES36a in LB+amp

Ran out of lysis solution ,unable to mini prep the remaining 25 mL

Yield:

5x50 mL = 250 mL

Concentration: 42 ug/mL, equivalent to 10.5 ng plasmid

Testing old Tubes and Preparing LB Agar (100mL)  
Wednesday, July 9, 2014

Testing Old Mini-prep Columns and Tubes:

Use 20 uL of 3.0 ug/mL DNA

Result

New tubes: 52 ug/mL

Old Tubes: 20 ug/mL

Conclusion: old tubes can no longer be used purchased new tubes

100 mL LB Agar and autoclaved:

1g tryptone

0.5g yeast extract

1g NaCl

1.5 g Agar (1.5%)

100 mL Water (ddH<sub>2</sub>O)

50 mL aliquots

## Ligating Hsp60 Inserts Into Vector Backbone

Monday, July 14, 2014

### 1. Digestion

Materials:

Reagent	Amount (ul)
DNA	40
Buffer	5
XbaI restriction enzyme	1.5
PstI restriction enzyme	1.5
ddH2O	2
Total	50

Procedures:

1. Mixed materials together in microcentrifuge tubes
2. Incubate 30 minutes at 37 degrees Celsius

### 2. PCR Purification (using Promega Wizard SV gel clean up kit)

### 3. Ligations

- a. Ligation 1: RBS and Vector
- b. Ligation 2: no RBS and Vector

Materials for Ligation 1:

Reagent	Amount (ul)
Ligase	1
Buffer T4	2
30 ng of RBS insert	4
60 ng of vector plasmid	3.5
ddH2O (Volume to 20)	9.5

	(ul)
Ligase	1
Buffer T4	2
30 ng of no RBS insert	3
60 ng of vector plasmid	3.5
ddH2O (Volume to 20)	10.5

Materials for Ligation 2

Reagent	Amount
---------	--------

Procedures:

1. 15 minutes at room temperature
2. Put on ice until transformation

### 4. Transformation (see transformation protocol)

Housekeeping

Monday, July 14, 2014

30 mL LB Agar and autoclaved:

0.3 g tryptone

0.18 g yeast extract

0.3 g NaCl

0.45 g Agar (1.5%)

30 mL Water (ddH<sub>2</sub>O)

Chloramphenicol (CAM)

10 mg in 1 mL of 95% EtOH

5ug/mL workign concentration is needed

$$C_1 V_1 = C_2 V_2$$
$$\left(10 \frac{mg}{mL}\right) V_1 = 5.0 \frac{ug}{mL} (30 mL)$$
$$V_1 = 15 mL$$

Used 35 ug/mL final

Prep 50 mL LB agar +CAM:

0.5 g tryptone

0.25 g yeast extract

0.5 g NaCl

0.75 g Agar (1.5%)

48 mL Water (ddH<sub>2</sub>O)

25 uL of 10 mg/mL CAM

Competent cells growing

Housekeeping and Ligation  
Wednesday, July 16, 2014

Prep 75 m

0.75 g tryptone

0.375 g yeast extract

0.5 g NaCl

1.125 g Agar (1.5%)

50 mL Water (ddH<sub>2</sub>O)

Ligations

Ligation 1: RBS and Vector

Ligation 2: no RBS and Vector

Materials for Ligation 1:

Reagent	Amount (ul)
Ligase	1
Buffer T4	2
30 ng of RBS insert	4
60 ng of vector plasmid	3.5
ddH <sub>2</sub> O (Volume to 20)	9.5

	(ul)
Ligase	1
Buffer T4	2
30 ng of no RBS insert	3
60 ng of vector plasmid	3.5
ddH <sub>2</sub> O (Volume to 20)	10.5

Materials for Ligation 2

Reagent	Amount
---------	--------

Transformation:

50 mL comp cells + 5 uL plasmid (shipment plasmid + RBS/no RBS)

On ice for 2 minutes

Heat shock 42 degree Celsius for 45 seconds

Ice 2 minutes

550 mL SOC media into mixture

Shake at 37 degree Celsius for 1 hour (recovery)

Plate onto CAM plate + incubate at 4 pm



Plan + CAM plate  
Friday, July 18, 2014

Plan of Attack with pBRES36a

1. Digest with individual restriction enzymes and a negative control
  - a. No enzyme
  - b. XbaI
  - c. PstI
  - d. SpeI
  - e. EcoRI
2. Visualize on Gel
3. Double Digest with Pair from Step 1
4. Visualize on gel
5. Map the sequence

CAM plates made with 35 ug/mL

Transformation with:

1. Ligation product (2 of them)
2. Plasmid that's known to express CAM (for control)

7/20/14

2<sup>nd</sup> Double Digestion: For the purpose of confirming that there was no mix up in the previous digestion.

Components	Amount (uL)
10x Fast Digest Buffer	2
pBRES36a (42ng/uL)	16
SpeI	1
PstI	1
Total	20



Lane # (from left)	Sample
1	Ladder (not visible)
2	No enzyme
3	PstI and SpeI
4	No enzyme
5	PstI and SpeI

Housekeeping

Wednesday, July 23, 2014

Cam plates: 35 ug/mL, dilute 1:1000

Used 75 mL, used 75 ul

For 4 plates

Transformation with RBS:

50 mL E. coli comp cells, 5 uL RBS ligation

50 mL E. coli comp cells, 5 uL no RBS ligation

Transformed and recovered for 1 hour

Plated 50 uL of transformed cells onto cam plates

# Restriction Digest of pBRES36a Plasmid

Monday, July 28, 2014

## 1. Digestion

**Purpose:** cut with restriction enzymes to provide a rough map for the pBRES36a plasmid and confirm the plasmid's identity

### Materials:

Reagents	Amount (ul)
pBRES36a plasmid	17
Buffer	2
Restriction Enzyme	1
Total Volume	20

### Procedures:

1. Mix reagents together and spin down gently.

## 2. Gel Electrophoresis

Run 1: Gel Set Up (7/18/14)

lane	Reaction/Reagent	Amount loaded (ul)
1	1 kb Marker	1
2	Plasmid + XbaI	5
3	Plasmid + PstI	5
4	Empty Lane	0
5	Plasmid + SpeI	5
6	Plasmid +EcoRI	5
7	Plasmid only	5

Run 1 Condition:

1. 2.0 hours
2. 80 Volts
3. 0.8% Agarose Gel
4. 1x TBE Buffer

2. Incubated for 10 minutes
3. Pipet up and down gently to mix samples
4. Load each sample onto the agarose gel
5. Run electrophoresis according to the conditions specified.

### Enzyme Used:

1. XbaI
2. PstI
3. SpeI
4. EcoRI

Run 2: Gel Set Up (7/28/14)

lane	Reaction/Reagent	Amount loaded (ul)
1	1 kb Marker	1
2	Plasmid + XbaI	5
3	Plasmid + PstI	5
4	Plasmid + SpeI	5
5	Plasmid +EcoRI	5
6	Plasmid +SpeI and EcoRI	5
7	Hsp60 + XbaI and PstI	5
8	mCherry + XbaI and PstI	5


Run 2 Condition:

1. 1.5 hours
2. 80 Volts
3. 0.8% Agarose Gel
4. 1x TBE Buffer

July 28, 2014 con't:

**Results**

Run 1 (7/18/14) Left and Run 2 (7/28/14) Right



Ligation, Transformation, Selection  
Thursday, July 24, 2014

Ligations

Ligation 1: RBS and Vector

Ligation 2: no RBS and Vector

Materials for Ligation 1:

Reagent	Amount (ul)
Ligase	1
Buffer T4	2
30 ng of RBS insert	4
60 ng of vector plasmid	3.5
ddH2O (Volume to 20)	9.5

Materials for Ligation 2

Reagent	Amount
---------	--------

Incubate at RT for 15 minutes

Transformation

1. 50 mL E. coli comp cells
2. 6 uL plasmid
3. 2 min on ice
4. 45s heat shock at 42 degrees Celsius
5. 2 min on ice
6. 950 mL SOC growth media
7. 1 hour recovery at 37 degrees Celsius

Selection:

plate 60 mL of transformed cells onto selection plate (LB + CAM plates)

	(ul)
Ligase	1
Buffer T4	2
30 ng of no RBS insert	3
60 ng of vector plasmid	3.5
ddH2O (Volume to 20)	10.5

8/1/14

1. Take out overnight liquid culture around 11:00 am
2. Mini prep iGEM plasmid parts

Part	Amount (ng/uL)
Desaturase	42.6
Cathelicidin	12.7
mCherry	12.7

8/4/14 Lab

**Make 1% Gel (30 mL):**

0.3 g Agarose  
30 mL 1x TBE  
1 uL EtBr

**Restriction Digests:**

hsp60 PCR products:

<b>RBS:</b>	<b>no RBS:</b>
6 uL ddH <sub>2</sub> O	6 uL ddH <sub>2</sub> O
2 uL Buffer	2 uL Buffer
1 uL XbaI	1 uL XbaI
1 uL PstI	1 uL PstI
10 uL PCR product	10 uL PCR product

iGEM Parts:

<b>mCherry:</b>	<b>Cathelicidin:</b>	<b>deSaturase:</b>
6 uL ddH <sub>2</sub> O	6 uL ddH <sub>2</sub> O	6 uL ddH <sub>2</sub> O
2 uL FD Green Buffer	2 uL FD Green Buffer	2 uL FD Green Buffer
1 uL XbaI	1 uL XbaI	1 uL XbaI
1 uL PstI	1 uL PstI	1 uL PstI
10 uL mCherry	10 uL Cathelicidin	10 uL deSaturase

**PCR Purification:**

Run PCR purification of hsp60 digests (RBS and no RBS)

- 1: RBS concentration-
- 2: no RBS concentration –

**Gel Electrophoresis:**

Ladder – mcherry – cath – desat - empty - cath – mcherry – desat

Failed. Will upload pic later

**Ligation:**

Ligate hsp60 RBS and no RBS into pSB1C3 backbones

<b>RBS:</b>	<b>no RBS:</b>
9.5 uL ddH <sub>2</sub> O	10.5 uL ddH <sub>2</sub> O
3.5 uL backbone	3.5 uL backbone
4 uL insert	3 uL insert
2 uL buffer	2 uL buffer



1 uL ligase

1 uL ligase

**Transformation:**

Transform RBS and no RBS ligations into *E. coli* competent cells

5 uL plasmid

100 uL competent cells

Ice 2 minutes

Heat shock (42C) 45s

Ice 2 minutes

Add 950 uL outgrowth media

Recover with shaking (37C) for 1 hour

Plate onto selection plate (chloramphenicol)

**Restriction Digest:**

Cut mCherry prep, Cath. Prep., deSat. Prep, 1 RBS, 1 no RBS, 1 mCherry with XbaI and PstI

6 uL ddH<sub>2</sub>O

2 uL FD Green Buffer


1 uL XbaI

1 uL PstI

10 uL plasmid DNA

**Gel Electrophoresis:**

1 kb+ ladder – mCherry prep – Cath. Prep – deSat prep – 1 RBS – 1 no RBS – 1 mCherry -



8/5/15 Lab

**Make 1% Gel (35 mL)**

0.35 g Agarose  
35 mL 1x TBE  
1 uL EtBR

**Restriction Digests:**

Master Mix (for 6 reactions): x6.5 =

6 uL ddH <sub>2</sub> O	39 uL
2 uL Green Buffer	13 uL
1 uL XbaI	6.5 uL
1 uL PstI	6.5 uL

Combine 10 uL of Master Mix with 10 uL of each plasmid:

mCherry Prep, Cath. Prep, deSat prep., 1 RBS, 1 no RBS, 1 mCherry


Incubate 45 mins at 37C

**Gel Electrophoresis:**

Run 1 hour 30 mins

Lanes:

1 kb+ ladder – mCherry mini – Desat – Cath – 1 RBS – 1 no RBS – 1 mCherry – 1 kb+ ladder




08/08/14

Used the culture collected from 8/1/14 mini-prep for the subsequent restriction digest

1. Restriction digest (45 minutes)
2. Run agarose gel

	Plasmid	Expected Band (kb)
1	Desaturase	~1
2		
3	Cathelicidin	~0.1
4		
5	RBS	~0.3
6		
7	No RBS	~0.3
8		



Digestion was incomplete and need to repeat experiment

Experiment was repeated on 8/9/14 and digestion was successful. The desire bands were excised and purified. Unfortunately the image file was lost due to a technical problem prior to routine file backup. However, further experiments continued to use parts isolated on this date and bands of correct size were present through all of these steps.

## Lab Work 8/12/14:

### **Gel Purification:**

Gel Fragment	Gel + Tube Mass (g)	Tube Mass (g)	Gel Mass (g)	Membrane Binding Solution Added (uL)	DNA Concentration (ng/uL)
1 Cath. Part	1.133	1.050	0.083	83	2.9
2 Cath. Part	1.146	1.050	0.096	96	3.9
3 Desat. Part	1.173	1.050	0.123	123	3.1
4 Desat. Part	1.192	1.050	0.142	142	3.1
1 Backbone	1.160	1.050	0.110	110	3.2
2 Backbone	1.152	1.050	0.102	102	3.2
3 Backbone	1.151	1.050	0.101	101	3.2
4 Backbone	1.200	1.050	0.150	150	4.3
5 Backbone	1.201	1.050	0.151	151	6.6
<del>6 Backbone</del>	<del>1.058</del>	<del>1.050</del>	<del>0.008</del>	<del>8</del>	N/A
7 Backbone	1.243	1.050	0.193	193	5.7
8 Backbone	1.159	1.050	0.109	109	5.7
Backbone	1.300	1.050	0.250	250	4.4

\*Add 10 uL of Membrane Binding Solution per 10 mg of Gel slice

- Vortex Gel and solution
- Incubate at 55C for 10 mins to melt gel
- Spin Down
- Gel Purify

### **Miniprep:**

Miniprepped mRFP1 and MelA liquid cultures

- used same elution buffer for two of the same sample to obtain double DNA concentration

Concentrations:

- 1 mRFP: 516 ng/uL
- 2 mRFP: 462.5 ng/uL
- 3 mRFP: 243.5 ng/uL

- 1 MelA: 117.9 ng/uL
- 2 MelA: 158.3 ng/uL
- 3 MelA: 146.1 ng/uL

## **Lab Work 8/13/14:**

### **Restriction Digest:**

Digesting 1 mRFP, 2 mRFP, 3 mRFP, 1 MelA, 2 MelA, 3 MelA with XbaI and PstI.

#### 20 uL Reactions:

5 uL Plasmid DNA  
2 uL FD Green Buffer  
1 uL XbaI  
1 uL PstI  
12 uL ddH<sub>2</sub>O

#### Master Mix:

13 uL FD Green Buffer  
6.5 uL XbaI  
6.5 uL PstI  
78 uL ddH<sub>2</sub>O  
104 uL Total

\*Reactions are 5 uL Plasmid + 15 uL Master Mix

\*There are 6 reactions. Master Mix is 6.5x

\*Incubate for 45 minutes at 37C

### **Gel Electrophoresis:**


Running a gel to check the validity of the mRFP1 and MelA mini-prepped DNA.

Expect to see mRFP1 band at 706 bp and MelA band at 1844 bp.

Lanes:

1 kb+ --- 1 mRFP --- 2 mRFP --- 3 mRFP --- 1 MelA --- 2 MelA --- 3 MelA --- EMPTY

Gel:



**Gel Purification:**

Gel Fragment	Gel + Tube Mass (g)	Tube Mass (g)	Gel Mass (g)	Membrane Binding Solution Added (uL)	DNA Concentration (ng/uL)
1 mRFP	1.142	1.017	0.125		
2 mRFP	1.143	1.009	0.134		
3 mRFP	1.149	1.012	0.137		

\*Did not complete Gel Purification

**Plan for Tomorrow:**

- 1) Finish Gel Purification (Possibly Re-Run. I believe I cut wrong fragment)
- 2) Make chloramphenicol plates
- 3) Resuspend YF1 & FixJ, Blue Light Sensor (Plate 1, Well 10N; 1:10N)
- 4) Resuspend FixK2 Promoter (Plate 1, Well 19G; 1:19G)
- 5) Transform and Plate (YF1 & FixJ) and FixK2 Promoter

If able to obtain Kanamycin:

- 1) Make Kanamycin plates (2-3)
- 2) Transform and Plate mCherry Bomb (3.4 uL of plasmid left)

## **Lab Work 8/14/14:**

### **Plan:**


- 1) Re-Run mRFP1 on Gel and then Gel Purify
- 2) Make Chloramphenicol Plates
- 3) Dilute mCherry Bomb and nanodrop
  - 3a) ~~If enough DNA is present, run PCR (1.4 ng/uL only)~~
  - 3b) If not enough DNA is present, make Kanamycin plates
- 4) Resuspend YF1 & FixJ, Blue Light Sensor (Plate 1, Well 10N; 1:10N)
- 5) Resuspend FixK2 Promoter (Plate 1, Well 19G; 1:19G)
- 6) Transform and Plate (YF1 & FixJ) and FixK2 Promoter (Chloramphenicol)
- 7) Transform and Plate mCherry Bomb plasmid (Kanamycin)

### **mRFP1 Gel Electrophoresis:**

- Create 25 mL or 1% agarose gel and let solidify
- Load digested mRFP1 DNA from 8/13/14

Lanes:

1 kb+ Ladder --- Empty --- 1 mRFP --- empty --- 2 mRFP --- empty --- 3 mRFP --- empty



### **Dilute mCherry Bomb and nanodrop:**

- 1 uL mCherry Bomb into 9 uL ddH<sub>2</sub>O
- Concentration of 1.4 ng/uL

**Make Plates:**

\*1 Kanamycin; 3 Chloramphenicol

-Kanamycin working concentration of 50 ng/mL

-CAM working concentration of 25 ng/mL

**Resuspend DNA:**

\*Added 10 uL ddH<sub>2</sub>O to each part to resuspend

**Transform:**

Transformed Blue Light Sensor, Blue Light Promoter, and mCherry Bomb

\*Plates put in 37C incubator at 4:16 PM

**Gel Purification:**

Gel Fragment	Gel + Tube Mass (g)	Tube Mass (g)	Gel Mass (g)	Membrane Binding Solution Added (uL)	DNA Concentration (ng/uL)
1 mRFP	1.195	1.014	0.181	181	3.6
2 mRFP	1.217	1.002	0.215	215	2.0
3 mRFP	1.158	1.012	0.146	146	4.6



## **Lab Work 8/15/15:**

### **Morning:**

Took out agar plates of Blue Light Sensor, Blue Light Promoter, and mCherry Bomb  
\*Colonies were spotted on each Plate

### **Evening:**

Made liquid cultures for Blue Light Sensor, Blue Light Promoter, and mCherry Bomb  
\*5 mL per liquid culture  
\*2 Liquid cultures per plasmid  
\*Added CAM to Blue Light Sensor and Blue Light Promoter cultures at working concentration of 25 ng/mL  
\*Added Kan to mCherry Bomb culture at working concentration of 50 ng/mL

**Lab Work 8/16/14:**

**Miniprep:**

-Miniprep liquid cultures of Blue Light Promoter, Blue Light Sensor, and mCherry Bomb

Concentrations:

Blue Light Promoter (1:19G): 44.5 ng/uL

Blue Light Sensor (1:10N): 36.6 ng/uL

mCherry Bomb: 35 ng/uL

## Lab Work 8/18/14:

### PCR:

-PCRing hsp60 out of plasmid mCherry Bomb

Master Mix			Primers and DNA			PCR Reaction	
Reagent	1x Rxn (uL)	2.5x Rxn (uL)	Primer	RBS (uL)	no RBS (uL)	Reagent/DNA	Volume (uL)
Buffer	10	25	Fwd Primer	2.5	2.5	Master Mix	11.5
dNTP	1	2.5	Rev Primer	2.5	2.5	Primers and DNA	7.5
Phusion	0.5	1.25	mCherry Bomb	2.5	2.5	ddH2O	31
Total	11.5	28.75	Total	7.5	7.5	Total	50

\*Two 20 uL PCR reactions of each (RBS and no RBS) were carried out

PCR conditions:

98C for 30s

**98C for 10s**

**66C for 30s x30 cycles**

**72C for 30s**

72C for 5 min

4C holding

### Gel Electrophoresis:

Ran two PCR products with a ladder on a gel. There were no bands shown. PCR was unsuccessful

## **Lab Work 8/19/14:**

### **Restriction Digest:**

Digesting Blue Light Sensor and Blue Light Promoter

#### 20 uL Reactions:

5 uL Plasmid DNA  
2 uL FD Green Buffer  
1 uL XbaI  
1 uL PstI  
12 uL ddH<sub>2</sub>O


\*Incubate for 45 minutes at 37C

### **Gel Electrophoresis:**

Running a gel to check the validity of the Blue Light Sensor (1796 bp), Blue Light Promoter (250 bp), and PCR products (~380 bp).

Lanes:

1 kb+ ladder --- RBS PCR --- no RBS PCR --- Blue Light Sensor --- Blue Light Sensor --- Blue Light Promoter --- Blue Light Promoter --- empty



9/12/14

Goal:

1. PCR purify Hsp60 with and without RBS
2. Digestion of purified product with XbaI and PstI
3. Gel purification of digestion products
  - a. Cut out slices and store at 4 °C

Experiments:

Tubes Labeled by MJ

Tube label	Content
1, PCR purification product, no RBS	1 <sup>st</sup> elution: PCR purification product, no RBS
2, PCR purification product, no RBS	2 <sup>nd</sup> elution: PCR purification product, no RBS
1, PCR purification product, RBS	1 <sup>st</sup> elution: PCR purification product, RBS
2, PCR purification product, RBS	2 <sup>nd</sup> elution: PCR purification product, RBS

1. PCR purification

Step	Hsp60 with RBS (uL)	Components	Hsp60 without RBS (uL)
1	34.5	Membrane binding solution	34.5
2	700	Membrane Wash solution	700
3	500	Membrane Wash solution	500
4	50	1 <sup>st</sup> elution with Nuclease free ddH <sub>2</sub> O	50
5	30	2 <sup>nd</sup> elution with Nuclease free ddH <sub>2</sub> O	30

2. Digestion of purified product

Components	Amount for Hsp60 with RBS (uL)	Amount for Hsp60 without RBS (uL)
10x Fast Digest Green Buffer	2	2
PCR Product of Hsp60	7	7
XbaI	1	1
PstI	1	1
ddH <sub>2</sub> O	13	13
Total	24	24

Incubated at 37 °C for 30 minutes

2 reactions each were digested for 1<sup>st</sup> elution of both with and without RBS; 1 reaction each for the 2<sup>nd</sup> elution for both with and without RBS


9/12/14 continued

3. Gel Set Up:

Lane	Sample	Amount (uL)	Expected band size (bp)
1	1kb DNA Ladder	2	None
2	Elution 1, no RBS digested	5	~ 300 bp
3	Elution 1, no RBS digested	5	
4	Elution 2, no RBS digested	5	
5	1kb DNA Ladder	2	None
6	Elution 1, RBS digested	5	~300 bp
7	Elution 1, RBS digested	5	
8	Elution 2, RBS digested	5	

4. Mass of tube and tube plus sample

Tube label	Sample (lane #)	Mass of tube (g)	Mass of Tube plus sample (g)
1	2	1.012	1.1
2	3	1.005	1.098
3	4	1.018	1.121
4	6	1.01	1.105
5	7	1.02	1.113
6	8	1.003	1.094



## **Lab Work 9/17/14:**

### **Goals:**

- 1) Digestion of Blue Light Parts w/ X+P
- 2) Run Parts on Gel
- 3) Gel Purify Blue Light Parts, hsp60 parts, vector backbone
- 4) Ligate hsp60 + backbone

### **Digestion:**

9 uL ddH<sub>2</sub>O  
2 uL Fast Digest Green Buffer  
1 uL XbaI  
1 uL PstI  
7 uL Plasmid (Blue Light Sensor and Blue Light Promoter)  
20 uL rxn

Incubate at 37 C for 30 min

### **Gel Electrophoresis:**

Lanes:

1 kb+ ladder – Blue Sensor – Blue Sensor – Blue Promoter – Blue Promoter – empty

Image at 30 mins:





Image at 1 hour, 15 mins:



\*Lanes 4/5 extracted at 30 mins

**Purification:**

	Tube Weight (g)	Tube + Gel Weight (g)	Gel Fragment Weight (g)	Membrane Binding Solution (uL)
Backbone	1.021	1.183	0.162	162
Blue Promoter	1.016	1.198	0.182	182
Backbone	1.016	1.187	0.171	171
Blue Sensor	1.004	1.126	0.122	122
no RBS			0.088	88
no RBS			0.093	93
RBS			0.095	95
RBS			0.093	93

**Ligation:**

RBS

1 uL ligase  
2 uL Buffer  
4 uL Insert  
3.5 uL Backbone  
9.5 uL ddH<sub>2</sub>O  
20 uL

no RBS

1 uL ligase  
2 uL Buffer  
3 uL Insert  
3.5 uL Backbone  
10.5 uL ddH<sub>2</sub>O  
20 uL



9/22/14

**Goals:**

1. Digestion of desaturase, and mRFP-1
2. Purification of desaturase and mRFP-1
3. Ligation into the vector
4. Digestion of 7 uL of plasmid DNA
5. Purification

**Experiment:**


Cut mRPF-1 with EcoRI and XbaI

Desaturase with EcoRI and SpeI

Components mRFP-1	Amount (uL)	Components desaturase
10x Fast Digest Buffer	2	10x Fast Digest Buffer
plasmid	7	plasmid
EcoRI	1	EcoRI
XbaI	1	SpeI
ddH <sub>2</sub> O	9	ddH <sub>2</sub> O
Total	19	Total

**Gel:**

Lane	1	2	3	4
Components	1kb DNA Ladder	Desaturase (EcoRI+SpeI)		
Amount (uL)	2	5	5	5
Lane	5	6	7	8
Components	1kb DNA Ladder	mRFP-1 (EcoRI+XbaI)		
Amount (uL)	2	5	5	5



9/24/14

mRFP-1 and desaturase plasmids were previously collected and this is the resulting digestion (9/22/14) followed by gel purification of the previous experiments.

1. Gel purification: and resulting concentrations:

Tube label	Part excised	Mass (g)	Concentration (ng/uL)
1	Desaturase Vector Backbone	Not used since this is the backbone	
2			
3			
4	Desaturase Part	0.124	3.9
5		0.123	4.7
6		0.095	4.2
7	mRFP-1 linearized plasmid	0.095	4.4
8		0.103	4.4
9		0.133	4.8

9/26/14

1. Ligation reaction BH

Components	Amount (uL)
Insert (desaturase)	9
Vector Backbone (single cut mRFP-1)	3.5
Ligase	1
Buffer	1.5
ddH <sub>2</sub> O	5
Total	20

Negative control: just VB (labeled A)

Ligation 1 (labeled B)

Ligation 2 (labeled C)

2. Transformation into E. coli competent cells BH

Incubation at 37 °C on LB + CAM from 2 pm

Note: the plate was examined once at 6 hours and another time at 8 hours—  
yielded no usable colonies

3. Mini-Prep (BH) of colonies culture grown by MJ

Labeling maintained the same (1-6)

Resuspended in 50 uL of elution buffer

Undetermined concentration

Ligation and Transformation were repeated 2 times yet the reaction never yielded colonies that survived in LB+CAM culture.

**Lab Work 9/29/14:**

**Diagnostic Digest:**

Hsp60(no RBS) – mRFP1 (~1090 bp) 1) X+P

Blue Promoter – mRFP1 (~1050 bp) 2) X+P, 3) S+P

Hsp60 (no RBS) – Blue Sensor (~2160 bp) 4) X+P

MelA (~1896 bp) 5) X+P

20 uL Rxns

9 uL ddH<sub>2</sub>O

2 uL FD Green

1 uL XbaI/SpeI

1 uL PstI


7 uL Plasmid

\*Incubate for 30 min @ 37C


**Gel Electrophoresis:**

1 kb+ ladder – 1 – 2 – 3 – 5 – 4 – empty

Gel @ 30 mins



Gel @ 1 hr, 30 mins



\*Note: Band 3 was extracted after 30 minutes. The lower of band 5 was extracted at 1 hr, 30 min

**Gel Extraction:**

Blue Promoter – mRFP1 S+P gel weight: 0.086 g

Hsp60 (no RBS) – Blue Sensor X+P gel weight: 0.062 g

**Made 6 CAM plates & 2 AMP plates**

**Transformation:**

Transformed RBS (4:1N) and Terminator (4:16G).  
Plate onto AMP plates

## **Lab Work 10/1/14:**

### **Mini Prep:**

Mini Prep RBS and Terminator Liquid Cultures

### **Gel Purify:**

Add 86 uL Membrane Binding Solution to Blue Promoter – mRFP1 (S+P) part

Add 62 uL Membrane Binding Solution to hsp60 (no RBS) – Blue Sensor (X+P) part

### **Restriction Digests:**

Terminator (E+X) (~2150 bp)

Cathelicidin (E+S) (~123 bp)

RBS (E+X) (~2082 bp)

Blue Promoter (E+S) (~250 bp)

#### 20 uL rxns:

9 uL ddH<sub>2</sub>O

2 uL FD Green

1 uL EcoRI

1 uL XbaI/SpeI

7 uL Plasmid


Incubate 30 mins @ 37C

### **Gel Electrophoresis:**

Lanes:

1 kb+ ladder – Terminator (E+X) – Cathelicidin (E+S) – RBS (E+X) – B. Promoter (E+S) – empty

\*Run Gel 15 minutes



**Gel Purify:**

RBS (E+X)

Blue Promoter (E+S)

Gel Fragment	Gel + Tube Mass (g)	Tube mass (g)	Gel Mass (g)	Membrane Binding Solution (uL)	DNA Concentration (ng/uL)
<del>Terminator (E+X)</del>	-	<del>1.106</del>	-	-	-
<del>Cathelicidin (E+S)</del>	-	<del>1.015</del>	-	-	-
RBS (E+X)	1.033	0.985	0.048	48	1.7
Blue Promoter (E+S)	1.082	1.005	0.077	77	1.8

**Ligation:**

Blue Promoter – RBS

15 uL Rxn:

3.5 uL Insert (Blue Promoter)

9 uL VB (RBS)

1.5 uL Buffer

1 uL Ligase

0 uL ddH<sub>2</sub>O**Transformation:**

Transform Blue Promoter – RBS onto AMP plate

## Lab Work 10/4/14:

Morning:

Made Liquid Cultures of Blue Promoter – RBS part in Turbo E.coli Comp Cells

Afternoon:

### Mini Prep:

Mini Prepped Blue Promoter – RBS : ??? ng/uL

### Restriction Digest:

Blue Promoter – RBS (X/P)

Blue Promoter – RBS (S/P)

Terminator (S/P)

20 uL rxns:

9 uL ddH<sub>2</sub>O

2 uL FD Green

1 uL XbaI/SpeI


1 uL PstI

7 uL Plasmid

### Gel Electrophoresis:

Lanes:

- 1) 1 kb+ ladder
- 2) Blue Promoter – RBS (X+P) (~262 bp)
- 3) Blue Promoter – RBS (S+P) (~2350 bp)
- 4) Terminator (S+P) (~2150 bp)



### Gel Purify:

Gel Fragment	Tube (g)	Gel + Tube (g)	Gel (g)	Membrane Binding Solution (uL)	Concentration (ng/uL)
Blue Promoter-RBS (S/P)	1.002	1.067	0.065	65	5.0
Terminator (S/P)	1.020	1.093	0.073	73	2.1

### Ligation:



- 1) Blue Promoter – mRFP1 (S/P) (~3000 bp) + Cathelicidin (X/P) (~120 bp)
- 2) Blue Promoter – RBS (S/P) (~2350 bp) + Cathelicidin (X/P) (~120 bp)
- 3) Terminator (S/P) (~2150 bp) + hsp60 (no RBS) – Blue Sensor (X/P) (~2160 bp)
- 4) Blue Promoter – mRFP1 (S/P) (~3000 bp) + hsp60 (no RBS) – Blue Sensor (X/P) (~2160 bp)

Ligation	VB bp	VB concentration (ng/uL)	Insert bp	Insert concentration (ng/uL)	VB (uL)	Insert (uL)
1	3000	2.9	120	3.9	4.00	8.50
2	2350	5.0	120	3.9	3.50	9.00
3	2150	2.1	2160	4.9	3.00	9.50
4	3000	2.9	2160	4.9	3.50	9.00

1.5 uL Buffer

1 uL Ligase

**Transformation:**

Blue Promoter – mRFP1 – Cathelicidin (CAM)

Blue Promoter – RBS – Cathelicidin (AMP)

Terminator – hsp60 (no RBS) – Blue Sensor (AMP)

Blue Promoter – mRFP1 – hsp60 (no RBS) – Blue Sensor (CAM)

## **Lab Work 10/5/14:**

### **Ligation:**

Ligation 3: Terminator (S/P) + hsp60 (no RBS) – Blue Sensor (X/P)

1.5 uL Buffer

1 uL Ligase

6 uL VB (Terminator)

6.5 uL Insert (hsp60-Blue Sensor)

### **Transformation:**

Transform 3: Terminator – hsp60 (no RBS) – Blue Sensor

AMP Plate

### **MiniPrep:**

Mini 1: Blue Promoter – mRFP1 – Cathelicidin      61.5 ng/uL

Mini 2: Blue Promoter – RBS – Cathelicidin      85.9 ng/uL

Mini 4: Blue Promoter – mRFP1 – hsp60 (no RBS) – Blue Sensor      81.0 ng/uL

### **Restriction Digest:**

- 1) Mini 1: Blue Promoter – mRFP1 – Cathelicidin (X/P)
- 2) Mini 1: Blue Promoter – mRFP1 – Cathelicidin (S/P)
- 3) Mini 2: Blue Promoter – RBS – Cathelicidin (X/P)
- 4) Mini 2: Blue Promoter – RBS – Cathelicidin (S/P)
- 5) Mini 4: Blue Promoter – mRFP1 – hsp60 (no RBS) – Blue Sensor (X/P)

20 uL rxns:

9 uL ddH<sub>2</sub>O

2 uL FD Green

1 uL XbaI/SpeI


1 uL PstI

7 uL Plasmid


## Gel Electrophoresis:

Lanes:

- 1) 1 kb+ ladder
- 2) Mini 1: Blue Promoter – mRFP1 – Cathelicidin (X/P) (~1000 bp) (Part)
- 3) Mini 1: Blue Promoter – mRFP1 – Cathelicidin (S/P) (~3100 bp) (VB)
- 4) Mini 2: Blue Promoter – RBS – Cathelicidin (X/P) (~370 bp) (Part)
- 5) Mini 2: Blue Promoter – RBS – Cathelicidin (S/P) (~2470 bp) (VB)
- 6) Mini 4: Blue Promoter – mRFP1 – hsp60 (no RBS) – Blue Sensor (X/P)(~3160 bp) (Part)



\*After 30 mins, hsp60 (no RBS) – Blue Sensor (X/P) (~2160 bp) added to lane 7. Ladder added lane 8.



**Gel Purify:**

Gel Fragment	Tube (g)	Gel + Tube (g)	Gel (g)	Membrane Binding Solution (uL)	Concentration (g/uL)
Mini 1: B.P. - mRFP1 - Cath (S/P)	1.013	1.072	0.059	59	
no RBS - Blue Sensor (X/P)	1.016	1.065	0.049	49	

**Ligation:**

1.5 uL Ligation Buffer

1 uL Ligase

Lig 1:

9 uL Blue Promoter (Insert) (~200 bp)

3.5 uL RBS (VB) (~2100 bp)

Lig 2:

6 uL Terminator (VB) (~2100 bp)

6.5 uL hsp60 (no RBS) – B. Sensor (Insert) (~2160)

Lig 3:

3.5 uL B. Prom – mRFP1 (VB) (~3100)

9 uL hsp60 (no RBS) – B. Sensor (Insert) (~2160 bp)

**Transformation:**

Transform 1: Blue Promoter – RBS onto AMP plate

Transform 2: Terminator – hsp60 (no RBS) – Blue Sensor onto AMP plate

Transform 3: Blue Promoter – mRFP1 – hsp60 (no RBS) – Blue Sensor onto CAM plate

**Make Liquid Cultures:**

Made 2 5mL Liquid cultures of LB + AMP, Blue Promoter – RBS

## **Lab Work 10/6/14:**

### **Make Liquid Cultures:**

(AMP) Transform 1: Blue Promoter – RBS

(AMP) Transform 2: Terminator – hsp60 (no RBS) – Blue Sensor

(CAM) Transform 3: Blue Promoter – mRFP1 – hsp60 (no RBS) – Blue Sensor

### **Mini prep:**

Mini 1: Blue Promoter – RBS

Mini 2: Terminator – hsp60 (no RBS) – Blue Sensor

Mini 3: Blue Promoter – mRFP1 – hsp60 (no RBS) – Blue Sensor

### **Restriction Digest:**

- 1) Blue Promoter – RBS (X/P)
- 2) Blue Promoter – RBS (S/P)
- 3) Terminator – hsp60 (no RBS) – Blue Sensor (X/P)
- 4) Blue Promoter – mRFP1 – hsp60 (no RBS) – Blue Sensor (X/P)

### **20 uL Rxns**

9 uL ddH<sub>2</sub>O

2 uL FD Green

1 uL XbaI/SpeI

1 uL PstI

7 uL Plasmid

### **Gel Electrophoresis:**



Lanes:

- 1) 1 kb+ ladder
- 2) Blue Promoter – RBS (X/P) (~262 bp)
- 3) Blue Promoter – RBS (S/P) (~2341)
- 4) Terminator – hsp60 (no RBS) – Blue Sensor (X/P) (~2236)
- 5) Blue Promoter – mRFP1 – hsp60 (no RBS) – Blue Sensor (X/P) (~3100)

\*Note: In lane four, the VB is ~2100, will need long separation to see. Part will be bigger than backbone

\*Stop at 15 mins and check lane 2, if correct band is present, extract lane 3 before continuing.

\*Run for another hour and fifteen before stopping a second time and extracting lane 4



**Gel Purify:**

**Ligation:**

Ligation 1: Blue Promoter – RBS (X/P) + Cathelicidin (S/P)

Ligation 2: Blue Promoter – mRFP1 – Cathelicidin (S/P) + Terminator – hsp60 (no RBS) – Blue Sensor (X/P)

Ligation	VB bp	VB Concentration (ng/uL)	Insert bp	Insert Concentration (ng/uL)	VB (uL)	Insert (uL)
1	2341	-	123	-	-	-
2	2500		2236			