

iGEM2014 – Microbiology – BMB – SDU	
Title:	Date issued: 2013.06.19
SOP number: SOP0015	Review date:
Version number: 02	Written by: ASF

1. Purpose

To ligate pieces of DNA

2. Area of application

Cloning

3. Apparatus and equipment

Apparatus/equipment	Location (Room number)	Check points	Criteria for approval/rejection
Vortex		•	
Pipettes (p20, p10)		•	
		•	
		•	
		•	
		•	
		•	

4. Materials and reagents – their shelf life and risk labelling

Name	Components	Supplier / Cat. #	Room (hallway storage)	Safety considerations
Purple pipette tips		Contact lab-manager	Micro storage	
Green pipette tips		Contact lab-manager	Micro storage	
Eppendorftubes		Contact lab-manager	Micro storage	

Distilled water		Contact lab-manager	Micro storage	
Ligasebuffer		Agilent Technologies	Freezer at 1. Floor	
Ligase				
DNA piece 1				
DNA piece 2				

5. QC – Quality Control**6. List of other SOPs relevant to this SOP****7. Environmental conditions required****8. Procedure**

1. Prepare the ligation mixture and mix by pipetting up and down
2. Leave the mixture overnight at 16°C
- 2a. If there is no time leave the ligation solution at 22.5°C for a minimum of 2 hours.
3. Use ligation solution for transformations but save ¼ of the solution in the fridge so that a second transformation can be done in case something goes wrong.

Reagents	Volume
10x T4 DNA ligase buffer	2 µL
T4 DNA ligase (add last!)	1 µL
PCR product (cut) of each brick which is to be ligated – or 1 part plasmid and 5 part bricks	5 µL or 10 fmol Plasmid, 0, 10 and 20 fmol PCR
H2O	to reach a total volume of 20µL

9. Waste handling

Chemical name	Concentration	Type of waste (C, Z...)	Remarks

10. Time consumption

11. Scheme of development

Date / Initials	Version No.	Description of changes
13.06.19 / ASF	01	The SOP has been written
13.06.26 / PRA	01	The SOP has been approved
14.07.09 / DWP	02	The SOP has been updated

12. Appendixes