

Reagents

Electrophoresis buffer (50x): 2M Tris-acetate, pH 8.5, 100 mM EDTA

KqGreen safe DNA dye (20,000x in water) (Shanghai Keqing Biotech)

DSTM5000 molecular weight marker (Dongsheng Biotech)

The sizes of the bands (base pairs) in the molecular weight marker sample are:

5,000

3,000

2,000

1,500

1,000 (index strip)

750

500

250

100

A. Preparation of the gel

1. 1.0g agarose was made up to 100mL with electrophoresis buffer.
2. The solution was heated in a microwave until it turned clear.
3. After the solution had cooled to about 50-60°C, 4uL of a concentrated stock of KqGreen safe DNA dye was added.
4. The mixture was gently swirled and poured into the gel tray containing the comb.
5. After the gel had set (approx 30 minutes) the comb was taken out.
6. The tape sealing the ends of the gel tray was removed and the gel tray placed in position in the tank.

B. Loading the samples and electrophoresis

1. Add sufficient electrophoresis buffer to completely cover the gel by about 2 mm.
2. Using an automatic pipette, add 3uL of each DNA sample or the molecular weight marker to separate wells. Use a fresh tip each time.
3. When all samples have been loaded, connect the leads to the power supply, making sure that the lead connected to the negative (cathode) is nearest the wells, and run at 160 V for 15 minutes.
4. After electrophoresis is complete, switch off the power pack. Take out the gel tray; take the gel to the UV gel documentation system, scan and save the picture. Analyze the positions of the bands in the different lanes.